A randomised comparison of healing response between the BuMA Supreme stent and the XIENCE stent at one-month and two-month follow-up: PIONEER-II OCT randomised controlled trial



Taku Asano^{1,2}, MD; Qinhua Jin³, MD; Yuki Katagiri¹, MD; Norihiro Kogame¹, MD; Kuniaki Takahashi¹, MD; Chun-Chin Chang⁴, MD; Ply Chichareon¹, MD; ChangQian Wang⁵, MD; Bei Shi⁶, MD; Xi Su⁻, MD; GuoSheng Fu⁶, MD; YanQing Wu⁶, MD; XuChen Zhou¹⁰, MD; ZuYi Yuan¹¹, MD; Joanna J. Wykrzykowska¹, MD, PhD; Jan J. Piek¹, MD, PhD; Patrick W. Serruys¹², MD, PhD; Yoshinobu Onuma⁴,¹³, MD, PhD; YunDai Chen³*, MD

1. Academic Medical Center, Amsterdam, the Netherlands; 2. St. Luke's International Hospital, Tokyo, Japan; 3. The General Hospital of Chinese People's Liberation Army, Beijing, China; 4. Erasmus Medical Center, Rotterdam, the Netherlands; 5. Shanghai 9th People's Hospital, Shanghai, China; 6. Affiliated Hospital of Zunyi Medical College, Guizhou Sheng, China; 7. Wuhan Asia Heart Hospital, Wuhan, China; 8. Sir Run Run Shaw Hospital, Zhejiang University School of Medicine, Zhejiang, China; 9. 2nd Affiliated Hospital of Nanchang University, Nanchang, China; 10. 1st Affiliated Hospital of Dalian Medical University, Liaoning Sheng, China; 11. 1st Affiliated Hospital of Xi'an Jiaotong University, Shaanxi Sheng, China; 12. Imperial

T. Asano and Q. Jin contributed equally to this article.

GUEST EDITOR: Antonio Colombo, MD, FACC, FESC, FSCAI; Department of Interventional Cardiology, EMO GVM Centro Cuore Columbus, Milan, Italy

This paper also includes supplementary data published online at: http://www.pcronline.com/eurointervention/145th_issue/236

KEYWORDS

- drug-eluting stentoptical coherence
- tomography
 stable angina

Abstract

College London, London, United Kingdom; 13. Cardialysis, Rotterdam, the Netherlands

Aims: The aim of this study was to compare the strut coverage of the XIENCE stent with that of the BuMA Supreme sirolimus-eluting cobalt-chromium stent, which has a shorter drug elution, on optical coherence tomography (OCT) one or two months after implantation.

Methods and results: The PIONEER-II OCT trial was a multicentre, two-arm randomised trial, which comprised two cohorts: cohort-1 underwent an OCT imaging one month after coronary intervention (BuMA: 16 patients with 18 lesions, XIENCE: 15 patients with 17 lesions), whereas cohort-2 underwent OCT at two months (BuMA: 21 patients with 21 lesions, XIENCE: 23 patients with 28 lesions). The primary hypotheses were non-inferiority of the BuMA stent to the XIENCE stent in percent strut coverage at one month (cohort-1) or two months (cohort-2). In cohort-1, the BuMA stent was non-inferior to the XIENCE stent in terms of the strut coverage (83.8 \pm 10.4% for BuMA vs. 73.0 \pm 17.5% for XIENCE, p_{for noninferiority} <0.001), and was also significantly higher than the XIENCE (p_{for superiority} 0.037). In cohort-2, the BuMA stent was non-inferior to the XIENCE stent in OCT strut coverage (80.3 \pm 18.3% vs. 73.3 \pm 21.3%, p_{for noninferiority} 0.006, p_{for superiority} 0.24). Healing scores showed better healing in the BuMA stent in cohort-1 (32.36 \pm 21.59 vs. 54.88 \pm 34.65, p=0.027), whereas there was comparable healing between the BuMA and XIENCE stents in cohort-2 (39.86 \pm 37.77 vs. 53.75 \pm 42.84, p=0.25).

Conclusions: The BuMA Supreme had a faster coverage than the XIENCE at one month, presumably due to faster and shorter sirolimus elution. The difference in tissue coverage became less evident at two months. ClinicalTrials.gov Identifier: NCT02747329.

^{*}Corresponding author: Department of Cardiology, The General Hospital of Chinese People's Liberation Army, Beijing, 100853, China. E-mail: cyundai@vip.163.com

Abbreviations

LLL

%DS percent diameter stenosis
 BP biodegradable polymer
 CI confidence interval
 DAPT dual antiplatelet therapy
 DES drug-eluting stent
 HS healing score

late lumen loss

MLD minimum lumen diameter

OCT optical coherence tomography

PBMA poly(n-butyl methacrylate)

PLGA poly(lactic-co-glycolic acid)

QCA quantitative coronary angiography

RVD reference vessel diameter
SES sirolimus-eluting stent
ZES zotarolimus-eluting stent

Introduction

The second generation of drug-eluting stents (DES) has emerged with more biocompatible durable polymer coatings on thinstrut stents that have shown a better safety than first-generation DES^{1,2}. Biodegradable polymer (BP)-DES were developed to prevent polymer-related adverse effects³. The BuMA Supreme sirolimus-eluting stent (SES) (SINOMED, Tianjin, China) consists of a thin-strut (80 µm) cobalt-chromium platform with a sub-micrometre-thin electro-grafted base layer, top-coated with a biodegradable coating that releases sirolimus relatively fast within 28 days. With the short drug elution time, the first-inman trial of the BuMA Supreme stent showed a greater ninemonth angiographic late loss of 0.29±0.33 mm as compared to 0.14±0.37 mm in the Resolute™ zotarolimus-eluting stent (ZES) (Medtronic, Minneapolis, MN, USA), although the additional analysis showed the absence of difference in quantitative angiography-derived FFR between the two devices^{4,5}. This moderate late loss value could be the result of an early coverage of the BuMA Supreme stent.

The degradation of the top-layer poly-glyco-lactide acid coating of the BuMA Supreme stent is complete in six weeks⁶. To date, the advantage of BP-DES has been debatable. Some studies have shown improvement of clinical outcomes with BP-DES as compared to durable polymer DES, but the number of those studies is limited^{3,7}. Recent studies have reported less advantage of BP-DES compared with second-generation DES with durable polymer^{8,9}.

We hypothesised that the BuMA Supreme stent with its short elution time could achieve earlier strut coverage than the XIENCE stent (Abbott Vascular, Santa Clara, CA, USA) which has a relatively longer antiproliferative drug elution and a durable polymer coating. In this PIONEER-II OCT randomised trial comparing BuMA and XIENCE, the strut coverage was assessed on optical coherence tomography (OCT) one and two months after implantation.

Editorial, see page 1247

Methods

The PIONEER-II OCT trial was a multicentre, two-arm randomised trial conducted in China, which comprised two cohorts: cohort-1 underwent an OCT imaging one month after coronary intervention, whereas cohort-2 underwent OCT at two months. In the current trial, 31 patients were assigned to cohort-1 (BuMA: n=16, XIENCE: n=15), whereas 44 patients were enrolled to cohort-2 (BuMA: n=21, XIENCE: n=23). The current OCT study was also designed to explore the early-phase strut coverage in the patients at high bleeding risk. Details of the eligibility criteria of the current study are presented in Supplementary Appendix 1. Briefly, the inclusion criteria were: (1) age 18-85; (2) presence of stable or unstable angina or silent ischaemia; (3) presence of one or two separate de novo target lesions (a single target lesion per major epicardial territory) with a 50-99% lumen diameter stenosis in a reference vessel of 2.5-4.5 mm (visually determined), and a target lesion length of <40 mm, and (4) presence of one or more high bleeding risks such as a history of prior intracerebral bleeding, presence of renal failure (calculated creatinine clearance <40 ml/min) or presence of non-skin cancer diagnosed or treated ≤3 years prior to the enrolment.

Eligible patients were randomly assigned to each cohort (cohort 1: one-month OCT follow-up or cohort 2: two-month OCT follow-up) and assigned to receive the BuMA Supreme SES or XIENCE everolimus-eluting stent family (the control device) in a 1:1 ratio in random blocks, stratified by centre using a web-based allocation system. The current study complied with the Declaration of Helsinki. Medical ethics committees at each participating institution approved the trial, and all patients provided a written informed consent.

STUDY DEVICES AND IMPLANTATION PROCEDURE

The BuMA Supreme SES is a DES with a rapid exchange balloon-expandable catheter delivery system. The device consists of an extremely thin (100-200 nm) electrografted (eG) base layer of poly(n-butyl methacrylate) (PBMA), covalently bound to an electronically polished thin-strut (80 μ) cobalt-chromium stent platform. The base layer's molecules are on their other side connected by interdigitation with a 3.8 to 10 μ m-thin top layer – a blend of a biodegradable poly(lactic-co-glycolic acid) (PLGA) polymer and sirolimus. This base layer prevents the bioactive coating from cracking and delamination during delivery and expansion of the stent⁶. Unlike other DES, the BuMA Supreme has a fast drug elution profile: it elutes 92% of sirolimus in the first 28 days and 100% at two months (**Figure 1**)⁷. The XIENCE family of stents was used in the control group. These DES are made from cobalt-chromium (81 μ m) with a permanent fluoropolymer eluting everolimus¹⁰.

The interventional procedure was performed and dual antiplatelet therapy (DAPT) was given according to current clinical guidelines¹¹. The patients underwent repeat angiography and OCT at either one month (cohort-1) or two months (cohort-2).

QUANTITATIVE CORONARY ANGIOGRAPHY (QCA)

Off-line QCA analyses were performed by an independent core laboratory (Cardialysis BV, Rotterdam, the Netherlands) with the




Figure 1. Drug release curves of the BuMA Supreme and XIENCE stents

CAAS system (Pie Medical Imaging, Maastricht, the Netherlands) according to standard operational protocols¹². The following QCA parameters were calculated both in-stent and in-segment (including the 5 mm proximal and distal stent margins): minimum lumen diameter (MLD), interpolated reference vessel diameter (RVD), percent diameter stenosis (%DS: [1-MLD/RVD]×100), acute lumen gain (difference in MLD between pre-procedure and post-procedure) and late lumen loss (LLL: difference in MLD between post-procedure and follow-up).

OCT ANALYSIS

OCT assessment of the stented coronary segment was performed using the final OCT recordings sent to an independent core laboratory (Cardialysis BV) for off-line analysis using QIvus software, version 3.0 (Medis, Leiden, the Netherlands). A covered strut was defined as having a neointimal thickness >0 µm¹³. For the secondary endpoint analysis, stent area and derived measures were based on the endoluminal stent contour¹⁴. To quantify the degree of vascular healing status, the healing score (HS) was calculated according to the previously published methods¹⁵. Briefly, the HS was calculated by summing up the points from four components of OCT findings with specific weights: presence of intraluminal mass (assigned a weight of 4); presence of both malapposed and uncovered struts (assigned a weight of 2); and presence of malapposed struts alone (assigned a weight of 1).

ENDPOINTS

The primary study endpoint of this imaging study was the percentage strut coverage on OCT. There were two independent hypotheses for each cohort, which were: 1) for cohort-1, non-inferiority of the BuMA Supreme stent to the XIENCE stent for the primary endpoint at one month; 2) for cohort-2, non-inferiority of the BuMA Supreme stent to the XIENCE stent for the primary endpoint at two months. OCT imaging endpoints included the mean and minimal scaffold/stent and diameter, area and volume,

the frequency of incomplete strut apposition including area and volume, the percentage of malapposed struts, the mean neointima thickness and neointimal hyperplasia area on top of the strut, and volume, the mean flow area and volume, the malapposition area and volume, and intraluminal defect area and volume.

STATISTICAL METHODS AND SAMPLE SIZE CALCULATION

For binary variables, counts, percentages, and 95% confidence intervals (CI) were calculated. Pearson's chi-squared test or Fisher's exact test was performed when appropriate. For continuous variables, means, standard deviations, and 95% CI for the mean were calculated and compared with a t-test. A two-sided p-value <0.05 was considered statistically significant. All statistical analyses were performed using SAS version 9.3 (SAS Institute Inc, Cary, NC, USA) and SPSS, Version 24.0.0 (IBM Corp., Armonk, NY, USA).

Multilevel linear regression analyses with random intercept and fixed slope model were used in order to investigate the impact of the usage of the BuMA stent on percentage of strut coverage while adjusting for baseline and procedural characteristics selected based on previously known determinants associated with neointimal hyperplasia or strut coverage (diabetes, unstable angina pectoris, post-dilatation)¹⁶⁻²⁰. In these multilevel models, lesions were clustered within a patient when the patient had multiple lesions treated.

For cohort-1, the sample size was based on the following assumptions: based on the previous literature, the coverage rate of the XIENCE at one month was 70% and the BuMA 71.6%²¹. With a pooled standard deviation of 19%, one-sided alpha of 0.05 and a non-inferiority margin of 15%, 20 OCT per arm were needed to prove non-inferiority of BuMA to XIENCE with an 85% power. For cohort-2, the coverage rate of XIENCE and BuMA at two months was assumed to be 83% and 86.4%, respectively^{22,23}. With a pooled standard deviation of 11%, one-sided alpha of 0.05 and a non-inferiority margin of 8.5%, 20 OCT per arm were required to prove non-inferiority of BuMA to XIENCE with a 95% power.

Results

From May 2016 to March 2017, 75 patients were recruited at eight national sites in China. The baseline characteristics of the two cohorts are presented in **Supplementary Table 1**. There were no significant differences in the two device groups, though there existed a non-significant trend towards higher incidence of diabetes in the XIENCE arm of cohort-2. The patient flow chart is shown in **Figure 2**. In cohort-1, all patients received DAPT at one-month follow-up in both arms. In cohort-2, 90.5% (19/21) of the patients in the BuMA arm received DAPT at two-month follow-up while all patients received DAPT in the XIENCE arm. All patients enrolled in the current study continued statin at follow-up. Details of the bleeding propensity of the enrolled patients are presented in **Supplementary Table 2**.

QCA AND OCT RESULTS AT ONE MONTH (COHORT-1)

The results of angiographic and OCT analysis are presented in **Supplementary Table 3** and **Supplementary Table 4**, respectively.




Figure 2. Study flow chart of the PIONEER-II OCT randomised controlled trial. L: number of lesions; N: number of patients

Post-procedural QCA results were comparable between the two groups. At one month, there was no binary restenosis and the in-stent LLL of the two devices was not statistically different (0.12 \pm 0.25 mm for BuMA vs. 0.03 \pm 0.20 mm for XIENCE, p=0.25).

On OCT at one month, the BuMA Supreme stent was non-inferior to the XIENCE stent in terms of strut coverage as assessed on OCT (83.8±10.4% for BuMA vs. 73.0±17.5% for XIENCE, $p_{for\ noninferiority}$ <0.001), whereas the one-month OCT coverage rate of the BuMA Supreme was significantly higher than the XIENCE ($p_{for\ superiority}$ 0.037) (**Figure 3A**). Mean neointimal hyperplasia area was greater in the BuMA than in the XIENCE (0.63±0.37 mm² vs. 0.36±0.12 mm², p=0.008) as was volume obstruction (7.5±4.9%

vs. 4.9±1.7%, p=0.045). The healing score was lower in the BuMA than in the XIENCE stent (32.36±21.59 vs. 54.88±34.65, p=0.027), suggesting earlier healing with the BuMA in terms of coverage and malapposition (Figure 4A). Representative cases in cohort-1 are shown in Figure 5A. The quantitative parameters did not differ between the two arms: minimum lumen area was 6.43±1.39 mm² and 5.80±2.12 mm² in the BuMA arm and the XIENCE arm, respectively (p=0.30).

QCA AND OCT AT TWO MONTHS (COHORT-2)

Quantitative assessment of angiography and OCT is shown in **Supplementary Table 3** and **Supplementary Table 4**, respectively. Angiographic LLL at two months was 0.06 ± 0.21 mm in the BuMA




Figure 3. Cumulative frequency distribution curves of strut coverage as assessed on OCT in cohort-1 (A) and cohort-2 (B). L: number of lesions




Figure 4. Cumulative frequency distribution curves of healing score as assessed on OCT in cohort-1 (A) and cohort-2 (B). L: number of lesions

arm and 0.06 ± 0.22 mm in the XIENCE arm. There was no binary restenosis. On OCT at two months, the BuMA Supreme stent was non-inferior to the XIENCE stent in OCT strut coverage ($80.3\pm18.3\%$ for BuMA vs. $73.3\pm21.3\%$ for XIENCE, $p_{for\ noninferiority}$ 0.006, $p_{for\ superiority}$ 0.24) (**Figure 3B**). Healing score was comparable between the two groups (39.86 ± 37.77 vs. 53.75 ± 42.84 , p=0.25) (**Figure 4B**). Representative cases in cohort-1 are presented in **Figure 5B**. The quantitative parameters did not differ between the two arms: minimum lumen area was 6.47 ± 2.51 mm² and 6.78 ± 2.26 mm² in the BuMA arm and the XIENCE arm, respectively (p=0.65).

MULTILEVEL LINEAR REGRESSION ANALYSIS INVESTIGATING THE IMPACT OF THE BUMA STENT ON STRUT COVERAGE

The results of the regression analyses are summarised in **Table 1**. In cohort-1, usage of the BuMA stent was independently associated with an increased percentage of strut coverage at one month (estimated coefficient: 14.55 [95% CI: 1.72 to 27.38]) while post-dilatation was not. In cohort-2, usage of the BuMA stent was not associated with the percentage of strut coverage. Although

post-dilatation has a tendency towards a negative effect on the percentage of strut coverage at two months, there was no statistical significance (estimated coefficient: -25.15 [95% CI: -50.31 to 0]).

Discussion

POTENTIAL ADVANTAGE OF EARLY STRUT COVERAGE

The PIONEER-II OCT randomised controlled trial results suggest that the BuMA Supreme stent has faster coverage at an earlier time point (one month) than the XIENCE stent, whereas the difference in tissue coverage becomes less evident at a later time point (two months).

In human post-mortem histological analysis, it was reported that strut coverage was an independent predictor for late (>30 days) stent thrombosis²⁴. In the current trial, most of the BuMA stent could reach sufficient strut coverage at the early time point. In cohort-1, the lesions with uncovered struts >30% were 16.7% (3/18) in the BuMA arm as opposed to 41.2% (7/17) in the XIENCE arm one month after implantation (p-value 0.15). In cohort-2, the lesions with uncovered struts >30% were 28.6% (6/21) in the BuMA arm as opposed to 42.3% (11/26) in the XIENCE arm two months after implantation (p-value 0.33).

Table 1. Results of multilevel linear regression analysis examining the correlation between strut coverage and baseline and procedural characteristics.

	Cohort	÷1	Cohort-2			
	Coefficient [95% CI]	<i>p</i> -value	Coefficient [95% CI]	<i>p</i> -value		
Age	-0.31 [-1.04, 0.41]	0.38	0.32 [-0.37, 1.01]	0.36		
Male gender	-4.69 [-18.16, 8.78]	0.48	4.44 [-8.64, 17.52]	0.49		
Diabetes	8.58 [-4.54, 21.69]	0.19	2.80 [-13.13, 18.74]	0.72		
Unstable angina pectoris	-9.26 [-21.75, 3.23]	0.14	-9.57 [-23.33, 4.19]	0.17		
BuMA usage	14.55 [1.72, 27.38]	0.03	11.08 [–2.74, 24.9]	0.11		
Post-dilatation performed	5.28 [-13.67, 24.22]	0.57	-25.15 [-50.31, 0]	0.05		




Figure 5. Representative cases in the PIONEER-II OCT randomised controlled trial. A) OCT in cohort-1. B) OCT in cohort-2. Three-dimensional (3D) reconstructions of OCT pullbacks of the representative cases in the BuMA and the XIENCE arms (upper). Cross-sectional images and 2D fold-out views of strut coverage of the same cases, in which the neointimal thickness on struts is presented with colour coding (lower).

THE ROLE OF FAST DRUG ELUTION AND POLYMER DEGRADATION

Generally, the technological determining factors of strut coverage after DES implantation have been related to the strut (its thickness), to the coating polymer (its biocompatibility or degradation process) and to the drug (release kinetics)²⁵⁻²⁷. The faster coverage

of the BuMA Supreme stent could be mainly due to the faster drug release of sirolimus. The strut thickness of the BuMA stent is comparable to that of the XIENCE stent (80 μ m for the BuMA stent and 81 μ m for the XIENCE stent)^{6,10}. The duration of complete polymer degradation of the BuMA stent is approximately eight weeks, which is shorter than other biodegradable

polymer-coated stents⁶. This technological aspect may play a role during the study phase (two months) but obviously acts longer than two months after implantation. The BuMA stent releases the drug relatively rapidly so that the release of sirolimus from the stent is nearly complete by 28 days after implantation, whereas the drug release of the XIENCE stent continues until 120 days after stenting (**Figure 1**)^{7,28} The drug concentration in the polymer coating of the BuMA stent is 120 μ m/cm² while that of the XIENCE stent is 100 μ m/cm², suggesting that the BuMA stent elutes more concentrated cytostatic drug in a shorter time^{4,10}. Nevertheless, this aspect probably did not impact on the early-phase neointimal healing in the current OCT study.

The fast strut coverage due to the modestly shorter drug elution of the BuMA stent has the potential to reduce the risk of stent thrombosis while lowering the risk of excessive luminal loss. The concept of fast drug elution was previously tested with the Endeavor® zotarolimus-eluting stent (E-ZES) (Medtronic). The E-ZES was designed to release the drug rapidly within 10 days and demonstrated a lower incidence of safety endpoints (stent thrombosis) and comparable efficacy endpoints (target lesion revascularisation) to first-generation DES (sirolimus-eluting stent and paclitaxel-eluting stent)^{29,30}. Regardless of the favourable clinical outcomes of E-ZES, this technology was replaced by the Resolute zotarolimus-eluting stent (R-ZES) with long-lasting drug elution, possibly because of the efficacy concern derived from angiographic results: the E-ZES generated a higher angiographic LLL and binary restenosis rate than other DES, even though the clinical outcomes were comparable²⁹. The BuMA Supreme has the same concept but the drug elution of BuMA is slower than that of E-ZES (28 days vs. 10 days). The difference in drug release between BuMA and R-ZES/E-ZES is in line with the results of angiographic in-stent LLL of each device in first-in-man trials (BuMA: 0.29±0.33 mm at nine months, R-ZES: 0.14±0.37 mm at nine months and E-ZES: 0.61±0.44 mm at 12 months), suggesting that the inhibition capability of neointimal proliferation of the BuMA stent is intermediate between R-ZES and E-ZES^{4,31}. The faster strut coverage could be an advantage of the BuMA stent unless the neointimal hyperplasia is excessive up to the point of creating luminal obstruction in a later phase. Further studies, such as an OCT study with late-phase quantitative and qualitative neointimal assessment after BuMA implantation, are warranted in a larger population to investigate further the benefit of their initial faster strut coverage.

POTENTIAL IMPACT OF BASELINE PROPENSITY ON NEOINTIMAL PROLIFERATION

In the current trial, the percentage of diabetic patients was numerically higher in the XIENCE arm in both cohorts (cohort-1: 18.8% [3/16] for BuMA vs. 46.7% [7/15] for XIENCE, p-value 0.14; cohort-2: 14.3% [3/21] vs. 39.1% [9/23], p-value 0.06) at baseline. Diabetes is known to be associated with greater and faster neointimal hyperplasia compared to non-diabetes 16,17. This is considered to be due to the hormonal and vascular abnormalities which

promote smooth muscle cell proliferation after vascular injury¹⁶. The difference in the percentage of diabetic patients between the two arms might have influenced the results of the OCT study, even though this difference affects strut coverage in favour of the XIENCE arm.

Implantation technique can influence the neointimal response and strut coverage. In the current study, the percentage of lesions with post-dilatation in the BuMA arm was higher than that in the XIENCE arm (cohort-1: 100% [18/18] for BuMA vs. 76.5% [13/18] for XIENCE, p=0.045; cohort-2: 95.2% [20/21] vs. 85.7% [24/28], p=0.38) (Supplementary Table 1). Post-dilatation can be associated with better strut apposition and, furthermore, with greater expansion of the vessel wall and deeper embedment of struts into the vessel wall, but can induce vessel injury followed by increased neointimal hyperplasia^{19,20}. Thus, post-dilatation may result in a good strut coverage in a late phase. More frequent postdilatation in the BuMA arm in the one-month cohort might have had a favourable influence for the BuMA arm in terms of neointimal coverage. In the current study, these factors were not independent predictors of strut coverage either one or two months after implantation, but a potential impact may not be completely excluded.

Furthermore, patients with high bleeding risk were included in the current trial. Mainly, the patients had a history of cerebrovascular disease (ischaemic and haemorrhagic) (Supplementary Table 2). To the best of our knowledge, there are no public domain data with regard to the relationship between bleeding propensity and early-phase neointimal proliferation after stenting.

Limitations

A limitation of the study is its small sample size and non-inferiority design. Strut coverage was assessed by visual assessment by experienced core lab analysts, which may have a limited reproducibility32. The current study did not have OCT assessments at baseline. Vessel healing with neointimal coverage depends on the underlying lesion morphology. Neointimal qualitative assessment was not performed because at one or two months very thin neointima precluded precise assessment of the neointimal quality, even using OCT-based quantitative tissue characterisation software³³. Non-significant trends in baseline characteristics could have confounded the results: for example, 1) diabetes rates were numerically higher in the XIENCE arm, which could have influenced the strut tissue coverage in favour of XIENCE; 2) post-dilatation rates were numerically higher in the BuMA arm, which could have affected the strut tissue coverage. Lastly, the non-serial design does not allow assessment of serial changes of coverage over time.

Conclusions

The results of the PIONEER-II OCT randomised controlled trial suggest that the BuMA Supreme stent has a faster coverage at a very early time point (one month) than the XIENCE stent, presumably due to faster sirolimus elution. The difference in tissue coverage became less evident at a later time point (two months).

Impact on daily practice

Fast strut tissue coverage after implantation has the potential to reduce the risk of stent thrombosis. In the current OCT study, the BuMA Supreme stent demonstrated better strut coverage compared with the XIENCE stent at one month. Taking advantage of the fast strut coverage, an abbreviated duration of dual antiplatelet therapy may be applied to patients treated with the BuMA Supreme stent.

Guest Editor

This paper was guest edited by Antonio Colombo, MD, FACC, FESC, FSCAI; Department of Interventional Cardiology, EMO GVM Centro Cuore Columbus, Milan, Italy.

Acknowledgements

This manuscript was generated mainly by the contribution of the two cooperative teams. The Chinese PLA general hospital team was responsible for the design and execution of the project. Cardialysis was responsible for the core lab analysis and its academic team drafted the manuscript.

Funding

This study was funded by SINOMED, Tianjin, China.

Conflict of interest statement

P.W. Serruys is a member of the Advisory Board for SINOMED. The other authors have no conflicts of interest to declare. The Guest Editor has no conflicts of interest to declare.

References

- 1. Stone GW, Rizvi A, Sudhir K, Newman W, Applegate RJ, Cannon LA, Maddux JT, Cutlip DE, Simonton CA, Sood P, Kereiakes DJ; SPIRIT IV Investigators. Randomized comparison of everolimus- and paclitaxel-eluting stents. 2-year follow-up from the SPIRIT (Clinical Evaluation of the XIENCE V Everolimus Eluting Coronary Stent System) IV trial. *J Am Coll Cardiol*. 2011; 58:19-25.
- 2. Räber L, Magro M, Stefanini GG, Kalesan B, van Domburg RT, Onuma Y, Wenaweser P, Daemen J, Meier B, Jüni P, Serruys PW, Windecker S. Very late coronary stent thrombosis of a newer-generation everolimus-eluting stent compared with early-generation drug-eluting stents: a prospective cohort study. *Circulation*. 2012;125:1110-21.
- 3. Stefanini GG, Byrne RA, Serruys PW, de Waha A, Meier B, Massberg S, Jüni P, Schömig A, Windecker S, Kastrati A. Biodegradable polymer drug-eluting stents reduce the risk of stent thrombosis at 4 years in patients undergoing percutaneous coronary intervention: a pooled analysis of individual patient data from the ISAR-TEST 3, ISAR-TEST 4, and LEADERS randomized trials. *Eur Heart J.* 2012;33:1214-22.
- 4. von Birgelen C, Asano T, Amoroso G, Aminian A, Brugaletta S, Vrolix M, Hernandez-Antolin R, van de Harst P,

- Iniguez A, Janssens L, Smits PC, Wykrzykowska JJ, Ribeiro VG, Pereira H, Canas da Silva P, Piek JJ, Onuma Y, Serruys PW, Sabaté M. First-in-man randomized comparison of BuMA Supreme biodegradable polymer sirolimus-eluting versus durable polymer zotarolimus-eluting coronary stents: the PIONEER trial. *EuroIntervention*. 2018;13:2026-35.
- 5. Asano T, Katagiri Y, Collet C, Tenekecioglu E, Miyazaki Y, Sotomi Y, Amoroso G, Aminian A, Brugaletta S, Vrolix M, Hernandez-Antolin R, van de Harst P, Iniguez A, Janssens L, Smits P, Wykrzykowska JJ, Ribeiro VG, Pereira H, da Silva PC, Piek JJ, Reiber JHC, von Birgelen C, Sabaté M, Onuma Y, Serruys PW. Functional comparison between BuMA Supreme biodegradable polymer sirolimus-eluting and durable polymer zotarolimus-eluting coronary stents using quantitative flow ratio: PIONEER QFR substudy. *EuroIntervention*. 2018;14:e570-9.
- 6. Qian J, Zhang YJ, Xu B, Yang YJ, Yan HB, Sun ZW, Zhao YL, Tang YD, Gao Z, Chen J, Cui JG, Mintz GS, Gao RL. Optical coherence tomography assessment of a PLGA-polymer with electro-grafting base layer versus a PLA-polymer sirolimus-eluting stent at three-month follow-up: the BuMA-OCT randomised trial. *EuroIntervention*. 2014;10:806-14.
- 7. Xu B, Gao R, Yang Y, Cao X, Qin L, Li Y, Li Z, Li X, Lin H, Guo Y, Ma Y, Wang J, Nie S, Xu L, Cao E, Guan C, Stone GW; PANDA III Investigators. Biodegradable Polymer-Based Sirolimus-Eluting Stents With Differing Elution and Absorption Kinetics: The PANDA III Trial. *J Am Coll Cardiol*. 2016;67:2249-58.
- 8. El-Hayek G, Bangalore S, Casso Dominguez A, Devireddy C, Jaber W, Kumar G, Mavromatis K, Tamis-Holland J, Samady H. Meta-Analysis of Randomized Clinical Trials Comparing Biodegradable Polymer Drug-Eluting Stent to Second-Generation Durable Polymer Drug-Eluting Stents. *JACC Cardiovasc Interv.* 2017;10:462-73.
- 9. Cassese S, Xhepa E, Ndrepepa G, Kufner S, Colleran R, Giacoppo D, Koppara T, Mankerious N, Byrne RA, Laugwitz KL, Schunkert H, Fusaro M, Kastrati A, Joner M. Vascular response to percutaneous coronary intervention with biodegradable-polymer vs. new-generation durable-polymer drug-eluting stents: a meta-analysis of optical coherence tomography imaging trials. *Eur Heart J Cardiovasc Imaging*. 2018;19:1294-301.
- 10. Ding NI, Pacetti SD, Tang FW, Gada M, Roorda W. XIENCE VTM Stent Design and Rationale. *Journal of Interventional Cardiology.* 2009;22:S18-27.
- 11. Authors/Task Force members, Windecker S, Kolh P, Alfonso F, Collet JP, Cremer J, Falk V, Filippatos G, Hamm C, Head SJ, Jüni P, Kappetein AP, Kastrati A, Knuuti J, Landmesser U, Laufer G, Neumann FJ, Richter DJ, Schauerte P, Sousa Uva M, Stefanini GG, Taggart DP, Torracca L, Valgimigli M, Wijns W, Witkowski A. 2014 ESC/EACTS Guidelines on myocardial revascularization: The Task Force on Myocardial Revascularization of the European Society of Cardiology (ESC) and the European Association for Cardio-Thoracic Surgery (EACTS)Developed with the special contribution of the European Association of Percutaneous Cardiovascular Interventions (EAPCI). *Eur Heart J.* 2014;35:2541-619.

- 12. Serruys PW, de Jaegere P, Kiemeneij F, Macaya C, Rutsch W, Heyndrickx G, Emanuelsson H, Marco J, Legrand V, Materne P, et al. A comparison of balloon-expandable-stent implantation with balloon angioplasty in patients with coronary artery disease. Benestent Study Group. *N Engl J Med.* 1994;331:489-95.
- 13. Suwannasom P, Onuma Y, Benit E, Gach O, von Birgelen C, Hofma SH, Sotomi Y, Bo X, Zhang YJ, Gao R, Garcia-Garcia HM, Wykrzykowska JJ, de Winter RJ, Serruys PW. Evaluation of vascular healing of polymer-free sirolimus-eluting stents in native coronary artery stenosis: a serial follow-up at three and six months with optical coherence tomography imaging. *EuroIntervention*. 2016;12: e574-83.
- 14. Räber L, Onuma Y, Brugaletta S, García-García H, Backx B, Iñiguez A, Okkels-Jensen L, Cequier A, Christiansen E, Suttorp M, Serruys P, Sabaté M, Windecker S. Arterial healing following primary PCI using the Absorb everolimus-eluting bioresorbable vascular scaffold (Absorb BVS) versus the durable polymer everolimus-eluting metallic stent (XIENCE) in patients with acute ST-elevation myocardial infarction: rationale and design of the randomized TROFI II study. *EuroIntervention*. 2016;12:482-9.
- 15. Sabaté M, Windecker S, Iniguez A, Okkels-Jensen L, Cequier A, Brugaletta S, Hofma SH, Räber L, Christiansen EH, Suttorp M, Pilgrim T, Anne van Es G, Sotomi Y, Garcia-Garcia HM, Onuma Y, Serruys PW. Everolimus-eluting bioresorbable stent vs. durable polymer everolimus-eluting metallic stent in patients with ST-segment elevation myocardial infarction: results of the randomized ABSORB ST-segment elevation myocardial infarction-TROFI II trial. *Eur Heart J.* 2016;37:229-40.
- 16. Aronson D, Bloomgarden Z, Rayfield EJ. Potential mechanisms promoting restenosis in diabetic patients. *J Am Coll Cardiol*. 1996;27:528-35.
- 17. Tanaka N, Terashima M, Rathore S, Itoh T, Habara M, Nasu K, Kimura M, Itoh T, Kinoshita Y, Ehara M, Tsuchikane E, Asakura K, Asakura Y, Katoh O, Suzuki T. Different patterns of vascular response between patients with or without diabetes mellitus after drug-eluting stent implantation: optical coherence tomographic analysis. *JACC Cardiovasc Interv.* 2010;3:1074-9.
- 18. Kim JS, Fan C, Choi D, Jang IK, Lee JM, Kim TH, Park SM, Paik SI, Ko YG, Hong MK, Jang Y, Chung N. Different patterns of neointimal coverage between acute coronary syndrome and stable angina after various types of drug-eluting stents implantation; 9-month follow-up optical coherence tomography study. *Int J Cardiol.* 2011;146:341-6.
- 19. van Beusekom H, Sorop O, Weymaere M, Duncker D, van der Giessen W. The neointimal response to stents eluting tacrolimus from a degradable coating depends on the balance between polymer degradation and drug release. *EuroIntervention*. 2008;4:139-47.
- 20. Suwannasom P, Sotomi Y, Corti R, Kurz DJ, Roffi M, von Birgelen C, Buzzi S, Zucker A, Dijkstra J, Wykrzykowska JJ, de Winter RJ, Windecker S, Onuma Y, Serruys PW, Daemen J, Räber L. First-in-man six-month results of a surface-modified coronary stent system in native coronary stenosis. *EuroIntervention*. 2017;12:2118-27.

- 21. Takahara M, Kitahara H, Nishi T, Miura K, Miyayama T, Sugimoto K, Nakayama T, Fujimoto Y, Kobayashi Y. Very early tissue coverage after drug-eluting stent implantation: an optical coherence tomography study. *Int J Cardiovasc Imaging*. 2017; 33:25-30.
- 22. Yano H, Horinaka S, Ishikawa M, Ishimitsu T. Early vascular responses after everolimus-eluting stent implantation assessed by serial observations of intracoronary optical coherence tomography. *Heart Vessels.* 2017;32:804-12.
- 23. Kim JS, Kim JH, Shin DH, Kim BK, Ko YG, Choi D, Jang Y, Hong MK. Serial randomized comparison of strut coverage of everolimus- and first-generation sirolimus-eluting stents. *Can J Cardiol.* 2015;31:723-30.
- 24. Finn AV, Joner M, Nakazawa G, Kolodgie F, Newell J, John MC, Gold HK, Virmani R. Pathological correlates of late drug-eluting stent thrombosis: strut coverage as a marker of endothelialization. *Circulation*. 2007;115:2435-41.
- 25. Bourantas CV, Räber L, Zaugg S, Sakellarios A, Taniwaki M, Heg D, Moschovitis A, Radu M, Papafaklis MI, Kalatzis F, Naka KK, Fotiadis DI, Michalis LK, Serruys PW, Garcia Garcia HM, Windecker S. Impact of local endothelial shear stress on neointima and plaque following stent implantation in patients with ST-elevation myocardial infarction: A subgroup-analysis of the COMFORTABLE AMI-IBIS 4 trial. *Int J Cardiol*. 2015;186: 178-85.
- 26. Virmani R, Guagliumi G, Farb A, Musumeci G, Grieco N, Motta T, Mihalcsik L, Tespili M, Valsecchi O, Kolodgie FD. Localized hypersensitivity and late coronary thrombosis secondary to a sirolimus-eluting stent: should we be cautious? *Circulation*. 2004;109:701-5.
- 27. Guagliumi G, Ikejima H, Sirbu V, Bezerra H, Musumeci G, Lortkipanidze N, Fiocca L, Tahara S, Vassileva A, Matiashvili A, Valsecchi O, Costa M. Impact of drug release kinetics on vascular response to different zotarolimus-eluting stents implanted in patients with long coronary stenoses: the LongOCT study (Optical Coherence Tomography in Long Lesions). *JACC Cardiovasc Interv.* 2011;4:778-85.
- 28. Serruys PW, Ruygrok P, Neuzner J, Piek JJ, Seth A, Schofer JJ, Richardt G, Wiemer M, Carrié D, Thuesen L, Boone E, Miquel-Herbert K, Daemen J. A randomised comparison of an everolimuseluting coronary stent with a paclitaxel-eluting coronary stent:the SPIRIT II trial. *EuroIntervention*. 2006;2:286-94.
- 29. Kandzari DE, Leon MB, Popma JJ, Fitzgerald PJ, O'Shaughnessy C, Ball MW, Turco M, Applegate RJ, Gurbel PA, Midei MG, Badre SS, Mauri L, Thompson KP, LeNarz LA, Kuntz RE; ENDEAVOR III Investigators. Comparison of zotarolimus-eluting and sirolimus-eluting stents in patients with native coronary artery disease: a randomized controlled trial. *J Am Coll Cardiol.* 2006;48:2440-7.
- 30. Kirtane AJ, Leon MB, Ball MW, Bajwa HS, Sketch MH Jr, Coleman PS, Stoler RC, Papadakos S, Cutlip DE, Mauri L, Kandzari DE; ENDEAVOR IV Investigators. The "final" 5-year follow-up from the ENDEAVOR IV trial comparing a zotarolimus-

eluting stent with a paclitaxel-eluting stent. *JACC Cardiovasc Interv.* 2013;6:325-33.

- 31. Meredith IT, Ormiston J, Whitbourn R, Kay IP, Muller D, Bonan R, Popma JJ, Cutlip DE, Fitzgerald P, Prpic R, Kuntz RE. First-in-human study of the Endeavor ABT-578-eluting phosphorylcholine-encapsulated stent system in de novo native coronary artery lesions: Endeavor I Trial. *EuroIntervention*. 2005;1:157-64.
- 32. Brugaletta S, Garcia-Garcia HM, Gomez-Lara J, Radu MD, Pawar R, Khachabi J, Bruining N, Sabaté M, Serruys PW. Reproducibility of qualitative assessment of stent struts coverage by optical coherence tomography. *Int J Cardiovasc Imaging*. 2013; 29:5-11.
- 33. Liu S, Sotomi Y, Eggermont J, Nakazawa G, Torii S, Ijichi T, Onuma Y, Serruys PW, Lelieveldt BPF, Dijkstra J. Tissue characterization with depth-resolved attenuation coefficient and backscatter term in intravascular optical coherence tomography images. *J Biomed Opt.* 2017;22:1-16.
- 34. Subherwal S, Bach RG, Chen AY, Gage BF, Rao SV, Newby LK, Wang TY, Gibler WB, Ohman EM, Roe MT, Pollack CV Jr, Peterson ED, Alexander KP. Baseline risk of major bleeding in non-ST-segment-elevation myocardial infarction: the

CRUSADE (Can Rapid risk stratification of Unstable angina patients Suppress ADverse outcomes with Early implementation of the ACC/AHA Guidelines) Bleeding Score. *Circulation*. 2009;119:1873-82.

35. Mehran R, Pocock SJ, Nikolsky E, Clayton T, Dangas GD, Kirtane AJ, Parise H, Fahy M, Manoukian SV, Feit F, Ohman ME, Witzenbichler B, Guagliumi G, Lansky AJ, Stone GW. A risk score to predict bleeding in patients with acute coronary syndromes. *J Am Coll Cardiol*. 2010;55:2556-66.

Supplementary data

Supplementary Appendix 1. Study eligibility criteria.

Supplementary Table 1. Baseline characteristics.

Supplementary Table 2. Bleeding propensity of the enrolled patients.

Supplementary Table 3. Results of quantitative coronary angiography.

Supplementary Table 4. Results of optical coherence tomography.

The supplementary data are published online at: http://www.pcronline.com/eurointervention/145th issue/236



Supplementary data

Supplementary Appendix 1. Study eligibility criteria

Inclusion criteria:

- 1) The patient is at least 18 to 85 years of age;
- 2) Evidence of myocardial ischaemia without raised troponin (e.g., stable or unstable angina, silent ischaemia demonstrated by positive territorial functional study);
- 3) CAD patients with a high risk of bleeding. Any one or more situation listed below can be considered a patient with high bleeding risk by the investigator:
 - a) Adjunctive oral anticoagulation treatment planned to continue after PCI
 - b) Baseline Hb ≤11 g/dl (or anaemia requiring transfusion during the prior 4 weeks)
 - c) Any prior intracerebral bleed
 - d) Any prior ischaemic stroke
 - e) Hospital admission for bleeding during the prior 12 months
 - f) Non-skin cancer diagnosed or treated ≤3 years
 - g) Planned daily NSAID or steroids for ≥30 days after PCI
 - h) Planned major surgery (within 1 year)
 - i) Renal failure (calculated creatinine clearance ≤40 ml/min) including the status with haemodialysis
 - j) Thrombocytopaenia (platelet count ≤100,000/mm³ or below normal level)
 - k) Severe chronic liver disease (variceal haemorrhage, ascites, hepatic encephalopathy or jaundice)
 - 1) Previous history of gastrointestinal bleeding
 - m) Other bleeding risk factors Note (1)
- 4) The patient has a planned intervention of up to two de novo lesions, in different epicardial vessels;
- 5) Lesion(s) must have a visually estimated diameter stenosis of \geq 70% and \leq 100%;
- Reference vessel diameter (RVD) must be visually estimated \geq 2.5 and \leq 4.0 mm, and the vessel length must be no more than 40 mm;
- 7) Written informed consent;
- 8) No coronary artery bypass graft surgery contraindications;
- 9) The patient and the patient's physician agree to the follow-up visits including angiographic follow-up and OCT controls at 1 or 2 months.

Exclusion criteria:

- 1) Evidence of ongoing acute myocardial infarction in ECG prior to procedure;
- 2) Left ventricular ejection fraction (LVEF) <30%;
- 3) The patient is a recipient of a heart transplant;
- 4) Known hypersensitivity or contraindication to aspirin, both heparin and bivalirudin, antiplatelet medication specified for use in the study (clopidogrel and ticlopidine), sirolimus, everolimus or cobalt metal. Sensitivity to contrast media, which cannot be adequately premedicated;
- 5) Other medical illness (e.g., skin cancer diagnosed or treated >3 years, neurological deficiency) or known history of substance abuse (alcohol, etc.) excluding those with relative risk as per investigator judgement;
- 6) Per investigator judgement that may cause non-compliance with the protocol or confound the data interpretation or is associated with a limited life expectancy (less than 1 year);

- 7) Pregnant or breastfeeding woman or woman in fertile period not taking adequate contraceptives;
- 8) Chronic total occlusion (TIMI 0), left main lesion, intervention-required three-vessel lesions, branch vessel diameter \geq 2.5 mm and bypass lesion;
- 9) Patients expected not to comply with 1-month DAPT;
- 10) Compliance with long-term single antiplatelet therapy unlikely;
- 11) Active bleeding at the time of inclusion;
- 12) Patients requiring a planned staged PCI procedure more than one week after the index procedure;
- 13) Procedure planned to require non-study stents, or stand-alone POBA or stand-alone atherectomy;
- 14) Reference vessel diameter <2.25 or >4.0 mm, vessel length >35 mm;
- 15) Cardiogenic shock;
- 16) Participation in another clinical trial (12 months after index procedure);
- 17) Those who are not suitable to attend this trial after the evaluation by the doctor.

OCT exclusion criteria:

- 1) Severe tortuous, calcified or angulated coronary anatomy of the study vessel that in the opinion of the investigator would result in suboptimal imaging or excessive risk of complication from placement of an OCT catheter;
- 2) Total occlusion or Thrombolysis In Myocardial Infarction, TIMI 0, prior to wire crossing.

Note (1):

Patient meets the following 2 or more factors from "aged female", "diabetes", "uncontrolled high blood pressure", "heart failure" and other potential bleeding risk factors^{34,35}.

Supplementary Table 1. Baseline characteristics.

	Coho	rt-1	Coho	Cohort-2		
Baseline demographics	BuMA Supreme	XIENCE	BuMA Supreme	XIENCE		
	N=16	N=15	N=21	N=23		
Male gender	62.5%	60.0%	61.9%	56.5%		
Age (years)	63.4±10.0	63.5±8.2	64.1±9.4	65.0±10.3		
Diabetes	18.8%	46.7%	14.3%	39.1%		
Dyslipidaemia	25.0%	13.3%	42.9%	21.7%		
Hypertension	68.8%	66.7%	76.2%	78.3%		
Renal failure	0.0%	0.0%	0.0%	0.0%		
Atrial fibrillation	0.0%	0.0%	0.0%	13.0%		
History of bleeding	12.5%	20.0%	9.5%	4.4%		
Previous cerebral vessel diseases	18.8%	20.0%	19.1%	26.1%		
Previous myocardial infarction	6.3%	6.7%	0.0%	0.0%		
History of PCI	0.0%	6.7%	0.0%	0.0%		
History of CABG	0.0%	0.0%	0.0%	0.0%		
Clinical presentation						
Silent ischaemia	12.5%	13.3%	10.0%	4.4%		
Stable angina pectoris	18.8%	13.3%	15.0%	26.1%		
Unstable angina pectoris	62.5%	66.7%	75.0%	65.2%		
Unknown	6.3%	6.7%	0.0%	4.4%		
Antiplatelet/anticoagulant drugs						
Aspirin	93.8%	100.0%	95.2%	100.0%		
Clopidogrel	87.5%	86.7%	81.0%	87.0%		
Ticagrelor	31.3%	40.0%	33.3%	21.7%		
Oral anticoagulant (at discharge)	0.0%	0.0%	0.0%	4.4%		
g.,						
Lesion characteristics	BuMA Supreme	XIENCE	BuMA Supreme	XIENCE		
	L=18	L=17	L=21	L=28		
Target vessel						
LM	0%	0%	0%	0%		
LAD	55.56%	35.29%	42.86%	50.00%		
LCX	11.11%	11.76%	23.81%	10.71%		
RCA	33.33%	52.94%	33.33%	39.29%		
Preprocedural QCA (core lab)						
Lesion length (mm)	14.23±7.60	15.20±9.81	15.88±6.11	16.60 ± 9.45		
Reference lumen diameter (mm)	2.83±0.49	2.84 ± 0.48	2.95±0.43	2.80 ± 0.54		
Minimum lumen diameter (mm)	1.14±0.42	1.06 ± 0.42	1.28±0.41	1.06 ± 0.41		
Diameter stenosis (%)	60.0±12.0	61.5±16.3	56.3±13.3	62.1±12.8		
Pre-balloon dilatation performed (%)	77.8%	82.4%	61.9%	85.7%		
Balloon nominal diameter (mm)	2.25±0.26	2.21±0.24	2.38±0.30	2.36±0.28		
Maximum pressure (atm)	10.6±2.95	11.9±2.77	11.4±1.89	11.8±3.3		
Stent implantation						
Stent nominal diameter (mm)	3.20±0.34	3.04±0.42	3.21±0.42	3.13±0.42		
Stent length (mm)	25.5±6.9	22.4±6.1	25.52±8.17	25.34±7.38		
Maximum dilation pressure (atm)	12.7±2.7	12.0±2.6	12.59±2.93	12.00±2.77		
Post-dilation performed	100.0%	76.5%	95.2%	85.7%		
Balloon nominal diameter (mm)	3.51±0.46	3.33±0.54	3.51±0.43	3.54±0.49		
Maximum dilation pressure (atm)	18.4±3.9	17.3±3.0	17.52±3.40	18.57±3.30		
Non-compliant balloon usage	83.3%	84.6%	85.0%	83.3%		
Device success (%)	100%	100%	100%	100%		

Non-compliant balloon usage 83.3% 84.6% 85.0% 83.3%

Device success (%) 100% 100% 100% 100%

CABG: coronary artery bypass graft; L: number of lesions; LAD: left anterior descending; LCX: left circumflex; N: number of patients; PCI: percutaneous coronary intervention; RCA: right coronary artery

Supplementary Table 2. Bleeding propensity of the enrolled patients.

Haemorrhagic disease

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	2	3
Cohort-2 (BuMA N=21, XIENCE N=23)	2	1

Cerebrovascular disease -- ischaemic

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	3	3
Cohort-2 (BuMA N=21, XIENCE N=23)	3	5

Cerebral vascular disease -- haemorrhagic

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	0	0
Cohort-2 (BuMA N=21, XIENCE N=23)	1	1

Non-skin cancer

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	0	0
Cohort-2 (BuMA N=21, XIENCE N=23)	1	2

Renal failure

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	0	0
Cohort-2 (BuMA N=21, XIENCE N=23)	0	0

Anaemia

	BuMA Supreme	XIENCE
Cohort-1 (BuMA N=16, XIENCE N=15)	0	0
Cohort-2 (BuMA N=21, XIENCE N=23)	1	1

Supplementary Table 3. Results of quantitative coronary angiography.

		Cohort-1 (1-month follow-up)						Cohort-2 (2-month follow-up)				
		MA Supreme =16 patients / 18 lesions)	(N	XIENCE =15 patients / 17 lesions)	Difference [95% CI]	<i>p</i> -value	(N	MA Supreme =21 patients / 21 lesions)	(N	XIENCE =23 patients / 28 lesions)	Difference [95% CI]	<i>p</i> -value
In-stent analysis						•			1			•
Pre-procedure												
Reference vessel diameter (mm)	18	2.83±0.49	16	2.84±0.48	-0.01 [-0.34, 0.33]	0.96	20	2.95±0.43	26	2.80±0.54	0.15 [-0.15, 0.44]	0.33
Minimum lumen diameter (mm)	18	1.14±0.42	17	1.06±0.42	0.08 [-0.21, 0.36]	0.59	20	1.28±0.41	27	1.06±0.41	0.22 [-0.02, 0.46]	0.07
Percentage diameter stenosis	18	60.0±12.0	17	61.5±16.3	-1.5 [-11.3, 8.3]	0.75	20	56.3±13.3	27	62.1±12.8	-5.8 [-13.5, 1.9]	0.14
Obstruction length (mm)	18	14.23±7.60	16	15.20±9.81	-0.97 [-7.06, 5.13]	0.75	20	15.88±6.11	26	16.60±9.45	-0.72 [-5.62, 4.18]	0.77
Post-procedure												
Reference vessel diameter (mm)	18	3.04±0.41	17	2.83±0.46	0.21 [-0.09, 0.51]	0.16	21	2.92±0.44	27	3.00±0.49	-0.08 [-0.35, 0.20]	0.58
Minimum lumen diameter (mm)	18	2.72±0.39	17	2.51±0.42	0.21 [-0.07, 0.49]	0.13	21	2.59±0.42	27	2.64±0.47	-0.05 [-0.31, 0.21]	0.71
Percentage diameter stenosis	18	10.4±5.5	17	11.1±5.6	-0.7 [-4.5, 3.1]	0.71	21	11.2±5.9	27	11.8±7.8	-0.6 [-4.7, 3.5]	0.77
Acute lumen gain (mm)	15	1.63±0.47	15	1.47±0.63	0.17 [-0.25, 0.59]	0.41	17	1.32±0.53	25	1.52±0.32	-0.20 [-0.47, 0.06]	0.17
Follow-up												
Reference vessel diameter (mm)	17	2.93±0.35	16	2.75±0.50	0.18 [-0.12, 0.48]	0.24	21	2.87±0.39	28	2.90±0.45	-0.03 [-0.28, 0.22]	0.81
Minimum lumen diameter (mm)	17	2.61±0.36	16	2.50±0.48	0.11 [-0.18, 0.41]	0.44	21	2.53±0.42	28	2.57±0.43	-0.04 [-0.29, 0.21]	0.75
Percentage diameter stenosis	17	10.9±6.6	16	9.2±5.7	1.7 [-2.7, 6.1]	0.44	21	12.1±6.2	28	11.5±5.2	0.5 [-2.8, 3.8]	0.74
Binary restenosis	17	0 (0%)	16	0 (0%)		NA	21	0 (0%)	28	0 (0%)		NA
Late lumen loss (mm)	17	0.12±0.25	16	0.03±0.20	0.09 [-0.07, 0.25]	0.25	21	0.06±0.21	27	0.06±0.22	0.00 [-0.12, 0.13]	0.96
In-segment analysis												
Post-procedure												
Reference vessel diameter (mm)	18	2.94±0.43	17	2.78±0.49	0.15 [-0.16, 0.47]	0.33	21	2.77±0.45	27	2.93±0.53	-0.17 [-0.46, 0.13]	0.26
Minimum lumen diameter (mm)	18	2.37±0.54	17	2.29±0.51	0.09 [-0.27, 0.44]	0.63	21	2.30±0.48	27	2.40±0.53	-0.10 [-0.39, 0.20]	0.51
Percentage diameter stenosis	18	19.8±10.2	17	18.0±9.4	1.8 [-5.0, 8.5]	0.59	21	17.0±7.5	27	18.3±8.8	-1.2 [-6.1, 3.6]	0.61
Follow-up												
Reference vessel diameter (mm)	17	2.88±0.32	16	2.63±0.52	0.25 [-0.06, 0.55]	0.11	21	2.76±0.43	28	2.81±0.49	-0.05 [-0.32, 0.22]	0.71
Minimum lumen diameter (mm)	17	2.35±0.38	16	2.23±0.61	0.12 [-0.23, 0.48]	0.48	21	2.30±0.46	28	2.38±0.48	-0.07 [-0.35, 0.20]	0.59
Percentage diameter stenosis	17	18.4±8.1	16	16.5±10.9	1.9 [-4.9, 8.7]	0.57	21	17.0±7.1	28	15.6±7.3	1.4 [-2.8, 5.6]	0.51
Binary restenosis	17	0 (0%)	16	0 (0%)		NA	21	0 (0%)	28	0 (0%)		NA
Late lumen loss (mm)	17	0.07±0.28	16	0.07±0.29	0.00 [-0.20, 0.21]	0.97	21	0.00±0.24	27	0.02±0.24	-0.02 [-0.16, 0.13]	0.81

CI: confidence interval; L: number of lesions; N: number of patients

Supplementary Table 4. Results of optical coherence tomography.

	Cohort-1 (1-month follow-up)				Cohort-2 (2-month follow-up)			
	BuMA Supreme (N=16 patients / 18 lesions)	XIENCE (N=15 patients / 17 lesions)	Difference [95% CI]	<i>p</i> -value	BuMA Supreme (N=21 patients / 21 lesions)	XIENCE (N=23 patients / 28 lesions*)	Difference [95% CI]	<i>p</i> -value
Follow-up	<u> </u>							
Strut neointimal coverage (%)	83.8±10.4	73.0±17.5	10.8 [1.0, 20.6]	0.037	80.3±18.3	73.3±21.3	7.0 [-4.8, 18.8]	0.24
Mean neointimal hyperplasia area (mm²)	0.63±0.37	0.36±0.12	0.27 [0.08, 0.46]	0.008	0.63±0.40	0.49±0.23	0.14 [-0.05, 0.32]	0.16
Percentage of neointimal hyperplasia volume obstruction	7.5±4.9	4.9±1.7	2.6 [0.1, 5.2]	0.045	7.5±4.8	5.7±2.4	1.8 [-0.4, 3.9]	0.14
Mean malapposition area (mm²)	0.02±0.04	0.05±0.08	-0.03 [-0.07, 0.01]	0.12	0.03±0.07	0.02±0.03	0.01 [-0.02, 0.04]	0.65
Healing score	32.4±21.59	54.9±34.65	-22.5 [-42.3, -2.8]	0.027	39.9±37.8	53.8±42.8	-13.9 [-37.9, 10.1]	0.25
Length stented region (mm)	29.4±12.0	24.8±13.2	4.6 [-4.0, 13.3]	0.29	33.9±15.4	28.0±9.4	6.0 [-1.4, 13.3]	0.13
Mean stent area (mm²)	8.56±1.95	7.52±2.06	1.04 [-0.33, 2.42]	0.13	8.63±2.54	8.84±2.52	-0.22 [-1.71, 1.28]	0.77
Minimum stent area (mm²)	7.08±1.38	6.03±2.24	1.05 [-0.22, 2.32]	0.10	7.13±2.64	7.16±2.22	-0.03 [-1.46, 1.40]	0.96
Mean lumen area (mm²)	8.11±2.04	7.35±2.08	0.76 [-0.66, 2.18]	0.28	8.13±2.47	8.55±2.43	-0.41 [-1.86, 1.03]	0.57
Minimum lumen area (mm²)	6.43±1.39	5.80±2.12	0.63 [-0.59, 1.86]	0.30	6.47±2.51	6.78±2.26	-0.31 [-1.72, 1.09]	0.65
Strut-level analysis								
Number of struts	349.1±149.4	269.7±157.5	79.3 [-26.2, 184.9]	0.14	407.3±198.3	293.3±123.0	114.0 [18.9, 209.0]	0.028
Percentage of uncovered struts	16.2±10.4	27.0±17.5	-10.8 [-20.6, -1.0]	0.037	19.7±18.3	26.7±21.3	-7.0 [-18.8, 4.8]	0.24
Percentage of apposed and covered struts	82.9±11.1	70.9±16.6	11.9 [2.3, 21.6]	0.017	79.0±19.5	72.2±21.4	6.8 [-5.4, 18.9]	0.27
Percentage of malapposed struts	0.7±1.2	1.5±1.8	-0.9 [-1.9, 0.2]	0.10	0.9±2.2	0.8±1.0	0.1 [-0.9, 1.1]	0.82
Percentage of malapposed and covered struts	0.6±1.0	1.4±1.9	-0.8 [-1.9, 0.2]	0.12	0.8±2.1	0.7±1.1	0.1 [-0.8, 1.1]	0.81
Percentage of malapposed and uncovered struts	0.1±0.2	0.1±0.4	-0.0 [-0.2, 0.2]	0.79	0.0±0.1	0.0±0.1	-0.0 [-0.1, 0.1]	0.91

^{*}Two lesions/patients were non-analysable. Therefore, data are based on 26 lesions. CI: confidence interval; L: number of lesions; N: number of patients