

Crossing chronic total occlusions with the Ocelot system: the initial European experience

Arne Schwindt^{1*}, MD; Bernhard Reimers², MD; Dierk Scheinert³, MD; Matthew Selmon⁴, MD; John P. Pigott⁵, MD; Jon C. George⁶, MD; Greg Robertson⁷, MD; Sean Janzer⁸, MD; Huey B. McDaniel⁹, MD; Gataum V. Shrikhande¹⁰, MD; Giovanni Torsello¹, MD; Johannes Schaefer¹, MD; Salvatore Saccà², MD; Francesco Versaci², MD

1. St. Franziskus-Hospital, Münster, Germany; 2. Mirano Hospital, Mirano, Italy; 3. Heart Center-University of Leipzig, Leipzig, Germany; 4. Heart Hospital of Austin, Austin, TX, USA; 5. Jobst Vascular Institute, Toledo, OH, USA; 6. Deborah Heart & Lung Center, Browns Mills, NJ, USA; 7. Emory Heart and Vascular Center, Atlanta, GA, USA; 8. Bryn Mawr Hospital, Bryn Mawr, PA, USA; 9. St. Dominic Hospital, Jackson, MS, USA; 10. NYP/Columbia U. Medical Center, New York, NY, USA

This paper also includes accompanying supplementary data published online at: http://www.pcronline.com/eurointervention/66th_issue/139

KEYWORDS

- chronic total occlusion
- CTO navigation
- optical coherence tomography
- peripheral arterial disease
- popliteal artery
- superficial femoral artery
- true lumen recanalisation

Abstract

Aims: The aim of the study was to determine the safety, efficacy and feasibility of a new chronic total occlusion (CTO) device using optical coherence tomography (OCT) technology, the Ocelot catheter (Avinger, Inc., Redwood City, CA, USA), for crossing of SFA CTOs following guidewire failure.

Methods and results: Prospective, multicentre, market preference testing. Thirty-three patients with confirmed CTO (99-100% stenosis by visual estimate) of their superficial femoral artery (SFA) were treated between September 28, 2011, and December 9, 2011, at three European centres. Ocelot crossed 94% (31/33) of CTOs, allowing guidewire placement in the distal true lumen. All (100%) lesions were treated without any major adverse safety events. Procedural time and contrast dose were significantly reduced ($p < 0.0001$) when compared with a similar, non-OCT-guided CTO crossing device (Wildcat catheter; Avinger, Inc.). Overall physician feedback on the catheter performance was positive with an 87% average rating of excellent or good across seven categories. Performance ratings of Ocelot's OCT imaging guidance were consistently positive with an 86% average rating of excellent or good across five OCT categories.

Conclusions: The Ocelot catheter combines advanced CTO crossing technology with real-time OCT guidance. When compared with a similar non-OCT-guided catheter, crossing efficacy and safety profile improved. Total procedure time and contrast volumes were significantly reduced. The Ocelot is a safe, efficient and effective tool for crossing CTOs.

* Corresponding author: St. Franziskus-Hospital, Hohenzollernring 72, 48145 Münster, Germany.
E-mail: arne.schwindt@sfh-muenster.de

Introduction

It is estimated that peripheral artery disease (PAD) affects approximately 12-14% of the general population, and approaches 20% in patients 75 years of age or older¹. Over 200,000 amputations are reported each year for patients with end-stage PAD in the United States and Europe². As the population continues to live longer, the prevalence of this disease is only expected to grow. Without early diagnosis and risk factor modification, approximately 1-2% of these patients will develop critical limb ischaemia (CLI)³, defined as the presence of lower extremity rest pain, with or without gangrene or non-healing ulcers, caused by inadequate blood supply⁴. Without successful revascularisation, CLI patients have significantly increased morbidity with a 95% amputation rate at one year⁵.

The Inter-Society consensus (TASC) II guidelines traditionally recommend surgical intervention for severe PAD including peripheral arterial chronic total occlusions (CTO), denoted as type C and D lesions. Bypass surgery has traditionally been the “gold standard”, with a five-year limb salvage rate of nearly 80%⁶⁻⁹. Recent advances in endovascular therapy enable patients with CLI to be treated using minimally invasive techniques. With revascularisation reported as high as 97%, endovascular crossing options are becoming first-line crossing treatments for patients across a wide health spectrum. Importantly, percutaneous therapeutics are a safe alternative for non-viable surgical candidates as a result of their comorbidities, lack of suitable target vessels or poor venous conduits¹⁰.

Percutaneous crossing of infra-inguinal arterial CTOs is dependent on emerging guidewire and catheter technology. Guidewire therapy rests on the premise of intentional extraluminal (subintimal) crossing, as originally described by Bolia et al¹¹. A recent large systematic review highlights one-year patency following subintimal angioplasty at 50%, with a limb salvage rate of 80-90%¹². Alternatively, targeting intraluminal crossing may alleviate subintimal restenotic factors while enabling a wider spectrum of therapeutic interventions, including stenting and atherectomy. Adam et al¹³ demonstrated statistically equivalent one-year clinical success rates (56% vs. 50%) between surgery and angioplasty. By contrast, the surgical arm was associated with higher complications and costs.

In recent years, medical device companies have continued to develop and improve catheters for crossing CTOs in peripheral arteries. These include the Crosser™ (Bard Peripheral Vascular, Inc., Tempe, AZ, USA), TruePath™ (Boston Scientific Corp., Natick, MA, USA), the Frontrunner® XP (Cordis Corp., Bridgewater, NJ, USA), the CrossBoss™ (BridgePoint, Plymouth, MN, USA) and the Wildcat catheter (Avinger, Inc., Redwood City, CA, USA).

While advancements in percutaneous CTO treatments are ongoing, a true paradigm shift is emerging with the advent of advanced intravascular imaging, including intravascular ultrasound (IVUS) and optical coherence tomography (OCT). These technologies provide another axis of image interpretation, an adjuvant to existing angiography, for successful CTO mapping and crossing.

IVUS, a well-established diagnostic and surveillance technique, has proven helpful in quantifying plaque morphology in coronary lesions¹⁴. Likewise, the optical properties of OCT enable a histologic

grade differentiation of atheromatous plaque across the wide spectrum of intimal disease^{15,16}. Rieber et al demonstrated the comparative advantages of OCT over IVUS in granular clarity when evaluating the lipid components of atheromas¹⁷. OCT consistently provided greater detail with respect to plaque composition, potentially leading to stronger indices for thin-cap fibroatheromas, the hallmark of vulnerable arterial disease leading to thrombosis, ischaemia and sudden cardiac death^{18,19}.

OCT, a relatively new technology, has already made large impacts in clinical fields²⁰ ranging from ophthalmology to gastroenterology²¹ and dermatology²². The discriminating optics range (10 µm) allows for an axial resolution in the order of 10-20 times greater than standard b-mode ultrasound imaging²³. Initial OCT applications in the endovascular arena have focused on diagnostic capabilities. Using OCT, clinicians are able to map and quantify arterial plaque structures *in vivo*. Regar et al showed the early usefulness of OCT in identifying a variety of plaque structures including thin-cap fibroatheromas²³. More recent studies have validated the superiority of OCT over IVUS for detection and characterisation of plaque structures including calcium, lipid pools and fibrin²⁴. As a result of the growing OCT momentum in cardiovascular medicine, the International Working Group for Intravenous OCT (IWG-IVOCT) has recently released a consensus standardisation and validation in an effort to improve further both diagnoses and treatment of atherosclerosis²⁵.

The Ocelot, the first-ever OCT-equipped CTO crossing catheter, utilises proprietary real-time OCT guidance and navigation for optimal intraluminal recanalisation. Ocelot's optical image resolution is 10 µm and incorporates 1,024 A-lines per image. These capabilities enable differentiation between various healthy arterial structures, including media and adventitia (layered structures), and diseased arterial walls, including plaque (non-layered structures). Either saline or contrast may be used via the local flush port to clear the imaging fibre when blood refracts the OCT image. Minimal flush (<1 cc) is adequate to displace blood artefact. The static flow environment of a CTO minimises the need for large flush volumes or proximal balloon dilation. Ocelot's working edge includes fixed spiral flutes with adjustable speed of rotation between 30, 45 and 60 rotations per minute and may be advantageous for crossing CTOs with live evaluation of luminal boundaries, improving recanalisation rates and ultimately limb salvage.

Methods

DEVICE DESCRIPTION

The Ocelot system consists of the Ocelot catheter and the Lightbox console (Avinger, Inc., Redwood City, CA, USA). The Ocelot system combines the use of advanced catheter design and real-time OCT for orientation and navigation of the catheter tip across chronic total occlusions.

OCELOT CATHETER

The Ocelot catheter is an over-the-wire device that is compatible with a 6 Fr sheath and a 0.014” guidewire, has a working length of

110 cm, and incorporates a 155 mm optical fibre. The catheter has a distal tip with spiral wedges (**Figure 1**) that can rotate clockwise and counterclockwise to facilitate a corkscrew-like action, allowing for advancement of the catheter across a CTO. Additionally, the distal tip is preshaped and contains radiopaque markers (**Figure 2**). The directional markers are used to orient the catheter during a CTO-crossing procedure. The middle directional marker (**Figure 3**) corresponds to the marker band located directly opposite the catheter's distal tip deflection. As such, the operator guides CTO crossing by positioning the middle marker over healthy tissue (layered structures) prior to distal advancement of the catheter. This orientation permits the catheter tip to engage the diseased segments, allowing for truer luminal crossing and a better chance of re-entry. Real-time, constant OCT imaging allows for reorientation of the

middle markers over layered structures for the duration of lesion crossing. This technique, referred to as OCT-assisted orientation, ensures that the deflected working edge of the Ocelot is pointing towards diseased tissue and away from healthy structures, as demarcated by the middle marker (**Figure 4**).

The Ocelot catheter handle (**Figure 5**) controls rotational orientation of the outer catheter shaft independent of the rotational movement of the distal tip. Rotating the distal tip also allows for real-time cross-sectional OCT of the surrounding vasculature during the procedure. For improved image quality, small amounts of saline can be flushed through the distal tip using the luer located at the proximal end of the catheter handle assembly. After the catheter enters the distal reconstitution site, it is removed, leaving the wire in place for therapeutic interventional treatment.




Figure 1. Distal tip of Ocelot catheter.




Figure 2. Pre-shaped distal tip deflection.




Figure 3. Directional markers on OCT system.

LIGHTBOX

The Ocelot catheter connects to the Lightbox console using an optical and electrical extension cable, which is used to provide an OCT image of the vasculature when coupled. The resulting image provides a qualitative intraluminal view of the vessel morphology as an adjunct to standard imaging techniques (e.g., fluoroscopy) (**Figure 6**). The Lightbox transmits light to the intraluminal environment




Figure 4. Ocelot OCT image with directional markers and layered structures.




Figure 5. Ocelot catheter handle.




Figure 6. Lightbox console (front and rear view).

through the optical fibre on the Ocelot catheter and receives and interprets the signal from the tissue using a PC-based processing system. The Lightbox processes the reflected light signal and displays a plot of tissue reflectivity versus depth, thereby creating an image of the surrounding vasculature. The Lightbox is supplied with two monitors, one monitor for physician viewing and the other used by the assisting Lightbox operator.

OCT images captured by the Lightbox are presented in two modes on the physician monitor: sector view and waterfall view. The sector view provides a cross-sectional intraluminal view during rotation of the distal catheter tip. The reflected OCT image is correlated to the distal tip's deflection angle using the middle marker, while the directional markers are fixed with respect to distal tip rotation. The resulting OCT image can then be used to help the physician orient the distal tip and navigate within the vessel as an adjunct to standard visualisation tools (e.g., fluoroscopy).

The Ocelot System (Avinger, Inc.) received CE mark approval in September 2011.

DATA COLLECTION AND DEFINITIONS

After informed consent, patients with CTO of their superficial femoral artery (SFA) were treated in the market preference testing (MPT) study between September 28, 2011, and December 9, 2011, at three European centres. Ocelot catheter and OCT image performance data were based on direct physician feedback during each case and were rated on a Likert scale from 1-5 where: 1=unacceptable; 2=poor; 3=acceptable; 4=good; 5=excellent. Ocelot catheter performance feedback included: catheter introduction into sheath, device trackability, device pushability, tip deflection, ability to stay in true lumen, visibility of radiopaque markers, visibility of catheter tip under fluoroscopy. OCT image quality feedback included: visualisation of vessel wall layered structures, visualisation of vessel wall non-layered structures, visualisation of side branches, ability of saline flush to displace blood adequately for visualisation, ability to orient the device using the middle marker.

Definitions for key safety and efficacy measures were based on the CONNECT study²⁶. The key safety measure was a composite endpoint capturing procedural major adverse events (MAE), clinically significant perforations requiring treatment, embolisations resulting in distal ischaemia and/or requiring rescue intervention, and grade C or greater dissections resulting in extraluminal, persisting extravasation. The key efficacy measure was defined as successful CTO crossing using the Ocelot catheter as identified by guidewire placement in the distal true lumen and confirmed by angiography (allowing use of other crossing/assist/re-entry devices). Secondary efficacy measures included: use of re-entry/assist devices in lesions that were successfully crossed, procedure time (defined as initial stick to last wire out), CTO crossing time, fluoroscopy time, contrast dose and amount of saline flush (through Ocelot for visualisation only).

STATISTICAL ANALYSES

The resulting values are presented as mean \pm standard deviation. Analyses comparing Ocelot MPT data to CONNECT study data were made using Fisher's exact test for categorical data and Student's t-test for independent groups for continuous data. Statistical analyses were performed using SAS for Windows software, version 9.1 or higher (SAS Institute Inc., Cary, NC, USA).

Results

Thirty-three patients with confirmed CTO (99-100% stenosed by visual estimate) of SFA were treated in this study. Mean patient age was 70 \pm 9 years old and 62% of patients were male. Eighty-one percent of patients were prior smokers, 81% had a history of hypertension, and 70% had previous non-coronary interventions. Complete patient demographics and comorbidities are provided in **Table 1**.

Baseline lesion characteristics are presented in **Table 2**. Of the 33 lesions, 18% were located at the ostium, 27% were proximal, 42% were mid, 6% were mid to distal and 6% were distal pertaining to the segment within the SFA. A mildly calcified lesion was identified in 67% of patients while 21% had severely calcified lesions. Mean lesion length and diameter by visual estimate were 205 \pm 124 mm and 5.5 \pm 0.7 mm, respectively. A *de novo* lesion was identified in 94% of patients.

Table 1. Baseline demographics and comorbidities.

Demographics and comorbidities	Ocelot MPT N=33 (% , m/n)
Age, yrs (mean \pm SD)	71 \pm 9 (33)
Male gender, %	62% (21/33)
History of smoking, %	63% (19/30)
History of hypertension, %	81% (27/33)
History of MI, %	12% (4/33)
History of diabetes, %	42% (14/33)
Prior CABG, %	15% (5/33)
Prior PCI, %	24% (8/33)
Prior non-coronary interventions, %	70% (23/33)

Table 2. Baseline general information and lesion characteristics.

Target lesion		Ocelot MPT N=33 (% , m/n)
Access site	Femoral	97% (32/33)
	Popliteal	3% (1/33)
Procedural approach	Contralateral	94% (31/33)
	Ipsilateral	6% (2/33)
Lesion location		
Target vessel laterality	Left	58% (19/33)
	Right	42% (14/33)
SFA		100% (33/33)
Ostial		18% (6/33)
Proximal		27% (9/33)
Mid		42% (14/33)
Mid to distal		6% (2/33)
Distal		6% (2/33)
Calcification		
Mild		67% (22/33)
Moderate		6% (2/33)
Moderate to severe		6% (2/33)
Severe		21% (7/33)
Other characteristics		
Avg. lesion length (mm)		205±124 (n=32)*
Avg. reference vessel diameter (mm)		5.5±0.7 (n=33)*
Lesion type	<i>de novo</i>	94% (31/33)
	restenotic	6% (2/33)
99-100% stenosis pre-procedure		100% (33/33)

*Data presented as mean±standard deviation (n)

PROCEDURAL RESULTS

All Ocelot MPT patients (100%) achieved the key safety measure, freedom from procedural MAE (clinically significant perforations, embolisations and grade C or greater dissections). Thirty-one out of 33 patients (94%) achieved the key efficacy measure, successful CTO crossing with the Ocelot system. In the two patients where the CTO was not successfully crossed, the lesions were identified as severely calcified along the entire SFA and subsequent attempts to cross the CTO with other crossing devices also failed (i.e., Wildcat; Avinger, Inc., Quick-Cross®; Spectranetics Corp., Colorado Springs, CO, USA, stiff guidewire, and Outback®; Cordis, Johnson & Johnson, Warren, NJ, USA). Both patients were subsequently referred for elective surgery. Key safety and efficacy results are provided in **Table 3**.

Ocelot MPT procedural results are provided in **Table 4**. Only five of 31 lesions successfully crossed (16.1%) required the use of a re-entry/assist device. The mean total procedure time was 83±28 minutes and the mean Ocelot catheter crossing time was 17±16 minutes for cases where the CTO was successfully crossed. While OCT visualisation was able to define the true lumen target site in these instances, the anatomic distortion of severely calcified distal caps proved prohibitive to catheter navigation, displacing the catheter into a subintimal plane

Table 3. Ocelot MPT procedural results.

Data parameter	Ocelot MPT N=33 (%, m/n or mean±SD)
Key safety measure	
Freedom from MAE*	100% (33/33)
Key efficacy measure	
Overall crossing success	94% (31/33)
Secondary efficacy measures	
Use of re-entry device (assisted crossing)#	16% (5/31)
Procedural time (min)	83±28 (n=31)
OCT crossing time (min)#	17±16 (n=31)
Ocelot alone crossing time (min)	13±10 (n=26)
Ocelot + assisted crossing time (min)	39±27 (n=5)
Fluoroscopy time (min)	28±12 (n=31)
Contrast dose (ml)	132±75 (n=31)
Flush amount (ml)	49±37 (n=30)
Additional data	
Post-procedure treatment type#	100% (31/31)
Angioplasty	29% (9/31)
Angioplasty + stent	65% (20/31)
Atherectomy + angioplasty + stent	3% (1/31)
Stent	3% (1/31)
Average post-procedure stenosis (%)#	4.9±6.5 (31)

*MAE defined as composite of clinically significant perforations, embolisations, and grade C or greater dissections; #Only includes patients where CTO was successfully crossed (31 of 33 patients)

requiring an assist device for re-entry. During two cases, the distal lesion cap was so severely calcified that true lumen re-entry was altogether unsuccessful. When no re-entry/assist devices were needed to cross the CTO (26/31), Ocelot crossing time decreased to 13±10 minutes. Mean contrast dose was 132±75 mL and mean flush amount was 49±37 mL. Subsequent therapies for those lesions crossed were consistent with standard EU practice and included: 29% angioplasty alone, 65% angioplasty and stent, 3% atherectomy, angioplasty and stent, and 3% treated with a stent alone. The average post-procedure stenosis was 4.9±6.5 percent, defined by angiographic imaging.

OCELOT CATHETER PERFORMANCE RATINGS

Overall physician feedback on the performance of the Ocelot catheter is presented in **Table 5** and was positive, with an 87% average rating of excellent or good across the seven catheter categories. Catheter performance was physician-rated as excellent or good in the following percentages: 100% visibility of catheter tip under fluoroscopy, 97% catheter introduction into sheath, 94% catheter trackability, 82% ability to stay in true lumen, 82% tip deflection, 76% visibility of radiopaque markers, and 73% device pushability. The ability to stay in the true lumen was rated as poor or unacceptable in 4 instances (poor, n=3, and unacceptable, n=1). In two of the cases, re-entry devices were used and failed to gain re-entry or cross the lesion

Table 4. Catheter feedback rating summary.

Ratings	Introduction into sheath N=33	Catheter trackability N=33	Catheter pushability N=33	Tip deflection N=33	Ability to stay in true lumen N=33	Visibility of radiopaque marker N=33	Catheter tip visibility N=33
5=excellent	70% (23/33)	55% (18/33)	24% (8/33)	12% (4/33)	61% (20/33)	55% (18/33)	79% (26/33)
4=good	27% (9/33)	39% (13/33)	49% (16/33)	70% (23/33)	21% (7/33)	21% (7/33)	21% (7/33)
3=acceptable	3% (1/23)	6% (2/33)	27% (9/33)	18% (6/33)	6% (2/33)	21% (7/33)	0% (0/33)
2=poor	0% (0/33)	0% (0/33)	0% (0/33)	0% (0/33)	9%* (3/33)	3%* (1/33)	0% (0/33)
1=unacceptable	0% (0/33)	0% (0/33)	0% (0/33)	0% (0/33)	3%* (1/33)	0% (0/33)	0% (0/33)
Combined 5=excellent / 4=good	97% (32/33)	94% (31/33)	73% (24/33)	82% (27/33)	82% (27/33)	76% (25/33)	100% (33/33)
Combined 2=poor / 1=unacceptable	0% (0/33)	0% (0/33)	0% (0/33)	0% (0/33)	12% (4/33)	3% (1/33)	0% (0/33)

*The ability to stay in true lumen was rated as poor or unacceptable in 4 instances (4/33, 12% [9% poor, 3% unacceptable]). In two cases, re-entry devices were used and failed to gain re-entry or cross the lesion successfully. In the remaining two cases, an upstream catheter and needle and a Pioneer device were used to gain re-entry back into the true lumen. *The ability to visualise the radiopaque marker was rated as poor in 1 instance (1/33, 3%). In this instance, the operator indicated it was difficult to determine the orientation of the marker band under fluoroscopy.

successfully and, in the remaining two cases, an upstream catheter and needle and a Pioneer® device (Medtronic CardioVascular, Santa Rosa, CA, USA) were used to gain re-entry back into the true lumen.

OCT PERFORMANCE RATINGS

Overall physician feedback on the performance of OCT is presented in **Table 5** and was positive with an 86% average rating of excellent or good across the five OCT categories. OCT performance was physician rated as excellent or good in the following percentages: 100% visualisation of side branches, 100% ability to orient the device using the middle marker, 97% ability of saline flush to displace blood adequately for visualisation, 87% visualisation of vessel wall layered structures, and 54% visualisation of vessel wall non-layered structures. The only poor rating occurred in the category of ability to visualise the radiopaque marker. In this instance, the operator indicated that it was difficult to determine the orientation of the marker band under fluoroscopy.

OCELOT COMPARED TO CONNECT (NON-OCT-GUIDED WILDCAT CATHETER)

Analyses comparing the OCT-guided Ocelot catheter to the non-OCT-guided Wildcat catheter (data collected in the CONNECT study¹⁶) are shown in **Table 6**. The key safety measure, freedom

from clinically significant perforations, embolisations or grade C or greater dissections, was improved with the Ocelot, although not significantly (100% vs. 95.2%, $p=0.58$). The key efficacy measure, ability to cross the CTO, was also improved with the Ocelot as compared to the Wildcat, although again not significant (94% vs. 83%, $p=0.72$). When the ability to cross the CTO was analysed in only those Ocelot patients who met the inclusion and exclusion criteria of the CONNECT study ($n=26$), a 100% CTO crossing rate was observed (100% vs. 83%, $p=0.11$). Seven patients treated with the Ocelot would have been excluded from the CONNECT study due to a severely calcified target vessel.

Procedure time and contrast dose were significantly reduced with the use of the OCT-guided Ocelot catheter (**Table 6**). Procedure time decreased from 124 ± 74 minutes with the Wildcat, to 83 ± 28 minutes with the Ocelot, $p<0.0001$, while contrast dose decreased from 248 ± 156 millilitres to 132 ± 75 millilitres, respectively, $p<0.001$. There was not a statistically significant change in fluoroscopy time between the Ocelot (28 ± 12 minutes) and Wildcat (30 ± 15 minutes), $p=0.37$.

Discussion

CTOs represent the most technically challenging subset of PAD lesions. The current shift from primary surgical bypass to

Table 5. OCT feedback rating summary*.

Ratings	Ability to visualise layered structures N=32	Ability to visualise non-layered structures N=32	Visualisation of side branches [#] N=25	Flush adequately displaces blood for visualisation [¶] N=32	Ability to orient using the middle marker N=32
5=excellent	56% (18/32)	13% (4/32)	32% (8/25)	63% (20/31)	28% (9/32)
4=good	25% (8/32)	41% (13/32)	68% (17/25)	28% (9/31)	59% (19/32)
3=acceptable	16% (5/32)	47% (15/32)	0% (0/25)	6% (2/31)	13% (4/32)
2=poor	3% (1/32) [‡]	0% (0/32)	0% (0/25)	0% (0/31)	0% (0/32)
1=unacceptable	0% (0/32)	0% (0/32)	0% (0/25)	0% (0/31)	0% (0/32)
Combined 5=excellent / 4=good	81% (26/32)	54% (17/32)	100% (32/32)	97% (29/31)	100% (32/32)
Combined 2=poor / 1=unacceptable	3% (1/32)	0% (0/32)	0% (0/32)	0% (0/31)	0% (0/32)

*OCT feedback data not captured for 1 patient; [#]N/A in 7 cases, no side branches present; [¶]Not collected in 1 patient; [‡]Ability to visualise layered structures rated poor in 1 instance

Table 6. Ocelot vs. Wildcat (CONNECT).

Data parameter	Ocelot MPT N=33 (%, m/n or mean±SD)	CONNECT (Wildcat) N=84 (%, m/n or mean±SD)	p-value
Key safety measure			
Freedom from MAE*	100% (33/33)	96.2% (80/84)	0.58
Key efficacy measure			
Overall crossing success	94% (31/33) [#]	89% (75/84)	0.72
Secondary efficacy measures			
Use of re-entry device (assisted crossing) [‡]	16% (5/31)	17.85% (15/84)	0.825
Procedural time (min)	83±28 (n=31)	124±74 (n=82)	<0.0001
OCT crossing time (min) [#]	17±16 (n=31)	N/A	N/A
Ocelot alone crossing time (min)	13±10 (n=26)		
Ocelot+assisted crossing time (min)	39±27 (n=5)		
Fluoroscopy time (min)	28±12 (n=31)	30±15 (n=78)	0.37
Contrast dose (ml)	132±75 (n=31)	248±156 (n=84)	<0.0001
Flush amount (ml)	49±37 (n=30)	N/A	N/A
*MAE defined as composite of clinically significant perforations, embolisations, and grade C or greater dissections; [#] 100% crossing rate if patients who met CONNECT study eligibility criteria are included (n=26); [‡] Only includes patients where CTO was successfully crossed (31 of 33 patients)			

endovascular intervention is the result of improved technologies and post-crossing patency rates^{27,28}. This may translate into shorter hospital stays, earlier patient recovery and ambulation, and reduced exposure to potential surgical and anaesthetic complications. Ultimately, the expanded availability of successful endovascular modalities for limb salvage and revascularisation has the potential to reduce associated healthcare costs, all the while improving patient satisfaction and outcomes³⁰.

Currently available CTO crossing catheters utilise a variety of technologies. The Crosser™ CTO Recanalisation System (FlowCardia Inc.[Bard], Sunnyvale, CA, USA) transmits high frequency ultrasound mechanical vibrations to the distal tip with a reported recanalisation success rate of 75% (n=12) with no clinically significant perforations²⁹. The rotational diamond-coated tip of the TruePath™ CTO device (Boston Scientific, Natick, MA, USA) demonstrated an 80% CTO crossing with 98.8% freedom from perforation (n=85)³⁰. The Cordis Frontrunner XP (Cordis Corp., Miami, FL, USA) has a distal tip with hinged jaws, creating micro-dissections as it traverses CTOs. In one report of 26 SFA CTOs, 65% were successfully crossed³². Mossop et al reported a 91% success rate in 36 patients following failed guidewire crossing³³. The CrossBoss™ (BridgePoint, Plymouth, MN, USA), intended for coronary procedures, employs an atraumatic rounded tip with fast spin technique for crossing. The Wildcat (Avinger, Inc.) contains a distal element housing spiral wedges which, when engaged, actively corkscrew through CTOs with reported 89% efficacy (n=84)¹⁶.

IVUS and OCT are established diagnostic imaging modalities for both the peripheral and coronary vascular beds³⁴. However, there are limited applications adapted for therapeutic intervention that would allow for comparison with the aforementioned non-imaging devices. Saket et al³⁵ highlighted the CrossPoint TransAccess catheter (Medtronic Inc., Santa Rosa, CA, USA), a re-entry device with an intravascular ultrasound (IVUS) transducer for true lumen

targeting. This technique was used successfully in seven patients following subintimal crossing and angioplasty.

When compared with IVUS, OCT's infrared spectrum consistently demonstrates more reliable histopathological analyses of plaque morphology³⁶. Coronary findings have been extrapolated to peripheral disease with OCT validations of histologic correlation via imagery. Therefore, OCT has been recognised as the technology that represents the most promising conduit for image-guided therapy with the advent of percutaneous therapies such as cutting balloon angioplasty, intravascular cryoablation, CTO crossing, and atherectomy³⁷.

The Ocelot catheter, a new device aimed at crossing complex SFA CTOs, represents the first time that OCT has been incorporated into the distal tip of a peripheral CTO crossing catheter. The proprietary common path, OCT, enables the length variance and tip flexibility necessary for successful navigation of CTOs. Accordingly, the rotating optical fibre provides continuous A-scans, allowing for an extrapolated image that lends itself to real-time clinical decision making.

This paper reports on a multicentre European experience, where the Ocelot catheter was used to treat 33 patients with confirmed SFA CTOs with 100% freedom from major adverse events and 94% successful guidewire placement in the distal true lumen. Important secondary measures demonstrated statistically significant reductions in procedural time and contrast dose when compared to the Wildcat catheter, a non-OCT-guided CTO crossing catheter. While further study is needed to validate these findings, initial outcomes suggest an added benefit of healthcare cost savings related to shortened procedure time and improved patient safety with a reduction in delivered contrast volumes.

In accordance with previous investigational trials, difficulty with CTO crossing appeared to correlate with lesion-specific calcium burdens²⁹. In the two patients where the CTO was not successfully crossed, the lesions were identified as severely calcified along the

entire SFA and subsequent attempts to cross the CTO with other crossing devices also failed.

The Ocelot is a safe and effective tool moving forward in the treatment of complex femoropopliteal CTOs. The easy adoption and use of both the catheter and OCT technology will likely change the landscape of treating PAD. With an additional axis of imaging, Ocelot and OCT provide interventionalists with a safe, reliable, and reproducible tool for intraluminal CTO revascularisation and limb salvage.

Limitations

A few limitations of this study are: 1) small sample size; 2) non-randomised patient population; and 3) selection bias for patients with failed initial guidewire crossing of proximal cap, albeit unfavourable for the Ocelot cohort as this probably resulted in longer, more complex, and calcified lesions. Additional experience is needed to validate the use of this device in endovascular therapy and will be forthcoming via the CONNECT II study.

Conflict of interest statement

J.C. George, J. Pigott and G. Robertson are consultants for Avinger, Inc., and have investments in Avinger, Inc. M. Selmon and B. Reimers are consultants for Avinger, Inc. G. Shrikhander has received an educational grant from Avinger, Inc., and has investments in Avinger, Inc. The other authors have no conflicts of interest to declare.

References

- Shammas NW. Epidemiology, classification, and modifiable risk factors of peripheral arterial disease. *Vasc Health Risk Manag.* 2007;3:229-34.
- Allie DE, Hebert CJ, Lirtzman MD, Wyatt CH, Keller VA, Khan MH, Khan MA, Fail PS, Vivekananthan K, Mitran EV, Allie SE, Chaisson G, Stagg SJ, Allie AA, McElderry MW, Walker CM. Critical limb ischemia: a global epidemic. A critical analysis of current treatment unmasks the clinical and economic costs of CLI. *EuroIntervention.* 2005;1:75-84.
- Weitz JI, Byme J, Clagett GP, Farkouh ME, Porter JM, Sackett DL, Strandness DE Jr, Taylor LM. Diagnosis and treatment of chronic arterial insufficiency of the lower extremities: a critical review. *Circulation.* 1996;94:3026-49.
- Jamieson C. The definition of critical limb ischemia of a limb. *Br J Surg.* 1982;69(suppl):S1.
- Wolfe JH, Wyatt MG. Critical and subcritical ischaemia. *Eur J Vasc Endovasc Surg.* 1997;13:578-82.
- Albers M, Romiti M, Brochado-Neto FC, Pereira CA. Meta-analysis of alternate autogenous vein bypass grafts to infrapopliteal arteries. *J Vasc Surg.* 2005;42:449-55.
- Wölfle KD, Bruijnen H, Loeprecht H, Rümnapf G, Schweiger H, Grabitz K, Sandmann W, Lauterjung L, Largiader J, Erasmi H, Kasprzak PM, Raithel D, Allenberg JR, Lauber A, Berlakovich GM, Kretschmer G, Hepp W, Becker HM, Schulz A. Graft patency and clinical outcome of femorodistal arterial reconstruction in diabetic and non-diabetic patients: results of a multi-centre comparative analysis. *Eur J Vasc Endovasc Surg.* 2003;25:229-34.
- Pomposelli FB, Kansal N, Hamdan AD, Belfield A, Sheahan M, Campbell DR, Skillman JJ, Logerfo FW. A decade of experience with dorsal pedis artery bypass: analysis of outcome in more than 1,000 cases. *J Vasc Surg.* 2003;37:165-71.
- Hunink MG, Wong JB, Donaldson MC, Meyerovitz MF, de Vries J, Harrington DP. Revascularization for femoropopliteal disease. A decision and cost-effective analysis. *JAMA.* 1995;274:165-71.
- Rabellino M, Zander T, Baldi S, Garcia Nielsen L, Aragon-Sanchez FJ, Zerolo I, Llorens R, Maynar M. Clinical follow-up in endovascular treatment for TASC C-D lesions in femoro-popliteal segment. *Catheter Cardiovasc Interv.* 2009;73:701-5.
- Bolia A, Miles KA, Brennan J, Bell PR. Percutaneous transluminal angioplasty of occlusions of the femoral and popliteal arteries by subintimal dissection. *Cardiovasc Intervent Radiol.* 1990;13:357-63.
- Met R, Van Lienden KP, Koelemay MJ, Bipat S, Legemate DA, Reekers JA. Subintimal angioplasty for peripheral arterial occlusive disease: a systemic review. *Cardiovasc Intervent Radiol.* 2008;31:687-97.
- Adam DJ, Beard JD, Cleveland T, Bell J, Bradbury AW, Forbes JF, Fowkes FG, Gillepsie I, Ruckley CV, Raab G, Storkey H; BASIL trial participants. Bypass versus angioplasty in severe ischemia of the leg (BASIL): multicentre, randomized controlled trial. *Lancet.* 2005;336:1925-34.
- Vaina S, Stefanadis C. Detection of the vulnerable coronary atheromatous plaque. Where are we now? *Int J Cardiovasc Intervent.* 2005;7:75-87.
- Tearney GJ, Brezinski ME, Bouma BE, Boppart SA, Pitris C, Southern JF, Fujimoto JG. In vivo endoscopic optical biopsy with optical coherence tomography. *Science.* 1997;276:2037-9.
- Fujimoto JG, Pitris C, Boppart SA, Brezinski ME. Optical coherence tomography: an emerging technology for biomedical imaging and optical biopsy. *Neoplasia.* 2000;2:9-25.
- Rieber J, Meissner O, Babaryka G, Reim S, Oswald M, Koenig A, Schiele TM, Shapiro M, Theisen K, Reiser MF, Klaus V, Hoffmann U. Diagnostic accuracy of optical coherence tomography and intravascular ultrasound for the detection and characterization of atherosclerotic plaque composition in ex-vivo coronary specimens: a comparison with histology. *Coron Artery Dis.* 2006;17:425-30.
- Jang IK, Bouma BE, Kang DH, Park SJ, Park SW, Seung KB, Choi KB, Shishkov M, Schlendorf K, Pomerantsev E, Houser SL, Aretz HT, Tearney GJ. Visualization of coronary atherosclerotic plaques in patients using optical coherence tomography: comparison with intravascular ultrasound. *J Am Coll Cardiol.* 2002;39:604-9.
- Sawada T, Shite J, Garcia-Garcia HM, Shinke T, Watanabe S, Otake H, Matsumoto D, Tanino Y, Ogasawara D, Kawamori H, Kato H, Miyoshi N, Yokoyama M, Serruys PW, Hirata K. Feasibility of combined use of intravascular ultrasound radiofrequency data analysis and optical coherence tomography for detecting thin-cap fibroatheroma. *Eur Heart J.* 2008;29:1136-46.
- Pinto TL, Waksman R. Clinical applications of optical coherence tomography. *J Interv Cardiol.* 2006;19:566-73.
- Li XD, Boppart SA, Van Dam J, Mashimo H, Mutinga M, Drexler W, Klein M, Pitris C, Krinsky ML, Brezinski ME, Fujimoto JG.

Optical coherence tomography: advanced technology for the endoscopic imaging of Barrett's esophagus. *Endoscopy*. 2000;32:921-30.

22. Welzel J, Lankenau E, Birngruber R, Engelhardt R. Optical coherence tomography of the human skin. *J Am Acad Dermatol*. 1997;37:958-63.

23. Farooq MU, Khasnis A, Majid A, Kassab MY. The role of optical coherence tomography in vascular medicine. *Vasc Med*. 2009;14:63-71.

24. Regar E, Schaar JA, Mont E, Virmani R, Serruys PW. Optical Coherence Tomography. *Cardiovasc Radiat Med*. 2003;4:189-204.

25. Tearney GJ, Regar E, Akasaka T, Adriaenssens T, Barlis P, Bezerra HG, Bouma B, Bruining N, Cho JM, Chowdhary S, Costa MA, de Silva R, Dijkstra J, Di Mario C, Dudek D, Falk E, Feldman MD, Fitzgerald P, Garcia-Garcia HM, Gonzalo N, Granada JF, Guagliumi G, Holm NR, Honda Y, Ikeno F, Kawasaki M, Kochman J, Koltowski L, Kubo T, Kume T, Kyono H, Lam CC, Lamouche G, Lee DP, Leon MB, Maehara A, Manfrini O, Mintz GS, Mizuno K, Morel MA, Nadkarni S, Okura H, Otake H, Pietrasik A, Prati F, Räber L, Radu MD, Rieber J, Riga M, Rollins A, Rosenberg M, Sirbu V, Serruys PW, Shimada K, Shinke T, Shite J, Siegel E, Sonoda S, Suter M, Takarada S, Tanaka A, Terashima M, Thim T, Uemura S, Ughi GJ, van Beusekom HM, van der Steen AF, van Es GA, van Soest G, Virmani R, Waxman S, Weissman NJ, Weisz G; International Working Group for Intravascular Optical Coherence Tomography (IWG-IVOC). Consensus standards for acquisition, measurement, and reporting of intravascular optical coherence tomography studies: a report from the international working group for intravascular optical coherence tomography standardization and validation. *J Am Coll Cardiol*. 2012;59:1058-72.

26. Pigott JP. Intraluminal Versus Subintimal Approach to Chronic Total Occlusions: the pros and cons of each approach for challenging femoropopliteal crossing. *Endovascular Today*. 2011;10:73-6.

27. Hong MS, Beck AW, Nelson PR. Emerging national trends in the management and outcomes of lower extremity peripheral arterial disease. *Ann Vasc Surg*. 2011;25:44-54.

28. DeRubertis BG, Faries PL, McKinsey JF, Chaer RA, Pierce M, Karwowski J, Weinberg A, Nowygrod R, Morrissey NJ, Bush HL, Kent KC. Shifting paradigms in the treatment of lower extremity vascular disease: a report of 1000 percutaneous interventions. *Ann Surg*. 2007;246:415-22; discussion 422-4.

29. Culverwell AD, Tapping CR, Ettles DF, Kessel D. Patient experience, pain, and quality of life after lower limb angioplasty: a multisite prospective cohort study. *Cardiovasc Intervent Radiol*. 2012;35:788-94.

30. Gandini R, Volpi T, Pipitone V, Simonetti G. Intraluminal recanalization of long infrainguinal chronic total occlusions using the Crosser System. *J Endovasc Ther*. 2009;16:23-7.

31. <http://www.bostonscientific.com/truepath/clinical-data.html>? Accessed 2/10/12

32. Charalambous N, Schäfer PJ, Trentmann J, Hümme TH, Stöhring C, Müller-Hülsbeck S, Heller M, Jahnke T. Percutaneous intraluminal recanalization of long, chronic superficial femoral and popliteal occlusions using the Frontrunner XP CTO device: a single-center experience. *Cardiovasc Intervent Radiol*. 2010;33:25-33.

33. Mossop PJ, Amukotuwa SA, Whitborn RJ. Controlled blunt microdissection for percutaneous recanalization of lower limb arterial chronic total occlusions: a single center experience. *Catheter Cardiovasc Interv*. 2006;68:304-10.

34. Fujimoto JG, Boppart SA, Tearney GJ, Bouma BE, Pitris C, Brezinski ME. High resolution in vivo intra-arterial imaging with optical coherence tomography. *Heart*. 1999;82:128-33.

35. Saket R, Razavi MK, Padidar A, Kee S, Sze DY, Dake MD. Novel intravascular ultrasound-guided method to create transintimal arterial communications: initial experience in peripheral occlusive disease and aortic dissection. *J Endovasc Ther*. 2004;11:274-80.

36. Brezinski ME, Tearney GJ, Weissman NJ, Boppart SA, Bouma BE, Hee MR, Weyman AE, Swanson EA, Southern JF, Fujimoto JG. Assessing atherosclerotic plaque morphology: comparison of optical coherence tomography and high frequency intravascular ultrasound. *Heart*. 1997;77:397-403.

37. Meissner OA, Rieber J, Babaryka G, Oswald M, Reim S, Siebert U, Redel T, Reiser M, Mueller-Lisse U. Intravascular optical coherence tomography: comparison with histopathology in atherosclerotic peripheral artery specimens. *J Vasc Interv Radiol*. 2006;17:343-9.

Online data supplement

Moving image 1. Animated picture of the Ocelot.