

Deformation of transcatheter heart valves with mitral valve-in-valve

Miho Fukui¹, MD, PhD; Atsushi Okada², MD, PhD; Marcus R. Burns³, DNP; Hiroto Sato², MD, PhD; Kiahltone R. Thao², BA; Cheng Wang², MD; Hideki Koike¹, MD, PhD; Nadira Hamid^{2,3}, MD; Maurice Enriquez-Sarano², MD; John R. Lesser^{2,3}, MD; João L. Cavalcante^{1,3}, MD; Paul Sorajja^{2,3}, MD; Vinayak N. Bapat^{2,3*}, MD

1. Cardiovascular Imaging Research Center and Core Lab, Minneapolis Heart Institute Foundation, Minneapolis, MN, USA; 2. Valve Science Center, Minneapolis Heart Institute Foundation, Minneapolis, MN, USA; 3. Allina Health Minneapolis Heart Institute – Minneapolis, Abbott Northwestern Hospital, Minneapolis, MN, USA

This paper also includes supplementary data published online at: <https://eurointervention.pronline.com/doi/10.4244/EIJ-D-23-00614>

KEYWORDS

- degenerative valve
- imaging modalities
- mitral valve disease
- valve-in-valve

Abstract

Background: The use of oversizing in mitral valve-in-valve (MViV) procedures can lead to non-uniform expansion of transcatheter heart valves (THV). This may have implications for THV durability.

Aims: The objective of this study was to assess the extent and predictors of THV deformation in MViV procedures.

Methods: We examined 33 patients who underwent MViV with SAPIEN prostheses. The extent of THV deformation (deformation index, eccentricity, neosinus volume, asymmetric leaflet expansion and vertical deformation) and hypoattenuating leaflet thickening (HALT) were assessed using cardiac computed tomography (CT), performed prospectively at 30 days post-procedure. For descriptive purposes, the THV deformation index was calculated, with values >1.00 representing a more hourglass shape.

Results: Non-uniform underexpansion of THV was common after MViV implantation, with a median expansion area of 74.0% (interquartile range 68.1-84.1) at the narrowest level and a THV deformation index of 1.21 (1.13-1.29), but circularity was maintained with eccentricity ranging from 0.24 to 0.28. The degree of oversizing was a key factor associated with greater underexpansion and a higher deformation index ($\beta=-0.634$; $p<0.001$; $\beta=0.594$; $p<0.001$, respectively). Overall, the incidence of HALT on the 30-day postprocedural CT was 27.3% (9 of 33). Most patients (32 of 33) were on anticoagulation therapy, but the prothrombin time and international normalised ratio (PT-INR) at the time of the CT scan was <2.5 in 23 of 32 patients. Among patients with a PT-INR of <2.5, HALT was predominantly observed with a high THV deformation index of ≥ 1.18 .

Conclusions: THV deformation, i.e., underexpansion and an hourglass shape, commonly occurs after MViV implantation and is negatively affected by excessive oversizing. Optimising THV expansion during MViV could potentially prevent HALT.

*Corresponding author: Minneapolis Heart Institute, Abbott Northwestern Hospital, Valve Science Center, Minneapolis Heart Institute Foundation, 800 E 28th Street, Suite 300, Minneapolis, MN, 55407, USA. E-mail: vnbapat@yahoo.com

Abbreviations

BVF	bioprosthetic valve fracture
CT	computed tomography
HALT	hypoattenuating leaflet thickening
MViV	mitral valve-in-valve
PT-INR	prothrombin time and international normalised ratio
SHV	surgical heart valve
THV	transcatheter heart valve
True ID	true internal dimension

Introduction

Transcatheter mitral valve-in-valve (MViV) implantation is an established therapeutic option for selected patients with surgical heart valve (SHV) degeneration, with demonstrated safety and effectiveness¹⁻⁴. While cases of MViV are commonly performed with balloon-expandable SAPIEN (Edwards Lifesciences) devices^{1,3} which were originally designed for the aortic position, the original intended design of this transcatheter heart valve (THV) prosthesis is to maintain a uniform, cylindrical shape for its ideal function.

During MViV procedures, oversizing the THV, with respect to the true internal diameter (True ID) of the SHV, and flaring of the outflow are used in order to ensure fixation and reduce the risk of embolisation^{1,5,6}. This can result in a non-uniform underexpansion of the THV, with potential consequences of pinwheeling, shear stress on the leaflets, increased blood stasis and flow stagnation as shown in previous *in vitro* testing⁷⁻¹⁰. This could affect THV durability. In a recent study, we observed that such prosthesis deformation was also related to an increased risk of hypoattenuating leaflet thickening (HALT) in patients treated with transcatheter aortic valve replacement (TAVR)¹¹. In addition, adverse haemodynamic effects due to excessive underexpansion of the functional portion of the THV were observed after aortic valve-in-valve procedures¹². Given these *in vitro* and *in vivo* findings, insight into the extent of THV deformation following MViV implantation and its relation to HALT and haemodynamics is important. This is particularly relevant as post-dilatation, with or without fracture, to optimise THV expansion is utilised less frequently in the current practice of MViV^{1,5,13,14}.

In our investigation, we performed the THV assessment with cardiac computed tomography (CT), with detailed examination of THV expansion and leaflet alignment to evaluate the extent and the predictors of THV deformation with MViV.

Methods

STUDY DESIGN AND PATIENT POPULATION

Starting on 1 July 2015, a prospective internal registry was initiated at the Minneapolis Heart Institute at Abbott Northwestern Hospital (Minneapolis, MN, USA) for patients who received THV implantation, whereby a retrospective electrocardiogram (ECG)-gated contrast-enhanced cardiac CT examination was performed to detect HALT and assess leaflet mobility at 30 days post-procedure. All patients were consented to participate unless there was a medical (e.g., high risk of contrast-induced nephropathy) or

social (e.g., unable to travel) contraindication to the performance of the study, as previously described¹⁵. This current study included a subset of patients with mitral SHV degeneration who had been successfully treated with a ViV procedure using either a balloon-expandable SAPIEN 3 or 3 Ultra (S3/S3 Ultra) prosthesis (sizes 23, 26, or 29 mm) and who had undergone such CT examinations between 1 July 2015 and 31 January 2023. All patients underwent a comprehensive clinical evaluation by a Heart Team, and the MViV procedure was performed according to established procedural strategies, as previously shown^{1,5,16}. The transeptal approach was the preferred approach, and the transapical approach was used only in the earlier part of our experience. Selection of the THV size and use of pre-/postdilatation were left to the operator's discretion. Patients who converted to surgery, had THV migration, needed a second THV implantation or who suffered perioperative death that precluded subsequent CT imaging were excluded.

The present study was conducted in accordance with the Declaration of Helsinki and approved by the Allina Institutional Review Board. All patients provided informed consent for the use of their medical records for research purposes.

INSTITUTIONAL ANTICOAGULATION PROTOCOL

Prior to the MViV procedure, patients on anticoagulation therapy discontinued anticoagulation for 4 days (for vitamin K antagonists [i.e., warfarin]) or for 48 hours (for direct oral anticoagulation). All patients routinely received a loading dose of aspirin (325 mg) on the day of the procedure. During the procedure, weight-adjusted unfractionated heparin was given to achieve an activated clotting time of ≥ 250 seconds. Following MViV implantation, patients received warfarin with an international normalised ratio (INR) goal of 2.5 – for the first 3 months with oral aspirin, followed by lifelong aspirin therapy. For patients with indications for long-term anticoagulation, a combination of warfarin or direct oral anticoagulation and single antiplatelet therapy for 3 to 6 months was administered, followed by lifelong anticoagulation therapy based on the risks of thromboembolism and bleeding.

CT IMAGE ACQUISITION AND ANALYSIS

Contrast-enhanced ECG-gated cardiac CT was performed at 30 days post-MViV procedure using previously published acquisition protocols (**Supplementary Appendix 1**). All CT imaging analyses were performed on a dedicated post-processing workstation equipped with Circle cvi42 (Version 5.13.7; Circle Cardiovascular Imaging) by a single imaging core lab expert (M. Fukui) in a dedicated cardiac CT core laboratory.

IN VITRO THV ASSESSMENT

To evaluate the extent of S3/S3 Ultra prosthesis deformation compared to the design intent, we performed CT examinations of uncrimped, fully expanded 23, 26 and 29 mm SAPIEN 3 prostheses with the same scanner used in this study cohort. The external stent frame areas for the SAPIEN THV prostheses were measured at 2 levels: leaflet outflow (i.e., the 3 commissural tabs of the

prosthetic leaflets) and leaflet inflow (i.e., the nadir of the prosthetic leaflets) (Figure 1) at nominal expansion. The nominal *in vitro* neosinus volume of each valve size was measured as the volume above the prosthesis leaflets within the prosthesis frame using a short stack of images reconstructed from the nadir to the tips of the THV, using a 1.0 mm slice thickness without a gap, following previously published protocols¹¹.

IN VIVO ASSESSMENT FOR THV DEFORMATION AFTER MVIV IMPLANTATION

Analyses of THV configuration after implantation were performed in a single phase between mid- to end-systole (20 to 40% of the R-R interval), where the best image quality without motion artefact was obtained.

1. Oversizing: oversizing was calculated based on the size of THV and the True ID of SHV as:
Oversizing (mm)=Size of THV–True ID of SHV
2. Implant depth: implant depth was measured as the average of the maximum and minimum distances between the inflow

aspect of THV and SHV, i.e., atrial projection of the THV above the sewing ring of the SHV (Figure 2).

3. Canting: canting was measured as the angle between the axes of the THV and SHV (Figure 2).
4. Deformation of the THV: the analysis of THV deformation comprised the degree of stent frame expansion, eccentricity, neosinus volume index and asymmetric leaflet expansion; this was performed as described in detail previously¹¹ and as summarised in **Supplementary Table 1**.
 - 4.1. THV configuration: The design intent of the S3/S3 Ultra prosthesis is a cylindrical shape with the same cross-sectional area throughout the frame. However, due to the features of the MVIV procedure, i.e., oversizing and desirable flaring at the outflow to prevent embolisation, a cylindrical shape may not be achievable. Hence, the THV configuration after implantation was studied comparing the *in vitro* and *in vivo* configurations.
 - 4.2. THV deformation index: The prosthesis deformation index was calculated as follows:




Figure 1. Assessment of transcatheter heart valve deformation. A,B) Transcatheter heart valve (THV) expansion area and eccentricity are assessed at 3 levels: leaflet inflow, the narrowest level and leaflet outflow. The THV deformation index is calculated using the expansion area at the 3 levels as $(\text{area at leaflet outflow} + \text{area at leaflet inflow}) / (2 \times \text{area at the narrowest level})$. The neosinus volume is measured as the volume above the prosthesis leaflet within the stent frame. C) The design intent of the SAPIEN prosthesis is a cylindrical shape with the same cross-sectional area throughout the frame (A, left). However, 2 distinct configurations of the S3/S3 Ultra frame were observed after MVIV implantation; an “hourglass” shape, where the 2 ends of the THV expand more than the mid-section (i.e., the narrowest level exists between the leaflet inflow and leaflet outflow and usually corresponds to the neo-annulus/sewing ring of the SHV) and a “cone” shape, where the expansion area gradually increases from the inflow to the outflow. MVIV: mitral valve-in-valve




Figure 2. Implant depth, canting and vertical deformation. *A)* To calculate the implant depth of the transcatheter heart valve (THV), the maximum and minimum depths from the inflow aspect of the surgical heart valve (SHV) to that of the THV are measured and averaged (red arrows). Canting is measured as the angle between the axes of the THV and the SHV. For vertical deformation of the THV, the maximum and minimum heights of the THV (orange dotted lines) were measured and we calculated the difference between the THV heights. *B)* The relative location of the minimum height of the transcatheter heart valve to the MV annulus was identified using multiplanar reconstruction. The figure is from the left ventricle view. Ao: aortic valve; LAA: left atrium appendage; max: maximum; min: minimum; MV: mitral valve; RA: right atrium.

THV deformation index=(area at leaflet outflow+area at leaflet inflow)/(2×area at the narrowest level)

When the THV expands to be perfectly cylindrical, this index would be 1.00, whereas an index of >1 represents varying degrees of THV deformation.

- 4.3. Vertical deformation: The minimum and maximum heights of the THV were measured, and vertical deformation was calculated as the difference between the two measurements in mm (**Figure 2**). The relative location of the minimum height of the THV to the mitral valve annulus was identified using multiplanar reconstruction.

HALT ASSESSMENT

HALT is visually identified increased thickness of the THV leaflet on the *in vivo* CT according to Valve Academic Research Consortium (VARC)-3 criteria. Using multiplanar reformats, each prosthetic leaflet was evaluated for the presence of HALT in at least 2 different diastolic phases, and the severity of HALT was semiquantitatively assessed on long-axis views, carefully aligned with the leaflet centre regarding involvement along the curvilinear leaflet beginning at the base, using a 5-tier grading scale: none, ≤25%, 26%-50%, 51%-75%, and >75%¹⁷.

ECHOCARDIOGRAPHIC ASSESSMENT

A non-invasive assessment of THV haemodynamics was performed using transthoracic echocardiography 30 days after MViv with standard assessments using several systems (Philips iE33, EpiQ, CX50; Philips Medical Systems). The mean gradient (MG) across the THV was measured with the modified Bernoulli equation.

STATISTICAL ANALYSIS

The electronic medical records were reviewed for collection of clinical characteristics, comorbidities, procedural characteristics,

and clinical outcomes. Categorical and continuous variables were expressed as frequencies (percentages) and medians (interquartile ranges [IQR]), respectively, and compared using the Kruskal-Wallis test. Due to the small cohort size, variables associated with THV deformation or haemodynamic status after the procedure were assessed using univariable linear regression analysis but not by multivariable model analysis. Correlations were examined with a calculation of Pearson coefficients to assess the relationship between the degree of oversizing and the THV deformation index. P-values<0.05 were considered *a priori* statistically significant. All analyses were performed using SPSS version 25 (IBM).

Results

IN VITRO MEASUREMENTS OF SAPIEN 3

The median external stent frame areas for SAPIEN 3 sizes 23, 26 and 29 mm were 401 (IQR 400-402) mm², 521 (IQR 520-522) mm² and 650 (IQR 650-650) mm², respectively. The neosinus volumes for valve sizes 23, 26 and 29 mm were 3.71 ml, 5.97 ml and 8.00 ml, respectively.

STUDY COHORT

Overall, 38 patients with SHV degeneration underwent MViv implantation with a balloon-expandable SAPIEN 3 or 3 Ultra device during the study period. Those without a postprocedural CT examination (n=5) were excluded, leaving a final study cohort of 33 patients. The post-procedure CT was performed at a median of 35 (IQR 31-47) days after MViv implantation, and adequate image quality for THV assessment was present in all 33 patients. Patients were predominantly elderly (median age: 81 [IQR 76-86] years), 69.7% were on anticoagulation at baseline and had high surgical risk as assessed by the Heart Team (Society of Thoracic Surgeons Predicted Risk of Mortality score: 7.1% [IQR 3.6-11.9%]) (**Table 1**). Mitral regurgitation was the predominant pathology

Table 1. Clinical and procedural characteristics (N=33).

Baseline clinical characteristics	
Age, years	81 [76-86]
Male	14 (42.4)
Duration between SMVR and mitral ViV, years	10 [7-12]
Diabetes mellitus	6 (18.2)
Hypertension	24 (72.7)
Atrial fibrillation/flutter	22 (66.7)
Coronary artery disease	13 (39.4)
eGFR, ml/min/1.73 m ²	48 [42-69]
Anticoagulant therapy	23 (69.7)
STS-PROM score, %	7.1 [3.6-11.9]
Baseline echocardiographic characteristics	
Reason for SHV dysfunction	
Mitral regurgitation	18 (54.5)
Mitral stenosis	8 (24.2)
Mixed mitral stenosis and regurgitation	7 (21.2)
LV ejection fraction, %	60 [56-65]
LV stroke volume index, ml/m ²	26 [22-36]
MVG, mmHg	10 [7-13]
Prosthesis valve characteristics	
SHV leaflet/stent frame type	
Porcine/polymer	8 (24.2)
Pericardial/metal+polymer	25 (75.8)
SHV True ID	
20.5/21/22 mm	4 (12.1)
23/24/25 mm	11 (33.3)
26/27/28/28.5 mm	18 (54.5)
THV size	
#23	1 (3.0)
#26	10 (30.3)
#29	22 (66.7)
Valve-in-valve procedural characteristics	
Transseptal access to mitral valve	30 (90.9)
Balloon-volume filling	
Nominal filling	30 (90.9)
Overfilling	3 (9.1)
Post-balloon dilatation	3 (9.1)
Medications at discharge	
Antiplatelet therapy	29 (87.9)
Anticoagulant therapy (warfarin)	32 (97.0)
At 1-month follow-up	
PT-INR (n=32)	2.2 [1.8-2.5]
≥2.5	9 (28.1)
2.0-2.4	13 (40.6)
<2.0	10 (31.3)
MVG, mmHg (on echocardiography)	9 [6-11]
Values are median [interquartile range] or n (%). eGFR: estimated glomerular filtration rate; ID: internal diameter; LV: left ventricular; MVG: mitral valve gradient; PT-INR: prothrombin time and international normalised ratio; SHV: surgical heart valve; SMVR: surgical mitral valve replacement; STS-PROM: Society of Thoracic Surgeons Predicted Risk of Mortality; THV: transcatheter heart valve; ViV: valve-in-valve	

of SHV degeneration, present in 18 (54.5%) patients. The nominal balloon inflation volume was used during the procedure in all but three patients in whom the S3 size used was either equal to or 0.5 mm larger than the True ID (i.e., 0 or 0.5 mm oversize). In these cases, an extra volume (1 or 2 cc) was added. Post-implant balloon dilatation was performed during the procedure in 3 patients (9.1%) using a True Balloon (Midwest Interventional Systems).

OVERSIZING, IMPLANT DEPTH AND CANTING

Oversizing THV with respect to the True ID of the SHV prosthesis was common, with a median value of 2.0 (1.0-3.5) mm (**Table 2**). The oversizing used was ≤1, 2, 2.5-3, and ≥4 mm in 10, 8, 7 and 8 patients, respectively. In 7 patients with borderline True ID of SHV (i.e., True ID of SHV=22 or 25 mm), a larger size of THV (i.e., 4 mm oversize) was chosen in 6 patients and a smaller one (i.e., 1 mm oversize) in 1 patient. An extreme oversize of 5 mm was used in 2 patients, only during the early experience, and was most likely due to use of the stent ID rather than True ID for sizing of the THV (**Supplementary Table 2**). The median implant depth (i.e., larger % of the THV towards the left atrium) was 6.2 (3.8-7.8) mm. Minimal or no canting (i.e., <10°) occurred in 30 (90.9%) patients, while 3 had canting >10°.

Table 2. One-month postprocedural cardiac computed tomography variables (N=33).

LV ejection fraction, %	48 [41-52]
LV stroke volume index, ml/m ²	36 [30-41]
Relationship between THV and SHV	
Oversizing, mm	2.0 [1.0-3.5]
Implant depth, mm	6.2 [3.8-7.8]
Canting, degree	4.7 [3.1-8.1]
THV deformation	
Expansion, %	
Leaflet inflow (LA side)	82.3 [75.4-88.8]
Narrowest level	74.0 [68.1-84.1]
Leaflet outflow (LV side)	99.8 [92.4-103.8]
Eccentricity	
Leaflet inflow (LA side)	0.25 [0.17-0.31]
Narrowest level	0.24 [0.17-0.33]
Leaflet outflow (LV side)	0.28 [0.20-0.33]
THV deformation index	1.21 [1.13-1.29]
THV configuration	
Hourglass shape	30 (90.9)
Cone shape	3 (9.1)
Vertical deformation, mm	1.7 [0.9-2.4]
Neosinus volume index	0.82 [0.78-0.88]
Asymmetric leaflet expansion, degree	12 [8-19]
Values are median [interquartile range] or n (%). LA: left atrium; LV: left ventricular; SHV: surgical heart valve; THV: transcatheter heart valve	

THV DEFORMATION

Relative to the nominal area of the prosthesis *in vitro*, underexpansion of the implanted THV was observed at all levels of the S3/S3 Ultra frame following the MViv procedure (Table 2), but circularity was maintained, with eccentricity ranging from 0.24 to 0.28, across all levels of the THV, suggesting a near-circular underexpansion throughout the THV.

Overall, THV stent frame deformation was common, with a median THV deformation index of 1.21 (IQR 1.13-1.29). Two distinct configurations of the S3/S3 Ultra frame were observed after MViv implantation: an “hourglass” shape, where the 2 ends of the THV expand more than the mid-section (i.e., the narrowest level exists between the leaflet inflow and leaflet outflow and usually corresponds to the neo-annulus/sewing ring of the SHV), and a “cone” shape, where the expansion area gradually increases from the inflow to the outflow. An hourglass configuration was common and was observed in 30 (90.9%) patients, whereas a cone shape was observed in only 3 (9.1%) patients (Figure 1). Vertical deformation, i.e., difference in THV height of ≥ 2 mm between the maximum and minimum, was observed in 11 (33.3%) patients. The minimum THV height more often occurred medially (near the interatrial septum) than laterally (near the left atrial appendage) (Figure 2).

The neosinus volume was smaller than the nominal volume *in vitro* in all but one patient. The range of neosinus volumes observed was 31-50% smaller in 3 patients, 10-30% in 25 and 0-9% in 4, compared to the nominal volume. In one patient where the neosinus volume was slightly greater than nominal (neosinus volume index, 1.03), the oversize was 0.5 mm for a 29 mm SAPIEN 3 prosthesis, and the THV deformation index was 1.09. In the majority of patients, some degree of leaflet asymmetry was observed with $<10^\circ$, 10° to 20° , and 21° to 30° in 9, 18 and 6 patients, respectively.

DETERMINANTS OF THV DEFORMATION

The degree of oversizing and final configuration of the THV within the SHV affected the extent of the THV deformation (Table 3). A greater degree of oversizing was associated with a greater degree of underexpansion at the narrowest level of the THV, with a higher THV deformation index and, consequently, a smaller neosinus volume index (Figure 3, Supplementary Table 3). A shallower implant depth (i.e., larger % of the THV towards the left atrium) was also

related to a greater degree of underexpansion and a smaller neosinus volume index. Canting was associated with worse eccentricity and greater vertical deformation.

HALT AND HAEMODYNAMICS

Overall, the incidence of HALT on the 30-day postprocedural CT was 27.3% (9 of 33 patients), of which all patients except one had HALT on multiple leaflets or >26 -50% HALT on at least one leaflet (Figure 4). Most patients (32 of 33) were on anticoagulation therapy with warfarin except one due to contraindication. However, the prothrombin time and INR (PT-INR) at the time of the CT scan was below the target range (i.e., <2.5) in 23 of 32 patients (Table 1). Only 2 of 9 patients with HALT had a PT-INR within the therapeutic range (i.e., ≥ 2.5). Notably, among patients with a PT-INR of <2.5 , HALT was predominantly observed with a high THV deformation index of ≥ 1.18 , while no HALT was observed with a THV deformation index of <1.18 (Central illustration).

The mean mitral valve gradient (MVG) on the 30-day postprocedural echocardiography was 9 (6-11) mmHg. This tended to be associated with the type and size of the SHV, oversizing and the degree of THV expansion at the leaflet inflow but was not associated with the presence of HALT (Table 4).

Of the 9 patients with HALT, 4 patients had a follow-up CT 3 months later, 3 of whom continued to demonstrate HALT despite being within the therapeutic range of anticoagulation, while one had resolution of HALT. The remaining 5 patients did not have follow-up CTs due to social reasons. All 9 remain on lifelong anticoagulation.

Discussion

The present investigation is the first study to systematically evaluate the extent and predictors of THV deformation following MViv with a balloon-expandable SAPIEN device. The key findings are as follows: 1) THV deformation, i.e., underexpansion and an hourglass shape, is common after implantation; 2) such THV deformation is negatively affected by an excessive oversizing and a shallower implant depth (i.e., larger % of the THV towards the left atrium), and 3) there is a higher incidence of HALT in patients with a higher degree of deformation and an INR <2.5 . Taken together, these findings might have implications in optimising MViv procedures and achieving better mid- and long-term clinical outcomes.

Table 3. Determinants of THV deformation.

	Expansion at the narrowest area		Eccentricity at the narrowest area		THV deformation index		Vertical deformation		Neosinus volume index		Asymmetric leaflet expansion	
	β	<i>p</i> -value	β	<i>p</i> -value	β	<i>p</i> -value	β	<i>p</i> -value	β	<i>p</i> -value	β	<i>p</i> -value
Oversizing, mm	-0.634	<0.001	-0.167	0.35	0.594	<0.001	0.071	0.69	-0.588	<0.001	0.158	0.38
Implant depth, mm	-0.413	0.02	-0.002	0.99	0.255	0.15	0.289	0.10	-0.425	0.01	0.241	0.18
Canting, degree	0.208	0.25	0.480	0.005	-0.370	0.03	0.411	0.02	-0.109	0.55	-0.104	0.57
\geq Moderate MS due to SHV degeneration	0.286	0.11	0.224	0.21	-0.430	0.01	-0.097	0.59	-0.135	0.45	-0.094	0.60

Univariable linear regression analysis was performed for each THV deformation. MS: mitral stenosis; SHV: surgical heart valve; THV: transcatheter heart valve




Figure 3. Oversizing and transcatheter heart valve deformation. As oversizing increases, the expansion area at the narrowest level of the stent frame decreases (A), THV deformation index increases (B) and neosinus volume index decreases (C). The horizontal lines within the boxes correspond to the median. The bottom of the boxes indicates the 25th percentile (quartile 1), and the top of the boxes represents the 75th percentile (quartile 3). The T-bars that extend from the boxes represent scores outside the middle 50% and the lowest and highest scores. The Xs in the boxes correspond to the mean, and the dots outside of the boxes and T-bars represent outliers. THV: transcatheter heart valve





Figure 4. Transcatheter heart valve deformation after MViV implantation. Example cases of 29 mm SAPIEN 3 (S3) in Epic #31 SHV (True ID=27 mm) with and without hypoattenuating leaflet thickening (HALT) on 30-day postprocedural computed tomography. Yellow lines correspond to the narrowest level of THV, and blue lines correspond to the leaflet inflow and leaflet outflow levels. ID: internal diameter; MViV: mitral valve-in-valve; PT-INR prothrombin time and international normalised ratio; SHV: surgical heart valve; THV: transcatheter heart valve

CENTRAL ILLUSTRATION Stent frame deformation of transcatheter heart valves (THV), i.e., underexpansion and an hourglass shape, is common and is influenced by the degree of oversizing.



Among patients with a prothrombin time international normalised ratio (PT-INR) of <2.5 , hypoattenuating leaflet thickening (HALT) was predominantly observed when the THV deformation index was high, i.e., ≥ 1.18 , while no HALT was observed when the THV deformation index <1.18 .

Prior reports have shown the safety and effectiveness of MViv therapy for patients with SHV degeneration¹⁻⁴. The number of patients who need this therapy is expected to grow with the increasing use of bioprosthetic surgical valves and approval of MViv as an indication for high-risk patients¹⁸. Nevertheless, the long-term durability of THV with MViv remains uncertain. Recent findings on TAVR, ViV and *in vitro* models have shown that the deformation of THV prostheses might increase the risk of THV dysfunction with the potential consequences of pinwheeling, leaflet shear stress and increased flow stagnation leading to HALT^{7-11,19-23}. In addition, there is higher risk of valve thrombosis in the mitral position, as evident from surgical and transcatheter mitral valve replacement experience²⁴. Hence, minimising THV deformation after an MViv procedure may be important for the long-term durability of THV.

Compared to the results of our previous study of TAVR in native aortic stenosis¹¹, we observed that the SAPIEN prosthesis after MViv implantation was more underexpanded (94% vs 82% at the leaflet inflow; 89% vs 74% at the waist level) and more deformed (THV deformation index: 1.09 vs 1.21). In this study,

excessive oversizing and a shallower implant depth (i.e., larger % of the THV towards the left atrium) were factors associated with a greater degree of underexpansion and THV deformation index.

As embolisation is a serious complication after MViv procedures, oversizing of at least 2 mm and the adoption of a flared or conical shape for THV after implantation are recommended due to early experience of embolisation with the SAPIEN XT⁶. However, this contributes to deformation, especially in borderline sizes. Newer-generation devices like the SAPIEN 3/3 Ultra have a different stent structure, with a higher stent height and an outer skirt, which offers an opportunity to reduce THV deformation by minimising oversizing. As SAPIEN devices are only available in 3 mm increments, the considerations above become relevant in patients with borderline SHV sizes where the current available sizes of the SAPIEN device could result in either 1 or 4 mm oversizing (**Supplementary Table 4**). Bioprosthetic valve fracture (BVF) is utilised increasingly during ViV procedures and can be also considered to optimise expansion during MViv procedures^{25,26}. When excessive oversizing is unavoidable and/or significant THV

Table 4. Determinants of postprocedural MVG at 30 days.

	Standardised β	p-value
Age, years	0.168	0.35
Male	-0.294	0.10
Body mass index, kg/m ²	0.032	0.88
≥moderate MS due to SHV degeneration	0.138	0.45
Porcine/polymer SHV (vs pericardial/polymer+metal)	0.321	0.07
SHV True ID ≤22 mm	0.330	0.06
THV #29 (vs #23 or #26)	-0.211	0.24
Post-balloon dilatation	0.051	0.78
Postprocedural CT-derived LV ejection fraction, %	0.149	0.41
Postprocedural CT-derived LV stroke volume index, ml/m ²	-0.148	0.41
Oversizing, mm	0.310	0.08
Implant depth, mm	0.040	0.83
Canting, degree	-0.146	0.42
Expansion at leaflet inflow (LA side), %	-0.328	0.06
Expansion at narrowest level, %	-0.199	0.27
Expansion at leaflet outflow (LV side), %	0.053	0.77
Eccentricity at narrowest level	-0.232	0.19
THV deformation index	0.121	0.50
Difference of THV height, mm	-0.081	0.65
Neosinus volume index	-0.177	0.33
Asymmetric leaflet expansion, degree	-0.376	0.03
HALT	-0.098	0.59

CT: computed tomography; HALT: hypoattenuating leaflet thickening; ID: internal diameter; LA: left atrium; LV: left ventricular; MS: mitral stenosis; MVG: mitral valve gradient; SHV: surgical heart valve; THV: transcatheter heart valve

deformation is observed on the operating table, BVF could be considered to optimise THV expansion after weighing the procedural risks, as it could offer superior THV expansion with lower residual transvalvular gradients^{26,27}. Vertical deformation was observed in 33.3% of patients and partly reflects the issues with the transseptal approach and orientation of the SHV in the mitral position as well as left ventricular cavity shape. Techniques to optimise trajectory can minimise this issue and also allow the operator to implant the valve in a shallow position, as higher implantations in the atrium were associated with more deformation. Further studies are needed.

Although HALT was observed in 27.3% of patients (9 of 33) on the 30-day postprocedural CT in this investigation, 24 of 33 patients were on insufficient anticoagulation therapy (i.e., 1 patient with no anticoagulation therapy and 23 with PT-INR <2.5). However, 2 of 9 patients with therapeutic PT-INR demonstrated HALT. This highlights the importance of adequate

anticoagulation and routine follow-up, as the risk of HALT may be higher in the mitral position, with patient-related risk factors such as atrial fibrillation and left ventricular dysfunction²⁴.

Based on the results of prior studies⁸⁻¹¹ and the present study, excessive THV deformation might be associated with an increased risk of HALT and, subsequently, concern about THV durability, but further studies in larger number of patients with sufficient anticoagulation therapy will be needed to understand factors contributing to the occurrence of HALT after MViv implantation. THV deformation and final configuration may be more important in the mitral position, as four of our patients with similar MViv combinations, i.e., 31 Epic SHV (Abbott) and 29 SAPIEN 3 (oversize of 2 mm), demonstrated varying degrees of HALT but also different final configurations (**Figure 4**).

Similar to what is seen with TAVR in native aortic stenosis, neither the THV deformation index nor the presence of HALT with MViv procedures were associated with an increase in MVG, suggesting it is subclinical and often unable to be detected by transthoracic echocardiography. This explains the lower incidence of clinical thrombosis (2.4%) after MViv therapy reported previously when assessed by transthoracic echocardiography¹⁴. Based on these findings, cardiac CT after MViv implantation might be considered in the short-term follow-up of these patients, as persistent HALT, if untreated, could affect THV degeneration via a process of leaflet fibrosis and then calcification^{26,28}.

Limitations

First, although this is the first study to comprehensively assess postprocedural CT, which was prospectively performed 30 days after MViv implantation, the number of patients is small, and the analysis is limited especially for factors associated with HALT. Second, despite our institutional anticoagulation protocol following MViv procedures, more than half of the patients were on insufficient anticoagulation therapy at the time of the CT scan, highlighting the real-world issues. Third, we did not assess pinwheeling of THV prosthesis leaflets; its assessment by CT remains challenging and can only be performed when excellent image quality is obtained throughout the entire diastolic phase. Lastly, the long-term effects of untreated HALT cannot be discerned from our study. As our protocol does not have longitudinal CT follow-up more than 30 days after the procedure if no HALT observed, our data cannot be used to evaluate later occurrences of HALT in these patients.

Conclusions

THV deformation, i.e., underexpansion and an hourglass shape, is common after MViv implantation with the currently available SAPIEN 3 sizes, and its key determinant is excessive oversizing. The relationship between deformation and HALT needs further study, although those with a higher deformation index and lower INR appeared to be at most risk. Optimisation of THV expansion and a higher and consistent INR target, i.e., >2.5, may be required to prevent HALT following MViv procedures.

Impact on daily practice

The deformation of transcatheter heart valves (THV), i.e., underexpansion and an hourglass shape, is common after mitral valve-in-valve implantation with a SAPIEN 3 prosthesis, and its key determinant is excessive oversizing and a shallower implant depth (i.e., larger % of the THV towards the left atrium). Since excessive THV deformation can cause THV dysfunction, avoiding excessive oversizing and a shallower implant depth and using post-dilatation should be considered to optimise THV expansion.

Conflict of interest statement

J.L. Cavalcante has received consulting fees from 4C Medical, Abbott Structural, Anteris, Aria CV, Boston Scientific, Edwards Lifesciences, HighLife, Medtronic, VDYne, W.L. Gore & Associates, and Xylocor; and has received research grant support from Abbott Northwestern Hospital Foundation. P. Sorajja has received consulting fees from 4C Medical, Abbott Structural, Anteris, Boston Scientific, Edwards Lifesciences, Evolution Medical, Foldax, GLG, Medtronic, Philips, Siemens, Shifamed, W.L. Gore & Associates, VDYne, and xDot Medical; and has received institutional research grant support from Abbott Structural, Medtronic, and Boston Scientific. M. Enriquez-Sarano has received consulting fees from CryoLife, Edwards Lifesciences, HighLife, and ChemImage. V.N. Bapat has received consulting fees from Abbott Structural, Medtronic, Boston Scientific, Edwards Lifesciences, 4C Medical, and Anteris. The other authors have no conflicts of interest to declare.

References

- Urena M, Vahanian A, Brochet E, Ducrocq G, Iung B, Himbert D. Current Indications for Transcatheter Mitral Valve Replacement Using Transcatheter Aortic Valves: Valve-in-Valve, Valve-in-Ring, and Valve-in-Mitral Annulus Calcification. *Circulation*. 2021;143:178-96.
- Bapat VN. Mitral Valve-in-Valve, Valve-in-Ring, and Valve-in-Mitral Annular Calcification: Are We There Yet? *Circulation*. 2021;143:117-9.
- Simonato M, Whisenant B, Ribeiro HB, Webb JG, Kornowski R, Guerrero M, Wijeyesundera H, Sondergaard L, De Backer O, Villablanca P, Rihal C, Eleid M, Kempfert J, Unbehaun A, Erlebach M, Casselman F, Adam M, Montorfano M, Ancona M, Saia F, Ubben T, Meincke F, Napodano M, Codner P, Schofer J, Pelletier M, Cheung A, Shuvy M, Palma JH, Gaia DF, Duncan A, Hildick-Smith D, Veulemans V, Sinning JM, Arbel Y, Testa L, de Weger A, Eltchaninoff H, Hemery T, Landes U, Tchetché D, Dumonteil N, Rodés-Cabau J, Kim WK, Spargias K, Kourkovieli P, Ben-Yehuda O, Teles RC, Barbanti M, Fiorina C, Thukkani A, Mackensen GB, Jones N, Presbítero P, Petronio AS, Allali A, Champagnac D, Bleiziffer S, Rudolph T, Iadanza A, Salizzoni S, Agrifoglio M, Nombela-Franco L, Bonaros N, Kass M, Bruschi G, Amabile N, Chhatriwalla A, Messina A, Hirji SA, Andreas M, Welsh R, Schoels W, Hellig F, Windecker S, Stortecky S, Maisano F, Stone GW, Dvir D. Transcatheter Mitral Valve Replacement After Surgical Repair or Replacement. *Circulation*. 2021;143:104-16.
- Guerrero M, Pursnani A, Narang A, Salinger M, Wang DD, Eleid M, Kodali SK, George I, Satler L, Waksman R, Meduri CU, Rajagopal V, Inglessis I, Palacios I, Reisman M, Eng MH, Russell HM, Pershad A, Fang K, Kar S, Makkar R, Saucedo J, Pearson P, Bokhary U, Kaptzan T, Lewis B, Tommaso C, Krause P, Thaden J, Oh J, Lang RM, Hahn RT, Leon MB, O'Neill WW, Feldman T, Rihal C. Prospective Evaluation of Transcatheter TMVR for Failed Surgical Bioprostheses: MITRAL Trial Valve-in-Valve Arm 1-Year Outcomes. *JACC Cardiovasc Interv*. 2021;14:859-72.
- Agarwal V, Kaple RK, Mehta HH, Singh P, Bapat VN. Current state of transcatheter mitral valve implantation in bioprosthetic mitral valve and in mitral ring as a treatment approach for failed mitral prosthesis. *Ann Cardiothorac Surg*. 2021;10:585-604.
- Bapat VV, Khaliel F, Ihleberg L. Delayed migration of Sapien valve following a transcatheter mitral valve-in-valve implantation. *Catheter Cardiovasc Interv*. 2014;83:E150-4.
- Sathananthan J, Sellers S, Barlow AM, Stanová V, Fraser R, Toggweiler S, Allen KB, Chhatriwalla A, Murdoch DJ, Hensey M, Lau K, Alkhourair A, Dvir D, Asgar AW, Cheung A, Blanke P, Ye J, Rieu R, Pibarot P, Wood D, Leipsic J, Webb JG. Valve-in-Valve Transcatheter Aortic Valve Replacement and Bioprosthetic Valve Fracture Comparing Different Transcatheter Heart Valve Designs: An Ex Vivo Bench Study. *JACC Cardiovasc Interv*. 2019;12:65-75.
- Khodae F, Barakat M, Abbasi M, Dvir D, Azadani AN. Incomplete expansion of transcatheter aortic valves is associated with propensity for valve thrombosis. *Interact Cardiovasc Thorac Surg*. 2020;30:39-46.
- Martin C, Sun W. Transcatheter Valve Underexpansion Limits Leaflet Durability: Implications for Valve-in-Valve Procedures. *Ann Biomed Eng*. 2017;45:394-404.
- Midha PA, Raghav V, Sharma R, Condado JF, Okafor IU, Rami T, Kumar G, Thourani VH, Jilaihawi H, Babaliaros V, Makkar RR, Yoganathan AP. The Fluid Mechanics of Transcatheter Heart Valve Leaflet Thrombosis in the Neosinus. *Circulation*. 2017;136:1598-609.
- Fukui M, Bapat VN, Garcia S, Dworak MW, Hashimoto G, Sato H, Gössl M, Enriquez-Sarano M, Lesser JR, Cavalcante JL, Sorajja P. Deformation of Transcatheter Aortic Valve Prostheses: Implications for Hypoattenuating Leaflet Thickening and Clinical Outcomes. *Circulation*. 2022;146:480-93.
- Fukui M, Sorajja P, Cavalcante JL, Thao KR, Okada A, Sato H, Wang C, Koike H, Hamid N, Enriquez-Sarano M, Lesser JR, Bapat VN. Deformation of Transcatheter Heart Valve Following Valve-in-Valve Transcatheter Aortic Valve Replacement: Implications for Hemodynamics. *JACC Cardiovasc Interv*. 2023;16:515-26.
- Guerrero M, Vemulapalli S, Xiang Q, Wang DD, Eleid M, Cabalka AK, Sandhu G, Salinger M, Russell H, Greenbaum A, Kodali S, George I, Dvir D, Whisenant B, Russo MJ, Pershad A, Fang K, Coylewright M, Shah P, Babaliaros V, Khan JM, Tommaso C, Saucedo J, Kar S, Makkar R, Mack M, Holmes D, Leon M, Bapat V, Thourani VH, Rihal C, O'Neill W, Feldman T. Thirty-Day Outcomes of Transcatheter Mitral Valve Replacement for Degenerated Mitral Bioprostheses (Valve-in-Valve), Failed Surgical Rings (Valve-in-Ring), and Native Valve With Severe Mitral Annular Calcification (Valve-in-Mitral Annular Calcification) in the United States: Data From the Society of Thoracic Surgeons/American College of Cardiology/Transcatheter Valve Therapy Registry. *Circ Cardiovasc Interv*. 2020;13:e008425.
- Yoon SH, Whisenant BK, Bleiziffer S, Delgado V, Dhoble A, Schofer N, Eschenbach L, Bansal E, Murdoch DJ, Ancona M, Schmidt T, Yzeiraj E, Vincent F, Niikura H, Kim WK, Asami M, Unbehaun A, Hirji S, Fujita B, Silaschi M, Tang GHL, Kuwata S, Wong SC, Frangieh AH, Barker CM, Davies JE, Lauten A, Deuschl F, Nombela-Franco L, Rampat R, Niez PFG, Masson JB, Wijeyesundera HC, Sievert H, Blackman DJ, Gutierrez-Ibanes E, Sugiyama D, Chakravarty T, Hildick-Smith D, de Brito FS Jr, Jensen C, Jung C, Smalling RW, Arnold M, Redwood S, Kasel AM, Maisano F, Treede H, Ensminger SM, Kar S, Kaneko T, Pilgrim T, Sorajja P, Van Belle E, Prendergast BD, Bapat V, Modine T, Schofer J, Frerker C, Kempfert J, Attizzani GF, Latib A, Schaefer U, Webb JG, Bax JJ, Makkar RR. Outcomes of transcatheter mitral valve replacement for degenerated bioprostheses, failed annuloplasty rings, and mitral annular calcification. *Eur Heart J*. 2019;40:441-51.
- Tang L, Lesser JR, Schneider LM, Burns MR, Gössl M, Garberich R, Niikura H, Witt D, Sorajja P. Prospective Evaluation for Hypoattenuated Leaflet Thickening Following Transcatheter Aortic Valve Implantation. *Am J Cardiol*. 2019;123:658-66.
- Harloff MT, Chowdhury M, Hirji SA, Percy ED, Yazdchi F, Shim H, Malarczyk AA, Sobieszczek PS, Sabe AA, Shah PB, Kaneko T. A step-by-step guide to transseptal valve-in-valve transcatheter mitral valve replacement. *Ann Cardiothorac Surg*. 2021;10:113-21.
- VARC-3 WRITING COMMITTEE; Gèneux P, Piazza N, Alu MC, Nazif T, Hahn RT, Pibarot P, Bax JJ, Leipsic JA, Blanke P, Blackstone EH, Finn MT, Kapadia S, Linke A, Mack MJ, Makkar R, Mehran R, Popma JJ, Reardon M, Rodes-Cabau Josep, Van Mieghem NM, Webb JG, Cohen DJ, Leon MB. Valve Academic Research Consortium 3: Updated Endpoint Definitions for Aortic Valve Clinical Research. *J Am Coll Cardiol*. 2021;77:2717-46.
- Goldstone AB, Chiu P, Baiocchi M, Lingala B, Patrick WL, Fischbein MP, Woo YJ. Mechanical or Biologic Prostheses for Aortic-Valve and Mitral-Valve Replacement. *N Engl J Med*. 2017;377:1847-57.
- Gunning PS, Saikrishnan N, McNamara LM, Yoganathan AP. An in vitro evaluation of the impact of eccentric deployment on transcatheter aortic valve hemodynamics. *Ann Biomed Eng*. 2014;42:1195-206.
- Vahidkhah K, Barakat M, Abbasi M, Javani S, Azadani PN, Tandar A, Dvir D, Azadani AN. Valve thrombosis following transcatheter aortic valve replacement: significance of blood stasis on the leaflets. *Eur J Cardiothorac Surg*. 2017;51:927-35.
- Lim WY, Lloyd G, Bhattacharyya S. Mechanical and surgical bioprosthetic valve thrombosis. *Heart*. 2017;103:1934-41.

22. Fuchs A, De Backer O, Brooks M, de Knecht MC, Bieliauskas G, Yamamoto M, Yanagisawa R, Hayashida K, Søndergaard L, Kofoed KF. Subclinical leaflet thickening and stent frame geometry in self-expanding transcatheter heart valves. *EuroIntervention*. 2017;13:e1067-75.
23. Hein M, Schoechlin S, Schulz U, Minners J, Breitbart P, Lehane C, Neumann FJ, Ruile P. Long-Term Follow-Up of Hypoattenuated Leaflet Thickening After Transcatheter Aortic Valve Replacement. *JACC Cardiovasc Interv*. 2022;15:1113-22.
24. Pagnesi M, Moroni F, Beneduce A, Giannini F, Colombo A, Weisz G, Latib A. Thrombotic Risk and Antithrombotic Strategies After Transcatheter Mitral Valve Replacement. *JACC Cardiovasc Interv*. 2019;12:2388-401.
25. Meier D, Payne GW, Mostaço-Guidolin LB, Bouchareb R, Rich C, Lai A, Chatfield AG, Akodad M, Salcudean H, Lutter G, Puehler T, Pibarot P, Allen KB, Chhatriwalla AK, Søndergaard L, Wood DA, Webb JG, Leipsic JA, Sathananthan J, Sellers SL. Timing of bioprosthetic valve fracture in transcatheter valve-in-valve intervention: impact on valve durability and leaflet integrity. *EuroIntervention*. 2023;18:1165-77.
26. Carlidge TRG, Doris MK, Sellers SL, Pawade TA, White AC, Pessotto R, Kwiecinski J, Fletcher A, Alcaide C, Lucatelli C, Densem C, Rudd JHF, van Beek EJR, Tavares A, Virmani R, Berman D, Leipsic JA, Newby DE, Dweck MR. Detection and Prediction of Bioprosthetic Aortic Valve Degeneration. *J Am Coll Cardiol*. 2019;73:1107-19.
27. Sathananthan J, Fraser R, Hatoum H, Barlow AM, Stanová V, Allen KB, Chhatriwalla AK, Rieu R, Pibarot P, Dasi LP, Søndergaard L, Wood DA, Webb JG. A bench test study of bioprosthetic valve fracture performed before versus after transcatheter valve-in-valve intervention. *EuroIntervention*. 2020;15:1409-16.
28. Sellers SL, Turner CT, Sathananthan J, Carlidge TRG, Sin F, Bouchareb R, Mooney J, Nørgaard BL, Bax JJ, Bernatchez PN, Dweck MR, Granville DJ, Newby DE, Lauck S, Webb JG, Payne GW, Pibarot P, Blanke P, Seidman MA, Leipsic JA. Transcatheter Aortic Heart Valves: Histological Analysis Providing Insight to Leaflet Thickening and Structural Valve Degeneration. *JACC Cardiovasc Imaging*. 2019;12:135-45.

Supplementary data

Supplementary Appendix 1. Methods: CT imaging acquisition.

Supplementary Table 1. Methodology for analysis of THV deformation.

Supplementary Table 2. Prosthesis valve characteristics.

Supplementary Table 3. THV deformation according to oversizing.

Supplementary Table 4. Size selection of SAPIEN and oversizing.

The supplementary data are published online at:

<https://eurointervention.pronline.com/>

[doi/10.4244/EIJ-D-23-00614](https://doi.org/10.4244/EIJ-D-23-00614)



Supplementary data

Supplementary Appendix 1. Methods: CT imaging acquisition.

Acquisition was performed with dual-source scanners (Siemens Healthineers, Germany) using either second-generation (Siemens SOMATOM Definition Flash 2×128×0.6 mm; temporal resolution = 75 msec, 2015 to 2016), or a third-generation (Siemens SOMATOM Definition Force 2×192×0.6 mm; temporal resolution = 66 msec, 2016 to 2022) scanners. Tube voltage ranged from 80 to 120 kV and tube current was adjusted according to Care kV and an in-house algorithm. A retrospective ECG-gating acquisition was used to cover the entire cardiac cycle, with dose modulation maximizing tube current during systolic phase. Intravenous contrast (Omnipaque 350, GE Healthcare Inc., Marlborough, MA) doses ranged from 60 to 120 ml at injection rates from 4 to 8 ml/sec depending upon the patient's renal function, body mass index, and kV used followed by 60 ml of normal saline flush at 4.5 ml/sec. Imaging acquisition was triggered using bolus tracking technique with region of interest at the ascending aorta after reaching a pre-defined threshold of 150 Hounsfield units followed breathing commands (averaging 3 seconds scan delay) and scan acquisition. Images were reconstructed at 1.5 mm with isotropic voxels throughout the entire cardiac cycle [minimal of 10 phases, maximal of 20 phases (5% of R-R interval for sinus rhythm, or in patients with atrial fibrillation, absolute milliseconds reconstruction at 50 msec intervals to the shortest R-R interval used during scan acquisition)].

Supplementary Table 1. Methodology for analysis of THV deformation.

Parameters of THV deformation	Methodology
<p>1. Stent frame expansion</p>	<p>$Expansion (\%) = (measured\ THV\ area / nominal\ THV\ area) \times 100$</p> <p>The external stent frame area was measured at 3 levels of SAPIEN prosthesis as described in Figure 1 (i.e., leaflet inflow, leaflet outflow and the narrowest level) and used to calculate expansion (%) at each level.</p>
<p>2. Eccentricity</p>	<p>$Eccentricity = \sqrt{1 - \frac{(minor\ diameter\ of\ THV\ area)^2}{(major\ diameter\ of\ THV\ area)^2}}$</p> <p>The major and minor diameters of the external stent frame were measured at the 3 levels of SAPIEN prosthesis and used to calculate eccentricity at each level.</p> <p>Eccentricity varies between 0 and 1, with 0 representing a perfect circle, and values approaching 1 representing greater degrees of ovality (Figure 1).</p>
<p>3. Neosinus volume index</p>	<p>$Neosinus\ volume\ index = THV\ neosinus\ volume\ in\ vivo / THV\ neosinus\ volume\ in\ vitro$</p> <p>The <i>THV neosinus volume</i> was measured as the volume above the prosthesis leaflets within the prosthesis frame and then standardized by nominal <i>in vitro</i> volume.</p> <p>A neosinus volume index of < 1.00 represents smaller volume relative to the nominal <i>in vitro</i> volume.</p>
<p>4. Asymmetric leaflet expansion</p>	<p>The expansion of each prosthesis leaflet was evaluated by measuring the angle formed by the border stent struts assigned to each prosthetic leaflet, and the THV center point at the coaptation level. Given that full leaflet expansion would be assumed to be 120 degrees (ideal), <i>asymmetric leaflet expansion</i> was calculated as the sum of the difference between each measured leaflet angle and 120 degrees.</p>

Supplementary Table 2. Prosthesis valve characteristics.

	SHV leaflet/stent frame	SHV	SHV label size	SHV True ID	THV type	THV size	Oversizing	SHV degeneration			Post-dilatation	Balloon type/size	HALT+	
								MR	MS	MSR				
1	Porcine / Polymer	Mosaic (n=2)	27	22.0	S3Ultra	26	+4mm	0	0	1	1	BLLN TRUE (24*4.5)	0	
2			29	24.0	S3	29	+5mm	1	0	0	0		1	
3		HancockII (n=2)	25	20.5	S3Ultra	23	+2.5mm	1	0	0	0		0	
4			27	22.0	S3	26	+4mm	0	0	1	0	0		
5		Biocor/Epic (n=21)		25	21.0	S3	26	+5mm	0	1	0	1	BLLN TRUE (24*4.5)	0
6				27	23.0	S3	26	+3mm	1	0	0	0		0
7				27	23.0	S3	26	+3mm	0	0	1	0	0	
8				27	23.0	S3Ultra	26	+3mm	1	0	0	0	0	
9				27	23.0	S3Ultra	26	+3mm	1	0	0	0	0	
10				27	23.0	S3Ultra	26	+3mm	0	0	1	0	0	
11				29	25.0	S3	29	+4mm	1	0	0	0	1	
12				29	25.0	S3	29	+4mm	1	0	0	0	1	
13				29	25.0	S3	29	+4mm	0	0	1	0	0	
14				31	27.0	S3	29	+2mm	1	0	0	0	1	
15				31	27.0	S3	29	+2mm	1	0	0	0	1	
16				31	27.0	S3	29	+2mm	1	0	0	0	1	
17				31	27.0	S3	29	+2mm	1	0	0	0	1	
18				31	27.0	S3	29	+2mm	1	0	0	0	0	
19				31	27.0	S3	29	+2mm	1	0	0	0	0	
20				31	27.0	S3	29	+2mm	1	0	0	0	0	
21				31	27.0	S3	29	+2mm	0	0	1	0	0	
22				33	28.5	S3	29	+0.5mm	1	0	0	0	1	
23				33	28.5	S3	29	+0.5mm	1	0	0	0	0	
24				33	28.5	S3	29	+0.5mm	1	0	0	0	0	
25				38	28.5	S3	29	+0.5mm	1	0	0	0	0	
26	Pericardial /Polymer+ Metal	Perimount (n=3)	27	25.0	S3	26	+1mm	0	1	0	1	BLLN TRUE (26*4.5)	0	
27			27	25.0	S3	29	+4mm	0	0	1	0		0	
28			33	28.5	S3	29	+0.5mm	0	1	0	0		0	
29		Magna (n=5)		27	26.0	S3Ultra	26	+0mm	0	1	0	0	0	
30				27	26.0	S3	29	+3mm	0	1	0	0	1	
31				29	28.0	S3	29	+1mm	0	1	0	0	0	
32				29	28.0	S3	29	+1mm	0	1	0	0	0	
33				33	28.5	S3	29	+0.5mm	0	1	0	0	0	

HALT, hypo-attenuating leaflet thickening; ID, internal diameter; MR, mitral regurgitation; MS, mitral stenosis; MSR, mitral stenosis regurgitation; SHV, surgical heart valve; THV, transcatheter heart valve. “n” represents number of patients.

Supplementary Table 3. THV deformation according to oversizing.

	Oversizing				<i>p</i>
	0-1 mm n=10	2 mm n=8	2.5- 3 mm n=7	4-5mm n=8	
Relationship between THV and SHV					
Implant depth (mm)	5.3 [2.2-7.8]	5.9 [3.4-7.9]	6.0 [5.3-7.1]	7.1 [2.4-8.5]	0.92
Canting (degree)	6.8 [3.0-9.9]	3.7 [3.0-7.4]	6.1 [4.5-8.9]	3.7 [0.4-6.2]	0.28
THV deformation					
Expansion (%)					
Leaflet inflow (LA side)	89.5 [86.8-93.8]	78.5 [73.5-83.4]	73.8 [72.4-79.0]	82.9 [75.8-83.9]	0.001
Narrowest level	87.5 [80.5-89.1]	72.9 [70.7-75.7]	67.6 [60.7-69.3]	68.0 [62.0-79.5]	<0.001
Leaflet outflow (LV side)	100.2 [98.4-105.7]	103.8 [97.9-108.1]	92.3 [92.1-99.0]	96.7 [91.1-101.5]	<0.05
Eccentricity					
Leaflet inflow (LA side)	0.28 [0.15-0.34]	0.26 [0.22-0.31]	0.17 [0.12-0.27]	0.25 [0.21-0.32]	0.61
Narrowest level	0.27 [0.20-0.33]	0.17 [0.15-0.32]	0.25 [0.21-0.41]	0.20 [0.11-0.31]	0.14
Leaflet outflow (LV side)	0.25 [0.19-0.31]	0.26 [0.11-0.34]	0.25 [0.18-0.30]	0.33 [0.30-0.39]	0.09
THV deformation index	1.09 [1.07-1.14]	1.25 [1.23-1.28]	1.28 [1.19-1.32]	1.27 [1.16-1.41]	0.001
THV hourglass shape – no. (%)	9 (90.0)	6 (75.0)	7 (100.0)	8 (100.0)	0.28
Vertical deformation (mm)	1.9 [1.2-3.3]	0.9 [0.5-1.7]	1.7 [1.1-2.00]	1.9 [1.1-3.9]	0.30
Neosinus volume index	0.89 [0.87-0.95]	0.84 [0.81-0.86]	0.74 [0.70-0.78]	0.80 [0.77-0.81]	<0.001
Asymmetric leaflet expansion (degree)	10 [8-18]	10 [5-18]	14 [10-26]	16 [8-20]	0.44
30-day outcome					
MVG (mmHg) on echocardiography	7 [6-8]	10 [8-13]	9 [5-11]	10 [7-11]	0.03
MVG ≥ 10mmHg – no. (%)	0 (0)	4 (50.0)	3 (42.9)	5 (62.5)	<0.04
HALT– no. (%)	1 (10.0)	4 (50.0)	1 (14.3)	3 (37.5)	0.22

HALT, hypo-attenuating leaflet thickening; LA, left atrium; LV left ventricle; MVG, transmitral valve gradient; SHV, surgical heart valve; THV, transcatheter heart valve.

Supplementary Table 4. Size selection of SAPIEN and oversizing.

SHV type and size		SHV True ID (mm)	SAPIEN size (mm)	Oversizing
Borderline size selection				
Biocor/Epic	#29	25	26	1mm
			29	4mm
Perimount	#27	25	26	1mm
			29	4mm
CE SAV Porcine	#25	22.5	23	0.5mm
			26	3.5mm
			29	4mm
CE Standard Porcine	#29	25	26	1mm
			29	4mm
Hancock II	#27	22	23	1mm
			26	4mm
Mosaic	#27	22	23	1mm
			26	4mm

SHV, surgical heart valve; ID, internal diameter.