

Sizing for self-expanding transcatheter aortic valve implantation



Hasan Jilaihawi, MD; Zhen-Gang Zhao, MD; Mathew Williams, MD

NYU Langone Health, New York, NY, USA

Computed tomography (CT) is the standard of care for device sizing for balloon-expandable (BE) and self-expanding (SE) transcatheter aortic valve implantation (TAVI) and has dramatically impacted outcomes in this field, particularly the frequency of paravalvular aortic regurgitation (PAR)^{1,2}. In this issue of EuroIntervention, Kim et al describe multicentric outcomes of 198 patients treated with a novel SE transcatheter heart valve (THV) device, the Edwards CENTERA (Edwards Lifesciences, Irvine, CA, USA) in the context of CT-based oversizing, as assessed prospectively by a dedicated CT core lab³.

Article, see page 511

The Edwards CENTERA is a relatively novel THV whose first-in-man (FIM) experience was described in 2010, but which has only recently completed a CE mark study⁴ and has yet to commence a US feasibility or pivotal study. The THV consists of a relatively short (in comparison to existing SE designs), contour-shaped, SE nitinol valve frame with bovine pericardial tissue leaflets. It is inserted by transfemoral access using an expandable 14 Fr introducer sheath and features a steerable delivery catheter to traverse the aortic arch and for coaxial alignment within the annulus, a loading capsule containing the pre-attached THV, and a tapered tip to facilitate valve crossing; a 6-V battery-powered

motorised handle enables valve loading and deployment, as well as repositioning with the possibility of recapturing up to 85% of the deployed valve⁴.

Cardiac CT annular perimeter sizing was performed at the discretion of the site with baseline guidance of optimal oversizing in the 10% to 20% range. Importantly, predilatation was performed in all cases and post-dilatation in one third. Favourable outcomes were achieved overall with this approach to sizing with low rates of new permanent pacemaker implantation (PPI) and PAR and low transvalvular gradients with high effective orifice areas (EOA). Each of these important outcomes are related to THV sizing but are multifactorial and also related to several other factors that would have been difficult to comprehensively address in the present manuscript.

PERMANENT PACEMAKER

The rate of PPI was only 4.7% (nine cases), 5.5% for pacemaker-naïve patients, extremely low for an SE THV. This number of cases makes it challenging to reliably investigate its association with the extent of oversizing or other potential predictors. The potential variation in pacemaker indication across participating centres may also be relevant. According to the original report

*Corresponding author: NYU Langone Health, 550 First Avenue, New York, NY 10016, USA.

E-mail: hasanjilaihawi@gmail.com

which described the 30-day results of the Safety and Performance Study of the Edwards CENTERA-EU study⁴, third-degree atrio-ventricular block (AVB) occurred in 14 patients, but only nine of them received new PPI. Moreover, the frequency of PPI must be interpreted in the context of pre-existing conduction disease in the form of right bundle branch block (RBBB), data that were not available in this focused CT-sizing paper.

PARAVALVULAR AORTIC REGURGITATION

PAR \geq mild was seen in 37.9% (with no cases of severe PAR and only 0.5% of PAR \geq moderate). This compares very favourably to both BE and contemporary SE TAVI data. However, in an additional three patients (1.5%), the THV was not successfully placed due to valve embolisation into the left ventricle in one patient and valve migration in one patient, with conversion to SAVR in both cases, as well as poor coaxiality in one case with subsequent implantation of a SAPIEN 3 THV.

THE RELEVANCE OF DEPTH OF IMPLANTATION

Although sizing for the BE SAPIEN 3 is relatively simple due to its cylindrical structure, self-expanding aortic THVs are contoured and non-cylindrical; this means that the interaction between device and annulus by degree of oversizing has a greater dependence on device positioning, data which are notably omitted from this paper and the original CENTERA CE mark study manuscript^{3,4}. Depth of implant is also a critical determinant of several outcomes after both BE and SE TAVI, including PPI and PAR^{5,6}. There are also historical data for the CoreValve[®] (Medtronic, Minneapolis,

MN, USA) with higher depth of implantation improving effective orifice area (EOA) and prosthesis-patient mismatch⁷.

The shorter design of the CENTERA frame combined with an aggressively demarcated waist are unique features for an SE THV (**Table 1**) that may mandate and help facilitate a more standardised depth of implantation but means that leaflets are intra-annular. In contrast, the longer frame of the Evolut[™] R/PRO (Medtronic), although permitting supra-annular leaflets and a greater degree of freedom in depth of implantation, may also allow a greater degree of variability in this parameter. Perhaps for this reason, a wide range of outcomes has been reported. Although instructions for use (IFU) for the Evolut R/PRO specify an implant depth of 3-5 mm, many operators, ourselves included, strive to routinely implant much higher at 1-3 mm. Interestingly, although SE THV data have generally reported much higher PPI rates than the CENTERA data, the ACURATE *neo*[™] device (Boston Scientific, Marlborough, MA, USA) that incorporates technology to standardise depth of implantation (**Table 1**) has a low reported rate of new PPI of 8%⁸. Both device technology and operator technique demand focus on implanting SE THV as high (aortic) as possible in the aortic annulus to mitigate PPI.

OVERSIZING, CENTERA AND OTHER SE TAVI DEVICES

The authors place considerable emphasis on higher rates of PPI and the greater influence of oversizing in SE TAVI devices other than the CENTERA. In this context, we examined data from our single-centre real-world experience of 175 consecutive patients undergoing TAVI between April 2017 and March 2018 with the latest

Table 1. Self-expanding TAVI devices: unique features pertaining to sizing in comparison to Evolut R/PRO.

	Edwards CENTERA	Medtronic Evolut PRO	Medtronic Evolut R (34 XL)	Abbott Portico	Boston Scientific ACURATE <i>neo</i>	Venus Medtech Venibri	MicroPort VitaFlow
Unique features compared with Evolut R/PRO	Contoured frame designed to anchor and seal in the annulus. Dry-tissue leaflets. Pre-crimped. Motorised deployment	NA	NA	Large cell geometry	Large cell geometry. Top-down deployment. Supra-annular anchoring mechanism	Increased radial strength. Dry-tissue leaflets. Pre-crimped	Large cell geometry. Motorised deployment
Retrievable after near complete deployment	Yes	Yes	Yes	Yes	No	No	No
Sizes (mm)	23, 26, 29	23, 26, 29	23, 26, 29 (34)	23, 25, 27, 29	23, 25, 27	23, 26, 29, 32	21, 24, 27, 30
Intended annular treatment range (diameter)	18-26 mm	18-26 mm	18-26 (26-30) mm	19-27 mm	21-27 mm	17-29 mm	17-29 mm
Frame length	18-23 mm	45 mm	45 (46) mm	49-53 mm	48-51 mm	41-51 mm	48-50 mm
Frame material	Nitinol	Nitinol	Nitinol	Nitinol	Nitinol	Nitinol	Nitinol
Covered external skirt/wrap	No	Porcine pericardial wrap	No	No	Porcine pericardial outer skirt	No	Polyethylene terephthalate (PET) skirt
Leaflet pericardium	Bovine	Porcine	Porcine	Bovine	Porcine	Porcine	Bovine
Leaflet position	Intra-annular	Supra-annular	Supra-annular	Intra-annular	Supra-annular	Supra-annular	Supra-annular
Sheath	14 Fr eSheath	16 Fr	14 Fr (16 Fr)	18/19 Fr	18 Fr	19 Fr	16/18 Fr

iteration Medtronic device (Evolut PRO) which treats a similar range of aortic annular dimensions to the CENTERA (Figure 1). CT data were measured prospectively by a dedicated CT core lab; pre-discharge echo parameters were reported clinically according to VARC-2 criteria by four dedicated valve centre expert echocardiographers. With the Evolut PRO, we had a mean implant depth of 3.0 ± 1.9 mm, PAR \geq mild of 29.9%, PAR \geq moderate of 0.6% and an EOA of 2.35 ± 0.45 cm² (Figure 1); excluding patients with prior PPI, new PPI was 5.7% (6.6% for pacemaker-naïve patients) at 30 days. Predilatation was performed in 27.0% and post-dilatation in 43.1%. In the setting of a systematically high depth of implantation and generally oversizing >10% (96% of cases) by perimeter (mean oversizing 17.6 ± 5.1 %), we saw a slightly greater influence of degree of oversizing and outcomes with our Evolut PRO experience than the CENTERA-EU data, but similarly low rates of PPI and favourable haemodynamics on echocardiography (Figure 1).

CALCIUM AND OVERSIZING

Aside from device sizing and depth of implantation, left ventricular outflow tract (LVOT) and aortic valve calcium have been identified as independent predictors of PAR and PPI in numerous studies of both BE and SE TAVI^{6,9,10}, but not in the present

study. Although data were not presented, it is conceivable that the extremes of aortic valve and subannular calcium severity may have been less frequently included in the setting of a clinical trial, where stricter anatomical selection criteria are generally adopted in comparison to a real-world experience. It is possible that the correlation of oversizing ratio and the incidence of PAR with the CENTERA could be potentially modified by aortic valve calcium severity and this merits further study in a real-world experience.

Conclusions

Although the data for the CENTERA device suggest a less dramatic interaction between degree of oversizing and THV outcomes than data previously published for other SE devices, this must be viewed in the context of depth of implantation and a systematic degree of oversizing $\geq 10\%$. Taking into account these factors, the favourable outcomes achieved with this device may also be achieved with other SE THV devices but may be more operator-dependent.

Conflict of interest statement

H. Jilaihawi is a consultant to Edwards Lifesciences and Venus Medtech and receives research funding from Edwards Lifesciences,




Figure 1. Comparison of oversizing in NYU Langone data with Evolut PRO (grey) with CENTERA-EU (red). Comparison of endpoints is shown for oversizing categories according to perimeter: A) Paravalvular regurgitation (PAR). B) Rate of permanent pacemaker implantation (PPI). C) Mean effective orifice area (EOA). To place % oversizing in context, the manufacturer-recommended Evolut PRO sizing table is shown with degree of oversizing for each valve size (D).

Medtronic and Abbott Vascular. M. Williams is a consultant to Medtronic and Abbott Vascular and receives research funding from LivaNova, Medtronic, Bracco, Edwards Lifesciences, NeoChord, Abbott Vascular. Z. Zhao has no conflicts of interest to declare.

References

1. Jilaihawi H, Kashif M, Fontana G, Furugen A, Shiota T, Friede G, Makhija R, Doctor N, Leon MB, Makkar RR. Cross-sectional computed tomographic assessment improves accuracy of aortic annular sizing for transcatheter aortic valve replacement and reduces the incidence of paravalvular aortic regurgitation. *J Am Coll Cardiol*. 2012;59:1275-86.
2. Popma JJ, Gleason TG, Yakubov SJ, Harrison JK, Forrest JK, Maini B, Ruiz CE, Pinto DS, Costa M, Resar J, Conte J, Crestanello J, Chang Y, Oh JK, Reardon MJ, Adams DH. Relationship of Annular Sizing Using Multidetector Computed Tomographic Imaging and Clinical Outcomes After Self-Expanding CoreValve Transcatheter Aortic Valve Replacement. *Circ Cardiovasc Interv*. 2016;9:e003282.
3. Kim U, Blanke P, Windecker S, Kasel AM, Schäfer U, Walters D, Linke A, Le Breton H, Schymik G, Spence MS, Sondergaard L, Abdel-Wahab M, Worthley S, Tchetché D, Reichenspurner H, Ohana M, Sellers SL, Leipsic JA. Computed tomography-based oversizing and incidence of paravalvular aortic regurgitation and permanent pacemaker implantation with a new-generation self-expanding transcatheter heart valve. *EuroIntervention*. 2018;14:e511-8.
4. Reichenspurner H, Schaefer A, Schäfer U, Tchetché D, Linke A, Spence MS, Sondergaard L, LeBreton H, Schymik G, Abdel-Wahab M, Leipsic J, Walters DL, Worthley S, Kasel M, Windecker S. Self-Expanding Transcatheter Aortic Valve System for Symptomatic High-Risk Patients With Severe Aortic Stenosis. *J Am Coll Cardiol*. 2017;70:3127-36.
5. Binder RK, Webb JG, Toggweiler S, Freeman M, Barbanti M, Willson AB, Alhassan D, Hague CJ, Wood DA, Leipsic J. Impact of post-implant SAPIEN XT geometry and position on conduction disturbances, hemodynamic performance, and paravalvular regurgitation. *JACC Cardiovasc Interv*. 2013;6:462-8.
6. Mauri V, Deuschl F, Frohn T, Schofer N, Linder M, Kuhn E, Schaefer A, Rudolph V, Madershahian N, Conradi L, Rudolph TK, Schäfer U. Predictors of paravalvular regurgitation and permanent pacemaker implantation after TAVR with a next-generation self-expanding device. *Clin Res Cardiol*. 2018 Apr 17. [Epub ahead of print].
7. Jilaihawi H, Chin D, Spyt T, Jeilan M, Vasa-Nicotera M, Bence J, Logtens E, Kovac J. Prosthesis-patient mismatch after transcatheter aortic valve implantation with the Medtronic-Corevalve bioprosthesis. *Eur Heart J*. 2010;31:857-64.
8. Möllmann H, Hengstenberg C, Hilker M, Kerber S, Schäfer U, Rudolph T, Linke A, Franz N, Kuntze T, Nef H, Kappert U, Walther T, Zembala MO, Toggweiler S, Kim WK. Real-world experience using the ACURATE neo prosthesis: 30-day outcomes of 1,000 patients enrolled in the SAVI TF registry. *EuroIntervention*. 2018;13:e1764-70.
9. Jilaihawi H, Makkar RR, Kashif M, Okuyama K, Chakravarty T, Shiota T, Friede G, Nakamura M, Doctor N, Rafique A, Shibayama K, Mihara H, Trento A, Cheng W, Friedman J, Berman D, Fontana GP. A revised methodology for aortic-valvar complex calcium quantification for transcatheter aortic valve implantation. *Eur Heart J Cardiovasc Imaging*. 2014;15:1324-32.
10. Maeno Y, Abramowitz Y, Kawamori H, Kazuno Y, Kubo S, Takahashi N, Mangat G, Okuyama K, Kashif M, Chakravarty T, Nakamura M, Cheng W, Friedman J, Berman D, Makkar RR, Jilaihawi H. A Highly Predictive Risk Model for Pacemaker Implantation After TAVR. *JACC Cardiovasc Imaging*. 2017;10:1139-47.