

The PICCOLETO study and beyond

Bernardo Cortese*, MD, FESC

Interventional Cardiology, Cliniche Humanitas Gavazzeni, Bergamo, Italy

KEYWORDS

- coronary angioplasty
- PCB
- DES
- restenosis

Abstract

Aims: After the excellent results of the PACCOCATH ISR and PEPCAD II trials in in-stent restenosis, paclitaxel-coated balloon (PCB) still has to prove its efficacy in native coronary disease.

Methods and results: In the PICCOLETO randomised trial, patients with stable or unstable angina undergoing PCI of small coronary vessels (≤ 2.75 mm) were randomised to the Dior PCB (28 patients) or the Taxus DES (29 patients). The primary study endpoint was percent diameter stenosis at 6-month angiographic follow-up (non-inferiority); secondary endpoints were angiographic binary restenosis and major adverse cardiac events (MACE: death, Q-wave myocardial infarction, TLR) at 9-month follow-up (non-inferiority). The two groups were not dissimilar for clinical and angiographic data. The study was interrupted after enrolment of two-thirds of the patients due to a clear superiority of one study group. The primary endpoint was not met: the PCB group had higher percent diameter stenosis (43.6% vs. 24.3%, $p=0.029$), angiographic restenosis (32.1 vs. 10.3%, $p=0.043$), and higher occurrence of MACE (35.7% vs. 13.8%, $p=0.054$).

Conclusions: Dior PCB failed to show equivalence to Taxus DES regarding angiographic endpoints during PCI of small coronary arteries. We hypothesise that it concerned the lack of efficacy of the device used (which has since been replaced by its second generation) rather than a class-effect in native coronary vessels.

* Corresponding author: *Interventional Cardiology, Cliniche Humanitas Gavazzeni, Via Mauro Gavazzeni, Bergamo, Italy*
E-mail: bcortese@gmail.com

Introduction

Since the results of studies from the group coordinated by B. Scheller et al became available, there has been great enthusiasm around the new technology of paclitaxel-coated balloons (PCB). The first study involved a PCB constructed with the PACCOATH technology,¹ that outperformed plain balloon angioplasty in in-stent restenosis, allowing for lower late lumen loss (LLL), 0.09 ± 0.49 vs. 0.76 ± 0.86 mm ($p=0.003$). Shortly afterwards, M. Unverdorben and colleagues published the results of the study Pepecad II², where a newer PCB (SeQuent Please; B. Braun, Melsungen, Germany) that utilised an evolution of the PACCOATH technology was compared with the Taxus Libertè drug-eluting stent (DES) (Boston Scientific, Natick, MA, USA) in a bare metal stent (BMS) restenosis population. One hundred and thirty-one patients were enrolled in this study, with a primary endpoint of in-segment LLL at 6-months angiography. The PCB was seen to outperform the DES: LLL rate was respectively 0.17 ± 0.42 mm (median 0.09 mm, interquartile range 0.42 mm) and 0.38 ± 0.61 mm (median 0.24 mm, interquartile range 0.55 mm), $p=0.03$. Moreover, at 12-months, the rate of MACE was 9% and 22%, respectively ($p=0.08$), mainly driven by a higher target vessel revascularisation (TVR) in the DES group (15% vs. 6% in the PCB group, $p=0.15$). A new revolution had begun. It was shortly after this that some researchers decided to test the efficacy and safety of this technology in native coronary artery disease.

While waiting for the results of the PEPCAD III study which had not yet been published, a few months earlier we published the results of the PICCOLETO trial³, a mechanistic study comparing the Dior I generation PCB (Eurocor GmbH, Bonn, Germany) with the Taxus Libertè stent (Boston Scientific, Natick, MA, USA) in small coronary vessels with diameter < 2.75 mm.

Methods

This single-centre randomised clinical trial aimed at enrolling 80 patients with complete 9-month follow-up. Patients were eligible if they had a clinical indication to PCI in a small vessel and stable or unstable angina at presentation. Myocardial infarction, unstable haemodynamics and serum creatinine > 2 mg/dl were the exclusion criteria.

The primary study endpoint was percent diameter stenosis at 6-month QCA (non-inferiority hypothesised). Secondary endpoints were binary angiographic restenosis (defined as restenosis of $> 50\%$) and MACE (death, myocardial infarction, target lesion revascularisation) at clinical follow-up (non-inferiority). In the case of an unsatisfying result in the PCB group (coronary dissection, TIMI flow < 3 and residual stenosis $> 50\%$), a bare metal stent had to be implanted. Predilatation was mandatory only in the DES group, but if gentle balloon deployment in the PCB group was unsuccessful, predilatation with an uncoated balloon was recommended. PCB inflation time was 45 seconds repeated twice, following the manufacturer's recommendations.

PCI was performed following standard protocols with unfractionated heparin or bivalirudin. Dual antiplatelet treatment with

aspirin and clopidogrel was mandatory for 12 months in the DES group, and one month in the PCB group (three months in case of BMS implantation).

The study protocol contemplated a 6-month angiographic follow-up, and a 9-month clinical follow-up.

Results

Pre-specified *ad interim* analysis performed after enrolment and analysis of two-thirds of the patients showed the superiority of one study group, therefore we decided to stop enrolment after having complete data on 28 patients in the PCB group and 29 in the DES group.

Table 1 shows the main characteristics of the patient population. Briefly, baseline and angiographic data showed no significant differences between the two study groups. An exception concerned the number of stents implanted (stents implanted per patient: 0.36 in the PCB group and 1.17 in the DES group, $p=0.024$) and lesion predilatation (25% and 86.2%, $p=0.001$). Angiographic and procedural success did not change across study groups.

The Taxus group was found superior to the PCB group regarding the primary efficacy endpoint ($24.3\pm 25.1\%$ vs. $43.6\pm 27.4\%$, $p=0.029$) (**Figure 1**). The Taxus stent outperformed the Dior balloon also concerning the secondary endpoints such as binary restenosis (10.3% vs. 32.1%, $p=0.043$) and MACE rates (13.8% vs. 35.7%, $p=0.054$), the latter mainly due to a higher clinically-driven TLR rate (**Table 2**). No safety endpoints were raised in this study, as myocardial infarction and death rates were similar between the two study groups.

Discussion

This was the first comparison of a PCB and a DES in a randomised clinical trial of *de novo* lesions of small coronary vessels, and the PCB failed to demonstrate its equivalence with the Taxus stent.

The analysis of our results astonished us, raising doubts about the overall efficacy of PCB in native vessels. However, due to the good results of previous trials (even in different lesions) and the clear advantages of the PCB technology over a stent strategy (**Table 3**),




Figure 1. Primary study endpoint. The Taxus group was found superior to the PCB group regarding percent diameter stenosis. PCB: paclitaxel-coated balloon

Table 1. Baseline clinical and angiographic characteristics.

	Dior I balloon first generation (n=28)	Taxus stent (n=29)	p value
Male	22 (78.6%)	22 (75.9%)	0.81
Age \pm SD (years)	68 \pm 9	67 \pm 10	0.67
Hypertension	21 (75.0%)	20 (70.8%)	0.61
Diabetes mellitus	13 (37.9%)	11 (46.4%)	0.51
Dyslipidaemia	17 (60.7%)	13 (54.2%)	0.23
Previous revascularisation with CABG	3 (10.7%)	4 (13.8%)	0.72
Previous revascularisation with PCI	3 (10.7%)	4 (13.8%)	0.72
Previous myocardial infarction	5 (17.9%)	6 (20.7%)	0.78
Stable angina	13 (46.4%)	13 (44.8%)	0.90
Type B2/C lesions*	17 (60.7%)	20 (69.0%)	0.51
Bifurcation lesion	6 (21.4%)	7 (24.1%)	0.81
Lesion length \pm SD (mm)	12.41 \pm 5.89	11.38 \pm 7.12	0.56
Reference vessel diameter \pm SD (mm)	2.45 \pm 0.28	2.36 \pm 0.25	0.20
MLD \pm SD (mm)	0.48 \pm 0.33	0.40 \pm 0.30	0.30
Lesion stenosis \pm SD (% LD)	86.0 \pm 12.1	89.14 \pm 10.6	0.31
Stents implanted per lesion	0.36	1.17	0.024
Stent or balloon diameter \pm SD (mm)	2.48 \pm 0.30	2.54 \pm 0.10	n.a.
Stent or balloon length \pm SD (mm)	18.6 \pm 5.24	18.9 \pm 7.63	n.a.
Lesion predilatation (uncoated balloon)	7 (25.0%)	25 (86.2%)	0.001
Maximal inflation pressure \pm SD (atmospheres)	7.71 \pm 2.2	13.41 \pm 1.38	0.02
Final MLD \pm SD (mm)	2.47 \pm 0.22	2.63 \pm 0.23	0.009
Final lesion stenosis \pm SD (% LD)	19.0 \pm 17.3	9.9 \pm 9.2	0.016
Angiographic success	96.4%	100%	0.30
Procedural success	96.4%	100%	0.30

LAD: left anterior descending artery; LCX: left circumflex artery; LD: lumen diameter; MLD: minimal lumen diameter; PCB: paclitaxel-coated balloon; PCI: percutaneous coronary intervention; RCA: right coronary artery; CABG: coronary artery bypass intervention * according to ACC/AHA classification, SG Ellis et al, *Circulation* 1990;82:1193-1202

Table 2. Angiographic and clinical results.

	Dior I balloon first generation (n=28)	Taxus stent (n=29)	p value
Reference vessel diameter, mm \pm SD	2.54 \pm 0.47	2.58 \pm 0.24	0.73
Minimal lumen diameter, mm \pm SD	1.11 \pm 0.65	1.94 \pm 0.72	0.0002
Angiographic binary restenosis, n (%)	9 (32.1)	3 (10.3)	0.043
Death, %	3.6	3.5	0.98
Cardiac death, %	0	0	0.97
Myocardial infarction, %	3.6	0	0.30
Target lesion revascularisation, %	32.1	10.3	0.15
Stent thrombosis/abrupt vessel closure (def/prob), %	0	0	0.97
MACE, %	35.7	13.8	0.054

MACE: major adverse clinical events

these results require further explanation. Particularly, the release kinetics of paclitaxel in PCB differs from that of currently available DES that elute this drug in weeks or months. In comparison, the PCB releases the drug almost entirely in the first 48 hours; however, its biological effect has been discovered to persist for the first 14 days (**Figure 2**). This is very important, because the process of

restenosis is seen in the first days after the baro-trauma induced by angioplasty, and paclitaxel primarily exerts its effect at that time. Moreover, this technology, that allows us to reserve stent implantation only in case of suboptimal immediate results, does not require the implantation of a permanent prosthesis and the use of a polymer, both potentially thrombogenic.

Table 3. Advantages of paclitaxel-coated balloon over drug-eluting stent strategies.

No polymer
No permanent prosthesis
No strut-related drug diffusion
Immediate, strong release
No drug found after the first days
Biological effect up to 14 days




Figure 2. Drug elution time (red dotted line) and duration of biological effect (yellow dotted line) by PCB, as well as main comparators among DES. PCB: paclitaxel-coated balloon

Therefore, it is possible that the true limit of this study is the device used itself, the Dior I generation PCB. It appears probable that this point has also been taken into account by the manufacturer as, in fact, a second generation of the Dior balloon has been recently released.

This second generation Dior balloon still elutes paclitaxel at the same dosage (3 mg/mm²) and has been prepared with a similar process including the charging of the balloon with the drug using a micropipetting process with the balloon folded. However, profound differences exist among the two devices. The first generation Dior balloon had a nanoporous surface containing microcrystals of pure paclitaxel that are then embedded on the vessel wall at the time of balloon inflation. The second generation Dior II balloon contains shellac as a paclitaxel carrier. Shellac is an inert substance that has already been approved by the FDA as a food additive. It is mostly composed of aleuritic acid, jalaric acid and shelloic acid, and has

a molecular weight of about 1000 u. The microporous balloon surface contains a 1:1 mixture of paclitaxel and shellac. Finally, this new technology allows complete paclitaxel deposition with a unique, shorter inflation (around 30 seconds) instead of the two 45-60 seconds inflations previously suggested.

The Dior II generation balloon has been recently tested in an animal model, where it showed a 20-fold higher drug concentration after 45 minutes when compared to its predecessor;⁴ moreover, tissue paclitaxel resulted in much lower concentration 12 hours after balloon inflation, a result comparable with other PCB in previous studies.

There are several ongoing studies that are testing this new device. Among them, we are about to start the enrolment for the Piccoletto II study, a trial comparing the Taxus Liberté DES with this new generation PCB for the treatment of small coronary vessels.

In conclusion, the PICCOLETO trial failed to demonstrate its equivalence to a DES for the treatment of small vessel disease, both in terms of angiographic and clinical restenosis. Further studies, thanks to newer technological advances, are needed to test the efficacy and safety of PCB in native coronary artery disease. We expect better results, both clinical and angiographic, from the next trial in the native coronary field.

Conflict of interest statement

B. Cortese is a consultant for Eurocor.

References

- Scheller B, Hehrlein C, Bocksch W, Rutsch W, Haghi D, Dietz U, Bohm M, Speck U. Treatment of coronary in-stent restenosis with a paclitaxel-coated balloon catheter. *N Engl J Med.* 2006;355:2113-2124.
- Unverdorben M, Vallbracht C, Cremers B, Heuer H, Hengstenberg C, Maikowski C, Werner GS, Antoni D, Kleber FX, Bocksch W, Leschke M, Ackermann H, Boxberger M, Speck U, Degenhardt R, Scheller B. Paclitaxel-Coated Balloon Catheter Versus Paclitaxel-Coated Stent for the Treatment of Coronary In-Stent Restenosis. *Circulation.* 2009;119:2986-94.
- Cortese B, Micheli A, Picchi A, Coppolaro A, Bandinelli L, Severi S, Limbruno U. Paclitaxel-coated balloon versus drug-eluting stent during PCI of small coronary vessels, a prospective randomised clinical trial. The PICCOLETO Study. *Heart.* 2010;96:1291-1296.
- Posa A, Nyolczas N, Hemetsberger R, Pavo N, Petnehazy O, Petrasi Z, Sangiorgi G, Gyongyosi M. Optimization of drug-eluting balloon use for safety and efficacy: Evaluation of the 2nd generation paclitaxel-eluting DIOR-balloon in porcine coronary arteries. *Catheter Cardiovasc Interv* 2010;76:395-403.