

The TandemHeart[®], percutaneous transseptal left ventricular assist device: a safeguard in high-risk percutaneous coronary interventions. The six-year Rotterdam experience

Pascal Vranckx^{1,3}, MD; Emanuele Meliga^{2,4}, MD; Peter PT De Jaegere², MD PhD; Martin Van den Ent², MD PhD; Evelyn S Regar², MD PhD; Patrick W Serruys^{2*}, MD PhD FESC FACC

1. Department of Cardiac Intensive Care & Interventional Cardiology, Virga Jesseziekenhuis, Hasselt, Belgium; 2. Department of Interventional Cardiology, Thoraxcenter, Erasmus Medical Centre, Rotterdam, The Netherlands; 3. Cardialysis B.V., Clinical Research Organisation & Core Laboratories, Rotterdam, The Netherlands; 4. Department of Interventional Cardiology, S. Giovanni Battista Hospital, University of Turin, Turin, Italy

None of the authors have a conflict of interest to declare.

KEYWORDS

Heart-assist device,
angioplasty,
TandemHeart[®]

Abstract

Aims: Percutaneous coronary interventions (PCI) in high-risk cardiac patients are preferentially referred to specialised myocardial intervention centres (MIC). Included in this group are patients with a haemodynamic collapse or high likelihood of haemodynamic collapse, either during balloon inflation or with acute vessel closure. The TandemHeart[®], a percutaneous transseptal left ventricular assist (PTVA[®]) that can be introduced using standard catheterisation laboratory techniques, offers interesting perspectives to reduce procedural risks.

Methods and results: Between September 2000 to July 2006, The TandemHeart[®], supported the circulation of 23 patients (age: range 46-74, mean 59) admitted to our centre for high risk, either emergency or elective, PCI. Successful implantation was achieved in 100% of patients. The mean time for implementation of circulatory support was 35 minutes (range 16-62). The index PCI was successful in all patients except two. A pump flow up to 4L/min was achieved with significant reduction of left ventricular filling pressures, pulmonary capillary wedge pressure and with significant increase of systemic arterial pressures. Duration of support ranged from 1-222 hours (mean 31±49.8 hours). Five patients died with the TandemHeart[®] in place, four of whom were in irreversible cardiogenic shock at admission. Mild to moderate access site bleeding was seen in 27% of patients. One patient experienced a loge syndrome of the leg. Core temperature (Ct) decreased to <36.5°C in six patients, profound hypothermia (Ct < 35°C) was observed in two patients. There was no technical device failure.

Conclusions: The TandemHeart[®] - PTVA[®] provides effective, total left ventricular support in very high risk PCI settings. The rate of device related cardiac and vascular complications was acceptable.

* Corresponding author: Dpt. of Interventional Cardiology, Thoraxcenter, Ba583a, Erasmus Medical Centre, Dr. Molewaterplein 40, 3015 GD Rotterdam, The Netherlands

E-mail: p.w.j.c.serruys@erasmusmc.nl

© Europa Edition 2008. All rights reserved.

Introduction

With increasing operator experience, refinement in technology and adjunctive pharmacological treatment, percutaneous coronary intervention (PCI) is now considered the treatment of choice for many high-risk subgroups in which PCI was previously contraindicated. PCI may even be a valuable option for those patients where coronary bypass grafting (CABG) is clinically contraindicated¹⁻³. Many of these procedures are elective.

In general the benefits of a specific (percutaneous) procedure should be weighed against the risks involved; taking into account alternative treatment strategies, the interventional and intensive care team experience⁴ and taking in consideration the individual patient risk scores, that may aid the operator in selecting or avoiding adjunctive pharmacotherapy or specific devices⁵⁻⁹.

Patients in whom it is considered that PCI poses a significantly high risk include these with high likelihood of haemodynamic collapse, either during balloon inflation or with acute vessel closure⁴. In general, these are the patients in whom a large amount of the viable myocardium is supplied directly or indirectly by the affected artery. The development of a percutaneous ventricular assist device (VAD) that could be easily and quickly inserted prophylactically or in case of an haemodynamic collapse proves invaluable in this setting.

The TandemHeart®, a percutaneous transeptal left ventricular assist device (PTVA®) system (CardiacAssist Inc., Pittsburgh, PA, USA), has been demonstrated to effectively reverse cardiogenic shock after myocardial infarction¹¹ and to play an important role for bridging to another definitive therapy¹². This system allows for rapid implementation of circulatory support using standard interventional techniques in the catheterisation laboratory and is designed to deliver up to 4.5 litres of blood flow per minute. In this report, we will report on a single centre six year clinical experience with the TandemHeart® in high risk PCI procedures.

Methods

Patients selection

Since 2000, the Thoraxcenter, Erasmus University Medical Centre, Rotterdam, The Netherlands started a program evaluating percutaneous left ventricular assist devices (LVAD) during high risk PCI. Between September 2000 to July 2006, twenty-three patients admitted to our centre for acute coronary syndromes (ACS)/ ST segment elevation myocardial infarction (STEMI) or for elective PCI, were treated with TandemHeart®.

In non-elective patients, the indication for TandemHeart® support was based on already established haemodynamic instability before the index PCI procedure, defined as typically with low cardiac output (cardiac index < 2.2 L/min/m²), peripheral signs of tissue hypoperfusion (decreased urine output and/or cold extremities), systemic hypotension (systolic blood pressure < 100 mmHg) despite vasopressor therapy and in the presence of appropriate left ventricular filling pressures (pulmonary capillary wedge pressure, PCWP ≥ 15 mmHg). In case of presence of an intra aortic balloon pump (IABP) haemodynamic measurements were done with the IABP paused for 60 seconds.

In patients who underwent elective PCI, circulatory assist was indicated in procedures which had a presumed high risk of ischaemic/haemodynamic complications, based on the presence of severely depressed left ventricular (LV) function and planned treatment of \geq one complex lesions in vessels supplying a large amount of myocardium. CABG was not considered a treatment option in any of these patients. Major exclusion criteria were predominant right ventricular failure requiring right ventricular support, severe aortic insufficiency, severe sepsis and anoxic brain damage.

System description

The TandemHeart® PTVA® (Figure 1-2) incorporates arterial perfusion cannula configurations ranging from 9 to 17 Fr, an unique 21 Fr venous transeptal cannula, and a centrifugal blood pump. Oxygenated blood from the patient's left atrium is supplied to the pump by the trans-septal cannula and then returned to the patient's systemic arterial circulation. The centrifugal pump contains a single moving part (rotor/impeller) that is suspended by a thin lubricating film of fluid to reduce heat and friction, thereby reducing the risk of thrombus formation. The pump connects to a microprocessor-based controller that displays PTVA speed and flow. These parameters are controlled by adjustment of a single knob. The controller also provides automatic system monitoring and alarms indicating conditions that require action. The system is designed to deliver up to 4.5-5 litres (L) of blood flow per minute, depending on the size of the arterial cannulation and the filling conditions of the left atrium, while operating at a relatively low speed (7500 revolutions per minute, RPM).




Figure 1. Schematic of the system deployed. The Tandemheart® removes oxygenated blood from the left atrium via a 21 Fr transeptal cannula that is advanced from the femoral vein through the inter-atrial septum into the left atrium (A). The oxygenated blood is returned into the arterial vascular system via a 15 up to 17 Fr arterial cannula by means of a centrifugal pump (B). The centrifugal pump contains a single moving part (rotor/impeller) that is suspended by magnetic force on a thin lubricating film of fluid to reduce heat and friction, thereby reducing the risk of thrombus formation.



Figure 2. Microprocessor-based controller used in this series (A) displaying PTVA speed and flow. These parameters are controlled by adjustment of a single knob. The controller includes extensive self-diagnostic and alarm features to ensure patient support without the need for constant operator surveillance. Panel (B) shows the recently introduced Escort Controller. This light (21 pound or 9.53 kg) weight controller can be mounted at the patient's bedside and includes an easy-load infusion system. The latter allows the lubrication fluid infusion line to be installed in seconds.

Insertion technique

A standard transseptal puncture technique, using an Inoue guidewire, was used to gain access into the left atrium from the right femoral vein. Transseptal puncture was carried out by an experienced operator only. The inter-atrial septum was dilated with an 2-stage 14/21 Fr dilator and the venous inflow cannula was inserted (proper positioning in the left atrium was checked by either angiography and/or transesophageal echocardiography). One 15 Fr or 17 Fr perfusion catheter (Bio-medicus cannula, Medtronic, Minneapolis, MN, USA) was inserted into the femoral artery and advanced to the common iliac artery. An ileac angiography was performed before the placement of the perfusion catheter to delineate the arterial anatomy and size and to disclose eventual luminal obstructions. In three patients a dual femoral approach with Y connection was used for arterial access. Both arterial and venous cannula were fixed to the skin with multiple sutures in order to secure their position. The pump was placed on the upper leg of the patient. After checking the central venous pressure, priming and cautious de-airing of the entire system, the arterial and venous cannulas sets were connected to the pump by a heparin-coated

tygon tubing (up to 30 cm length) and the pump was activated to its maximum rotation speed. Output pump flow was measured by an external electromagnetic flow meter (HT 311, transonic).

Device weaning and removal

A predefined weaning protocol was initiated at the moment vasoactive drugs were reduced to a stable minimal level (dobutamine up to 4 $\mu\text{g}/\text{kg}/\text{min}$., norepinephrine up to 0.1 4 $\mu\text{g}/\text{kg}/\text{min}$.) provided stable haemodynamic parameters (mixed venous or central venous saturation > 65%). A stepwise reduction of pump assist was performed, by reducing pump speed from 7500 to 3500 RPM (steps of 500-1,000 mL/min), adapted to the medical condition of each patient individually. The final removal decision was based on medical judgement. A weaning period up to six months was respected for most patients. The pump was not stopped until immediately before removal. No specific instructions on access site closure were issued.

Concomitant treatment

The index coronary intervention was performed after TandemHeart® implantation and functioning. Arterial access was obtained via the contra-lateral femoral artery and the interventional strategy was left to the operator, according to standard techniques. Standard intensive care was provided to each patient. Patient sedation and intubation was considered for clinical and/or comfort reasons. A balloon tipped pulmonary artery (PA) catheter was placed via the contra lateral femoral vein for haemodynamic monitoring purposes if clinically indicated. Cardiac output determinations were made using standard thermodilution techniques.

During the index procedure, a constant flow (10 ml/h) of heparinised infusate is maintained providing a localised concentration of heparin in the interior of the pump in order to obtain localised anti-coagulation thereby minimising systemic heparinisation, the risk of bleeding and thrombus formation. Additional boluses of heparin were administered peripherally to maintain the activated coagulation time to approximately 200 seconds during routine support (400 seconds during insertion) or an activated partial thromboplastin time of between 65 and 80 seconds.

Predefined safety endpoints and time periods

The following predefined clinical events related to the use of the TandemHeart® were assessed in each patient: any minor or major TIMI bleeding⁹, laceration of the insertion vessel/limb ischaemia, device related thromboembolic events, infective complications, residual atrial septum defect after device removal, device malfunction or failure.

The insertion time was defined as the time from providing of the Brockenbrough needle to the 'full' heparinisation following connection of the femoral cannula to the pump. The duration of support was defined as the time from 'full' heparinisation at time of connection to the pump till the removal of the femoral cannula.

Statistical analysis

Variables with normal distribution were analysed using parametric tests while variables with a non-normal distribution were analysed with non-parametric tests. Continuous variables are expressed as mean \pm SD or median \pm SD and differences are compared using

Student t test or Mann Whitney test. Categorical variables are expressed as counts and percentages. Differences were assessed by Fisher exact test or chi-square test, as appropriate. All analyses were performed using SPSS version 12 statistical software (SPSS Inc., Chicago, IL, USA). A two-tailed p value < 0.05 was considered significant for hypothesis testing.

Results

Baseline characteristics

Baseline clinical and procedural characteristics are listed in Table 1. Mean age was 59 ± 9.4 years, 19 patients (83%) were men. The ejection fraction was <30% in 16 patients (70%). Fifteen patients (65%) suffered a three vessel disease, six of whom with left main involvement. The indications for TandemHeart® support were elective/high risk PCI in 15 patients and non-elective (emergency) in eight patients with ACS and acute heart failure. Eight patients experienced an evolving (< 36 hours) myocardial infarction, five of which were complicated by acute heart failure/cardiogenic shock (CS). Eight patients needed intubation and mechanical ventilation as part of their medical treatment and 11 were on inotropic support at the start of the implant procedure. Standard EuroScore in our patients ranged between 2 and 14 (mean 6.5) providing a surgical risk of mortality between 1.5% and 46.6% (mean: 11.3±11.6%). Very high risk patients (EuroScore > 9) were nine¹³.

Haemodynamic effects

The TandemHeart® showed good performance, achieving flow rates up to 4.0 L/min. Mean systemic arterial pressure was 74.8±18 mmHg at baseline and 85.6±19 mmHg after pump functioning (p= 0.023).

Table 1. Baseline characteristics; n (%)

Age (mean±SD)		59±9.4
Sex	M	19 (82.6)
	F	4 (17.4)
Previous MI		11 (47.8)
Previous CABG		4 (17.4)
Clinical presentation	STEMI&CS	5 (21.7)
	STEMI	3 (13)
	NSTEMI	1 (4.3)
	UA	5 (21.7)
	Angina/dyspnea	8 (34.7)
	Scheduled PCI	1 (4.3)
CAD extension (23 pts)	1VD	3 (13)
	2VD	5 (21.7)
	3VD	9 (39.1)
	3VD + LM	6 (26)
EF	<20%	11 (47.8)
	20-30%	5 (21.7)
	30-50%	3 (13)
	>50%	4 (17.4)
Indication for PTVA	elective	15 (65.2)
	emergency	8 (34.7)
Intubation		8 (34.7)
Inotropic support		11 (47.8)
EuroScore (mean±SD)		6.5±3.4
	HRP (>6)	3 (13)
	VHRP (>9)	9 (39.1)

Pulmonary wedge pressure was 16.8±5.6 mmHg at baseline and 13.6±6 after pump functioning (P= 0.002). The pulse pressure was reduced on support from 39.7±17.7 to 31.5±17 (p<0.001) (Table 2). In two patients the initial period after implantation of the TandemHeart® was characterised by a complete non-pulsatile arterial blood pressure (pulse pressure < 7 mmHg). In one patient with the recovery of the heart function < 24 hours after implantation, increasing pulse pressures became evident. Modulation of pulsatility with low amplitude to the non-pulsatile blood flow produced by the PTVA could be achieved by varying pump flow. One patient did not show recovery of pulsatile blood flow and he died with the TandemHeart® *in situ*.

Table 2. Hemodynamic effects of PTVA.

	Baseline	PTVA	P value
MSBP (mean±SD)	74.8±18	85.6±19	0,023
Pulse P (mean±SD)	39.7±17.7	31.5±17	0,001
PCWP (mean±SD)	16.8±5.6	13.6±6	0,002

A pulmonary artery catheter was inserted in 20 out of 23 patients. The hemodynamic measurements reported were immediately before Tandem Heart® insertion and at the end of the PCI procedure (before the patient was transferred to the cardiac intensive care department)

Procedural and clinical outcomes

Procedural and clinical outcomes are summarised in Table 3. Index-PCI procedural success was achieved in 42 lesions (96%). A total of 44 lesions were treated (one up to four lesions per patient), all including stent implantation (mean stent per patient rate: 2.3). The mean time for insertion of the TandemHeart® was 34±12 minutes (range 16-62), the procedural success was 100%. The circulatory assist period ranged from 1h to 222h, mean 31±49.8h. One patient with a severe femoral artery stenosis needed PTA and stenting of the iliac artery vessel prior to cannulation.

Seventeen patients were successfully weaned from the TandemHeart®. In-hospital death occurred in six patients. In total five patients died with the device *in situ*, four of which were considered in irreversible CS at admission. Another patient admitted for CS was initially stabilised after 222 hours on TandemHeart®

Table 3.

PCI Results; n (%)		
PCI success rate		96%
Number of treated lesions		44
Stent per patient		2,3
PTVA Insertion success rate		100%
PTVA Insertion Time (mean±SD)		34±12
Assist period, hours (range/mean±SD)		1-222/31±49.8
Death	In hospital	6 (26)
	with PTVA in place	5* (22)
	at 6 months FU	9 (39)
	at 12 months FU	10 (43)
	at 2 years	11 (47)
	at 5 years	11 (47)

* 4/5 patients suffered cardiogenic shock

support. However, this patient developed abrupt circulatory failure after removal of the pump and died a few hours later. No further invasive treatment was possible due to the abrupt worsening of his clinical condition. During the follow-up period, three patients were lost and one patient died due to progressive heart failure < 30 days. An uneventful six months follow-up was observed in 13 patients.

Complications

Potential complications related with the use of TandemHeart® are summarised in Table 4. In this series of patients, no procedural complications during device insertion were noted and no cardiac tamponade, thromboembolic events nor device failure occurred. In two elective cases (8.7%) the weaning needed to be postponed due to the occurrence of profound hypothermia (core temperature 32.7 °C and 34.2 °C) and the need for active re-warming. Minor to mild access site bleeding complications occurred in 26% of

Table 4. Potential complications of PTVA; n(%)

Puncture of aortic root, coronary sinus, posterior free wall	0 (0)
TandemHeart pVAD system failure	0 (0)
Thromboembolism	0 (0)
Neurologic dysfunction	0 (0)
Hemolysis	0 (0)
Cardiac tamponade	0 (0)
Deep venous thrombosis	0 (0)
Arrhythmias	0 (0)
Systemic hypothermia (35-36.5)	4 (17.4)
Profound hypothermia (<35)	2 (8.7)
Bleeding	6 (26)
Cannulation site infection	0
Cannula dislodgment	0
Distal leg ischaemia	1 (4.3)
Total no. of events	13

Table 5. Exclusion criteria for this study were:

- Participation in another trial with an investigational drug/ device during the last 60 days
- Concomitant disease that interferes with the prognosis.
- Contra indications to standard drugs for coronary intervention and coronary heart disease.
- Coagulopathy-chronic anticoagulant therapy / uncontrolled active bleeding.
- CNS damage resulting in fixed dilated pupils not related to pharmacologics.
- Severe aortic/mitral valve stenosis, significant aortic valve insufficiency.
- Rupture of the ventricular wall
- Right heart failure, defined as the need of a right ventricular assist device.
- Severe peripheral vascular disease (e.g. aorto-iliac occlusive disease).
- Stroke of haemorrhagic or unknown origin
- Non-haemorrhagic stroke, that took place within the past 30 days
- A transient ischaemic attack, that took place within the past 30 days.

patients (n=6), none required specific treatment. Lower limb ischaemia complicated by a loge syndrome occurred in one patient with pre-existing severe peripheral vascular disease. This patient needed a fasciotomy. Dislocation of the venous cannula occurred in one patient and required pump removal.

Most pumps were surgically removed in the operation theatre (n=12, 52%). A closure device (Prostar) failed to close completely the puncture site in two out of three attempts. A local compression device (Femostop, Radi) was used successfully in these patients as well as in three others without any complication.

Discussion

PCI in high-risk coronary patients with (potential) haemodynamic compromise is challenging and may be considered a special dimension of MIC¹⁴. Active circulatory support using VAD constitutes a valuable safeguard to reduce potential fatal intervention-related complications in this setting and is nowadays considered part of the current therapeutic armamentarium.

IABP is currently the most widespread mechanical device used to support and to improve coronary perfusion in patients with acute heart failure. IABP implantation is easy and its use is safe and associated with a low incidence of serious complications. IABP may provide a lot of diastolic augmentation of coronary flow in patients in shock, however only limited forward output augmentation (up to about 0.5 L/min). This level of support has been shown to be scarcely able to improve clinical outcomes in the setting of complete hemodynamic collapse¹⁵.

Active circulatory support using VAD combines the beneficial effects of the myocardial unloading and an increase in tissue perfusion pressure. Theoretically, this technique should also allow very early decompression (“unloading”) of the left ventricle in massive myocardial infarction, thus enhancing the chances of recovery of the jeopardised ischemic non-infarcted areas¹⁶. In animal models, LV unloading prior to revascularisation of the infarcted artery led to significant myocardial salvage in comparison to the implementation of unloading after reperfusion¹⁷. A trial testing this concept, including one patient reported in this series, was stopped early due to poor recruitment. However, the concept could nicely be illustrated by echo contrast imaging in one of our patients (Figure 3). Our data confirmed that TandemHeart® support was associated with significant decreases in left ventricular filling pressures, mean systolic blood pressure and arterial pressures. The blood flow, up to 4.5 L/min, provided by the device may be sufficient to both optimally unload the LV and to prevent and even reverse organ dysfunction in cardiogenic shock patients^{11,18}.

Patient selection is the single most crucial factor in determining a successful outcome in patients who receive temporary mechanical circulatory support. Rapid diagnosis and treatment are essential but are often based on limited information. The patient's history and overall clinical setting should be considered in the decision process to initiate circulatory support. A pulmonary-catheter was used in the majority of the patients as adequate left ventricular filling conditions are essential for proper functioning of (preload to) the device. Haemodynamic data, including mixed (or central) venous oxygen saturation, may help to guide the overall patient selection, the




Figure 3. Echocardiographic contrast perfusion images in a patient on Tandemheart(r) circulatory support. Panel A: low circulatory support (high wall tension, poor myocardial perfusion, contrast remains in the ventricle cavity). Panel B: Full circulatory support with effective unloading of the left ventricle (low wall tension, good myocardial perfusion).

management and potential weaning decisions. In our series five patients died while on TandemHeart® support, four were in profound CS at the time of treatment initiation. Patients beyond the 'irreversible edge' of CS prove bad candidates for this technique in this and other series¹¹. Of interest, only three out of the successfully weaned patients (n=17) deceased during the FU period, suggesting this device may have a tangible impact on the long-term survival as well.

The TandemHeart® can be used prophylactically and in bail-out situations. In experienced hands cardiovascular support can be implemented in less than 35 minutes. The most important drawback of the technique may be the difficulty of cannula insertion. Substantial atherosclerotic changes in the iliac artery and the consequent peripheral vascular disease are most common in the group of coronary patients who most benefit from the effects of the pump. In this series, a total number of 13 device-related in-hospital adverse events rate were observed, 0.6 events per patient. Many of them were access site complications. One patient developed a loge syndrome requesting a fasciotomy, but no other iatrogenic arterial insufficiency was noticed. Mild to moderate groin bleeding was encountered in six out of 23 patients (27%), none needing a specific treatment. Most pumps were surgically removed in the operation theatre in order to allow optimal local haemostasis and proper inspection of the cannulation site. A local compression or closure device may provide a valuable alternative given a clean puncture technique. Access site complications may be minimised by prior iliac angiography and by implementing an ultrasound guided puncture technique of the (common) femoral artery¹³. When necessary, the use of distal perfusion and antegrade cannulation of the superficial femoral artery may obviate peripheral vascular complications as well.

As opposed to other modalities of percutaneous cardiopulmonary support (CPS) used in this indication^{19,20}, the TandemHeart® system keeps the patient's lungs as its own ventilator and may be used to support patients for a longer period of time without major

haematological or pulmonary complications (up to 10 days in our series). No thrombi formation was noticed either during the support or after the removal of the device and no thromboembolic event occurred. Evidence of clinical relevant haemolysis, which was a particular concern, was looked for but absent in our series (data not shown). Also this observation is in line with the Leipzig series¹¹.

Systemic hypothermia (<36,5°C) was observed in 25% of patients (n=6) but the temperature decreased below 35°C only in two cases. Contact of the system circuit with room temperature may have contributed to a cooling effect on the blood flowing through the pump. Other factors that may have contributed to heat loss are the use of general anaesthetic agents and neuromuscular blockers leading to vasodilatation and the loss of muscular tonus. Correcting measures may be straightforwardly introduced to avoid or rapidly detect this complication in the future.

Conclusion

Our results create optimism for more widespread introduction of the TandemHeart® in the MIC setting. The concept of closed chest left heart bypass, with left atrial to femoral artery bypass, constitutes a relatively safe and reliable safeguard in high risk PCI and acute heart failure. TandemHeart® PTLA® can be rapidly and percutaneously implanted in the catheterisation suite using standard interventional techniques in both a prophylactic or emergency settings. The TandemHeart® provides more haemodynamic support than the IABP and is effective in stabilising haemodynamics. Access site complications remain the Achilles' heel of this technique.

References

1. Serruys PW, Unger F, Sousa JE, Jatene A, Bonnier HJ, Schonberger JP, Buller N, Bonser R, van den Brand MJ, van Herwerden LA, Morel MA, van Hout BA. Comparison of coronary-artery bypass surgery and stenting for the treatment of multivessel disease. *N Engl J Med.* 2001; 344: 1117-1124.

2. Serruys PW, Ong AT, van Herwerden LA, Sousa JE, Jatene AD, Bonnier JJ, Schonberger JP, Buller N, Bonser R, Disco C, Backx B, Firth B, Unger F. Five year outcomes after coronary stenting versus bypass surgery for the treatment of multivessel disease: the final analysis of the Arterial Revascularization Therapies Study (ARTS) randomized trial. *J Am Coll Cardiol*. 2005; 46: 575-581
3. Topol EJ, Serruys PW. Frontiers in interventional cardiology. *Circulation* 1998;98:1802-1820.
4. de Feyter PJ, McFadden E. Risk score for percutaneous coronary intervention: forewarned is forearmed. *J Am Coll Cardiol* 2003;42(10):1729-30.
5. Singh M, Lennon RJ, Holmes DR, Jr., Bell MR, Rihal CS. Correlates of procedural complications and a simple integer risk score for percutaneous coronary intervention. *J Am Coll Cardiol* 2002;40(3):387-93.
6. Singh M, Rihal CS, Selzer F, Kip KE, Detre K, Holmes DR. Validation of Mayo Clinic risk adjustment model for in-hospital complications after percutaneous coronary interventions, using the National Heart, Lung, and Blood Institute dynamic registry. *J Am Coll Cardiol* 2003;42(10):1722-8.
7. Vanagas G, Kinduris S, Buivydaite K. Assessment of validity for EuroSCORE risk stratification system. *Scand Cardiovasc J* 2005;39(1-2):67-70.
8. Berman M, Stamler A, Sahar G, Georghiou GP, Sharoni E, Brauner R, Medalion B, Vidne BA, Kogan A. Validation of the 2000 Bernstein-Parsonnet score versus the EuroSCORE as a prognostic tool in cardiac surgery. *Ann Thorac Surg* 2006;81(2):537-40.
9. Rao AK, Pratt C, Berke A, Jaffe A, Ockene I, Schreiber TL, Bell WR, Knatterud G, Robertson TL, Terrin ML. Thrombolysis in Myocardial Infarction (TIMI) Trial: phase I: hemorrhagic manifestations and changes in plasma fibrinogen and the fibrinolytic system in patients treated with recombinant tissue plasminogen activator and streptokinase. *J Am Coll Cardiol*. 1988; 11: 1-11.
10. Kawachi Y, Nakashima A, Toshima Y, Arinaga K, Kawano H. Risk stratification analysis of operative mortality in heart and thoracic aorta surgery: comparison between Parsonnet and EuroSCORE additive model. *Eur J Cardiothorac Surg* 2001;20(5):961-6.
11. Thiele H, Lauer B, Hambrecht R, Boudriot E, Cohen HA, Schuler G. Reversal of cardiogenic shock by percutaneous left atrial-to-femoral arterial bypass assistance. *Circulation* 2001;104(24):2917-22.
12. Burkhoff D, O'Neill W, Brunckhorst C, Letts D, Lasorda D, Cohen HA. Feasibility study of the use of the TandemHeart percutaneous ventricular assist device for treatment of cardiogenic shock. *Catheter Cardiovasc Interv* 2006;68(2):211-7.
13. Toumpoulis IK, Anagnostopoulos CE. Can EuroSCORE accurately predict long-term outcome after cardiac surgery? *Nat Clin Pract Cardiovasc Med* 2005;2(12):620-1.
14. Vaina S, Stefanadis C. Treatment of multi-vessel coronary artery disease. What is the optimal revascularisation approach? What do we know, what will we learn? *Hellenic J Cardiol* 2007;48(1):1-4.
15. DeWood MA, Notske RN, Hensley GR, et al. Intraaortic balloon counterpulsation with and without reperfusion for myocardial infarction shock. *Circulation* 1980;61:1105-1112.
16. Laks H, Rosenbranz ER, Buckberg GD. Surgical treatment of cardiogenic shock after myocardial infarction. *Circulation* 1986; 74 (suppl III): 15-22.
17. Meyns B, Stolinski J, Leunens V, et al. Left ventricular support by catheter-mounted axial flow pump reduces infarct size. *J Am Coll Cardiol* 2003;41:1087-1095.
18. Thiele H, Sick P, Boudriot E, Diederich KW, Hambrecht R, Niebauer J, Schuler G. Randomized comparison of intra-aortic balloon support with a percutaneous left ventricular assist device in patients with revascularized acute myocardial infarction complicated by cardiogenic shock. *Eur Heart J* 2005;26(13):1276-83.
19. Shawl FA, Domanski MJ, Wish MH, Davis M. Percutaneous cardiopulmonary bypass support in the catheterization laboratory: technique and complications. *Am Heart J* 1990;120(1):195-203.
20. Vogel RA, Shawl F, Tommaso C, O'Neill W, Overlie P, O'Toole J, Vandormael M, Topol E, Tabari KK, Vogel J and others. Initial report of the National Registry of Elective Cardiopulmonary Bypass Supported Coronary Angioplasty. *J Am Coll Cardiol* 1990;15(1):23-9.